

## Product datasheet for **RC236746**

### **C12orf43 (NM\_001286197) Human Tagged ORF Clone**

#### Product data:

**Product Type:** Expression Plasmids  
**Product Name:** C12orf43 (NM\_001286197) Human Tagged ORF Clone  
**Tag:** Myc-DDK  
**Symbol:** C12orf43  
**Synonyms:** Custos  
**Vector:** pCMV6-Entry (PS100001)  
**E. coli Selection:** Kanamycin (25 ug/mL)  
**Cell Selection:** Neomycin  
**ORF Nucleotide Sequence:** >RC236746 representing NM\_001286197  
**Red=Cloning site Blue=ORF Green=Tags(s)**

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC  
GCC**CGATCGCC**

ATGGCGGCGCCCAGTGGCACAGTGAAGCGATTCCGAAAGTAGTAACAGCAGTAGCGATGCGGAGGAGCTGG  
AGCGGTGCCGCGAGGCGGCAATGCCGGCTTGGGGCTTGGAGCAACGCCCGCACGTGGCAGGGAAGCCAAG  
AGCCGGTCTGCAAATAGCCAGTTGTCAACCTCCAACCGAGCCTCAGCTTCATTACCATCTCAGAAGCA  
GCAAAGGAGCCAGCAAAAGCTAAGGTACAGAAAGTCGCTTGGAGGATGATGGTTTCCGCCCTTTCTTCA  
CATCTGTCCCTGGAGGCCGTGAGAAGGAAGAGTCTCCCAACCCCGCGAAAGCGACAGCCCTCCAGCTC  
CAGCAGTGAAGCAGTGACGAGGAGTGGCGGGTCCCGGAGGCGAGCTGTGTCGGCTCCGACATCCTA  
CAGGAGTCAGCCATCCACAGCCCTGGAACAGTGGAGAAGGAGGCAAGAAGAAAAGGAAGTTGAAAAAGA  
AAGCCAAGAAGGTGGCCAGTGTGACTCGGCTGTGCTGCCACCACCCACCAGCATGGCCACAGTCCA  
GAAGCAGAAGTCAGGTGAGCTCAACGGGGACCAGGTGTGCTTGGGACCAAAAAGAAAAGGCAAAAG  
AAGGCCAGCGAGACCTCTCCATTCCACCAGCAAAGAGTGCTACAGCTATACCTGCAAAC

**ACGCGT**ACGCGGCGCTCGAGCAGAAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT  
ACAAGGATGACGACGATAAGGTTTAA



[View online »](#)

**Protein Sequence:** >RC236746 representing NM\_001286197  
 Red=Cloning site Green=Tags(s)

MAAPSGTVSDSESSNSSSDAEELERCREAAMPWGLEQRPHVAGKPRAGAANSQ LSTSQP SLSFITISEA  
 AKEPAKAKVQKVALEDDGFRLFFTSVPGGREKEESPQPRRKRQPS SSSSESDSEWRRCREAAVSASDIL  
 QESAIHSPTVEKEAKKRKLKKKAKKVASVDSAVAATPTSMATVQKQKSGELNGDQVSLGTKKKKKKAK  
 KASETSPFPKASATAIPAN


TRTRPLEQKLISEEDLAANDILDYKDDDDKV

**Restriction Sites:** Sgfl-MluI

**Cloning Scheme:**



**Plasmid Map:**



**ACCN:** NM\_001286197

**ORF Size:** 690 bp

<b>OTI Disclaimer:</b>	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. <a href="#">More info</a>
<b>OTI Annotation:</b>	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
<b>Components:</b>	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
<b>Reconstitution Method:</b>	<ol style="list-style-type: none"><li>1. Centrifuge at 5,000xg for 5min.</li><li>2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.</li><li>3. Close the tube and incubate for 10 minutes at room temperature.</li><li>4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.</li><li>5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.</li></ol>
<b>RefSeq:</b>	<a href="#">NM_001286197.2</a>
<b>RefSeq Size:</b>	2452 bp
<b>RefSeq ORF:</b>	693 bp
<b>Locus ID:</b>	64897
<b>UniProt ID:</b>	<a href="#">Q96C57</a>
<b>Cytogenetics:</b>	12q24.31
<b>MW:</b>	24.9 kDa
<b>Gene Summary:</b>	Plays a role in the regulation of Wnt signaling pathway during early development. [UniProtKB/Swiss-Prot Function]