

Product datasheet for **RC236784**

PAQR4 (NM_001284511) Human Tagged ORF Clone

Product data:

Product Type: Expression Plasmids
Product Name: PAQR4 (NM_001284511) Human Tagged ORF Clone
Tag: Myc-DDK
Symbol: PAQR4
Vector: pCMV6-Entry (PS100001)
E. coli Selection: Kanamycin (25 ug/mL)
Cell Selection: Neomycin
ORF Nucleotide Sequence: >RC236784 representing NM_001284511
Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCCCGATCGCC

ATGGCGTTCTGGCCGGCCGCGCCTGCTGGACTGGCCAGCTCGCCGCCACCTGCAGTTCAATAAGT
TCGTGCTGACCGGGTACCGGCCCGCCAGCAGCGGCTCGGGCTGCCTGCGCAGCCTTTCTACCTGCACAA
CGAACTGGCAACATCTACACGACGGCTCCGTGCTCTACACCTTTTATGTGCCACCAAGGGGCGAGC
GCTGTGTACGCCCGGCTCCTCGCCCTGGACATGTGTGGGGTCTGCCTTGTCAACACCTTGGGGCCCTGC
CCATCATCCACTGCACCCTGGCCTGCAGGCCCTGGCTGCGCCCGGCTGCCCTGGTGGGCTACACTGTGT
GTCGGGTGTGGCCGGCTGGCGTCTCTACCGCCCTCCACCAGTGTGCGCTCCGGGCATTTGGATGG
CAGGCTGTGCCCGCTACTGGTATTTGGGGCCCGGGAGTGGGTCTGGGTTCCAGGGCTCCAGGCTCCC
TGCCCTGCTACCTGCGCATGGACGCACTGGCGTGTGGGGGACTGGTAAATGTAGCCCGTCTGCCCGA
GCGCTGGGACCTGGCCGCTTTGACTACTGGGGCACTCCACCAGATCATGCACCTGCTGAGCGTGGGC
TCCATCTGCAGCTGCACGCCGGCTCGTGCCCGACCTGCTGCGGCTGCCACCACGCTGTCCCCGGG
AC

ACGCGTACGCGGCCGCTCGAGCAGAAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT
ACAAGGATGACGACGATAAGGTTTAA

Protein Sequence: >RC236784 representing NM_001284511
Red=Cloning site Green=Tags(s)

MAFLAGPRLLDWASSPPLQFNKFVLTGYRPASSGSGCLRSLFYLHNELGNIYTHGSVLVYHLMCHQGGG
AVYARLLALDMCGVCLVNTL GALPIIHCTLACRPWLRPAALVGYTVL SGVAGWRALTAPSTSARLRAFGW
QAAARLLVFGARGVGLGSGAPGSLPCYL RMDALALLGGLV NVARLPERWGPGRFDYWGNSHQIMHLLSVG
SILQLHAGVVPDLLWAAHHACPRD


TRTRPLEQKLI SEEDLAANDILDYKDDDDKV




[View online »](#)

Restriction Sites: SgfI-MluI

Cloning Scheme:



Plasmid Map:



ACCN: NM_001284511

ORF Size: 702 bp

OTI Disclaimer: The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

OTI Annotation: This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.

Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
Reconstitution Method:	<ol style="list-style-type: none">1. Centrifuge at 5,000xg for 5min.2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.3. Close the tube and incubate for 10 minutes at room temperature.4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
RefSeq:	<u>NM_001284511.2</u>
RefSeq Size:	2694 bp
RefSeq ORF:	705 bp
Locus ID:	124222
UniProt ID:	<u>Q8N4S7</u>
Cytogenetics:	16p13.3
Protein Families:	Druggable Genome, Transmembrane
MW:	25.6 kDa