

Product datasheet for **RC236791**

FRA1 (FOSL1) (NM_001300844) Human Tagged ORF Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	FRA1 (FOSL1) (NM_001300844) Human Tagged ORF Clone
Tag:	Myc-DDK
Symbol:	FOSL1
Synonyms:	FRA; fra-1; FRA1
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)
Cell Selection:	Neomycin
ORF Nucleotide Sequence:	>RC236791 representing NM_001300844 Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCC**GCGATCGCC**

ATGTTCCGAGACTTCGGGAACCCGGCCCGAGCTCCGGAAACGGCGGGGTACGGCGCCCGCGCAGC
CCCCGGCCGAGCGCAGGCAGCCAGCAGAAGTCCACCTGGTGCCAAGCATCAACACCATGAGTGGCAG
TCAGGAGCTGCAGTGGATGGTACAGCCTCATTTCTGGGGCCAGCAGTTACCCAGGCCTCTGACCTAC
CCTCAGTACAGCCCCCAACCCGGCCAGGAGTCATCCGGCCCTGGGGCCGCCTCCAGGGTACGTC
GAAGGCCTTGTAACAGGAGACTGACAACTGGAAGATGAGAAATCTGGGCTGCAGCGAGAGATTGAGGA
GCTGCAGAAGCAGAAGGAGCGCCTAGAGCTGGTGTGGAAGCCACCGACCCATCTGCAAAATCCCGGAA
GGAGCCAAGGAGGGGACACAGGCAGTACCAGTGGCACCAGCAGCCACCAGCCCTGCCGCCCTGTAC
CTTGATCTCCCTTTCCCAGGGCCTGTGCTTGAACCTGAGGCACTGCACACCCACACTCATGACCAC
ACCCTCCCTAACTCCTTTACCCCCAGCCTGGTCTTCACTACCCAGCACTCCTGAGCCTGTGCTCA
GCTCATCGCAAGAGTAGCAGCAGCAGCGGAGACCCATCCTCTGACCCCTTGGCTCTCAACCCCTCTCG
CTTTG

ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT
ACAAGGATGACGACGATAAGGTTAA



[View online »](#)

OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
Reconstitution Method:	<ol style="list-style-type: none">1. Centrifuge at 5,000xg for 5min.2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.3. Close the tube and incubate for 10 minutes at room temperature.4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
RefSeq:	NM_001300844.2
RefSeq Size:	1651 bp
RefSeq ORF:	708 bp
Locus ID:	8061
UniProt ID:	P15407
Cytogenetics:	11q13.1
Protein Families:	Druggable Genome, Transcription Factors
Protein Pathways:	Wnt signaling pathway
MW:	25.5 kDa
Gene Summary:	The Fos gene family consists of 4 members: FOS, FOSB, FOSL1, and FOSL2. These genes encode leucine zipper proteins that can dimerize with proteins of the JUN family, thereby forming the transcription factor complex AP-1. As such, the FOS proteins have been implicated as regulators of cell proliferation, differentiation, and transformation. Several transcript variants encoding different isoforms have been found for this gene. [provided by RefSeq, Jul 2014]