

## Product datasheet for **RC236866**

### Cyclin D3 (CCND3) (NM\_001287427) Human Tagged ORF Clone

#### Product data:

**Product Type:** Expression Plasmids  
**Product Name:** Cyclin D3 (CCND3) (NM\_001287427) Human Tagged ORF Clone  
**Tag:** Myc-DDK  
**Symbol:** CCND3  
**Vector:** pCMV6-Entry (PS100001)  
**E. coli Selection:** Kanamycin (25 ug/mL)  
**Cell Selection:** Neomycin  
**ORF Nucleotide Sequence:** >RC236866 representing NM\_001287427  
**Red=Cloning site Blue=ORF Green=Tags(s)**

TTTTGTAATACGACTCACTATAGGGCGGCCGGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC  
GCC**CGATCGCC**

ATGAAAGTGCAGGAGGCAGCGGCAGGAGCCTCCTCCCTCCCTCGGGTATGTGAGGAGCAGCGCTGTG  
AGGAGGAAGTCTCCCCCTGGCCATGAACTACCTGGATCGTACCTGTCTTGCCTCCCAACCGAAAGGC  
GCAGTTGCAGCTCTGGGTGCGGTCTGCATGCTGCTGGCCTCCAAGCTGCGCGAGACCACGCCCTGACC  
ATCGAAAACTGTGCATCTACACCGACCAGCTGTCTCTCCCGCCAGTTGCGGGACTGGGAGGTGCTGG  
TCCTAGGGAAGCTCAAGTGGGACCTGGCTGCTGTGATTGCACATGATTTCTGGCCTTATTCTGCACCG  
GCTCTCTCTGCCCGTGACCGACAGGCCTTGGTCAAAAAGCATGCCAGACCTTTTTGGCCCTCTGTGCT  
ACAGATTACCTTTGCCATGTACCCGCCATCCATGATCGCCACGGGCAGCATTGGGGCTGCAGTCAAG  
GCCTGGGTGCCTGCTCCATGTCCGGGATGAGCTCACAGAGCTGCTGGCAGGGATCACTGGCACTGAAGT  
GGACTGCCTGCGGCCTGTCAGGAGCAGATCGAAGCTGCACTCAGGGAGAGCCTCAGGGAAGCCTCTCAG  
ACCAGTCCAGCCAGCGCCAAAGCCCCCGGGGCTCCAGCAGCAAGGGCCAGCCAGACCAGCACTC  
CTACAGATGTCACAGCCATACACCTG

**ACGCGT**ACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT  
ACAAGGATGACGACGATAAGGTTTAA

**Protein Sequence:** >RC236866 representing NM\_001287427  
**Red=Cloning site Green=Tags(s)**

MKSAGGSGRSLLPSPRVCEEQRCEEEVFLAMNYLDRYLSCVPTRKAQLQLLGAVCMLLASKLRETTPLT  
IEKLCIYTDHAVSPRQLRDWEVLVLGKWKDLAAVIAHDFLAFILHRLSLPRDRQALVKKHAQTFLALCA  
TDYTFAMYPPSMIATGSIGAAVQGLGACSMGDELTELLAGITGTEVDCLRACQEQIEAALRESLREASQ  
TSSSPAPKAPRGSSSQGPSQSTPTDVTAIHL


**TR**TRPLEQKLISEEDLAANDILDYKDDDDKV




[View online »](#)

Restriction Sites: SgfI-MluI

Cloning Scheme:



Plasmid Map:



ACCN: NM\_001287427

ORF Size: 726 bp

OTI Disclaimer: The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

OTI Annotation: This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.

<b>Components:</b>	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
<b>Reconstitution Method:</b>	<ol style="list-style-type: none"><li>1. Centrifuge at 5,000xg for 5min.</li><li>2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.</li><li>3. Close the tube and incubate for 10 minutes at room temperature.</li><li>4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.</li><li>5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.</li></ol>
<b>RefSeq:</b>	<u><a href="#">NM_001287427.1</a></u> , <u><a href="#">NP_001274356.1</a></u>
<b>RefSeq Size:</b>	2585 bp
<b>RefSeq ORF:</b>	729 bp
<b>Locus ID:</b>	896
<b>UniProt ID:</b>	<u><a href="#">P30281</a></u>
<b>Cytogenetics:</b>	6p21.1
<b>Protein Families:</b>	Druggable Genome
<b>Protein Pathways:</b>	Cell cycle, Focal adhesion, Jak-STAT signaling pathway, p53 signaling pathway, Wnt signaling pathway
<b>MW:</b>	26.7 kDa
<b>Gene Summary:</b>	The protein encoded by this gene belongs to the highly conserved cyclin family, whose members are characterized by a dramatic periodicity in protein abundance through the cell cycle. Cyclins function as regulators of CDK kinases. Different cyclins exhibit distinct expression and degradation patterns which contribute to the temporal coordination of each mitotic event. This cyclin forms a complex with and functions as a regulatory subunit of CDK4 or CDK6, whose activity is required for cell cycle G1/S transition. This protein has been shown to interact with and be involved in the phosphorylation of tumor suppressor protein Rb. The CDK4 activity associated with this cyclin was reported to be necessary for cell cycle progression through G2 phase into mitosis after UV radiation. Several transcript variants encoding different isoforms have been found for this gene. [provided by RefSeq, Oct 2008]