

Product datasheet for **RC236897**

TMEM116 (NM_001294314) Human Tagged ORF Clone

Product data:

Product Type: Expression Plasmids
Product Name: TMEM116 (NM_001294314) Human Tagged ORF Clone
Tag: Myc-DDK
Symbol: TMEM116
Vector: pCMV6-Entry (PS100001)
E. coli Selection: Kanamycin (25 ug/mL)
Cell Selection: Neomycin
ORF Nucleotide Sequence: >RC236897 representing NM_001294314
Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCC**CGATCGCC**

ATGAAACACACCCAGAGTGGACAGAGCACATCTCCACTGGTGATAGATTATACTTGTGCGAGTTTGTCAA
TGGCCTTTGTTTTCTCAAGCCTGATACCTCTGCTATTGATGACACCTGTATTCTGTCTGGGAAATACTAG
TGAATGTTTCCAAAACCTCAGTCAGAGCCACAAGTGTATCTTGATGCACTCACCACCATCAGCCATGGCT
GAACTTCCACCTTCTGCCAACACATCTGTCTGTAGCACACTTTATTTTTATGGTATCGCCATTTTCTGG
GCAGCTTTGTAAGTCAAGTCAAGTCAAGTCAAGTCAAGTCAAGTCAAGTCAAGTCAAGTCAAGTCAAGTCAAGT
TGTGAAGTCAACTGGCTTTCTGGGGAGTGAACAGTGGGAGTGGGAGTGGGAGTGGGAGTGGGAGTGGGAGT
TTCTACCCAGTGGCCTTCTTTGCTGCTGGGGCCAGCTGTCATTCTAATGATCATAAAGCTGACTAAGC
CACAGGACACCAAGCTTACATGGCCCTTTATGTTCTCCAGGCTCTAACGGCAACATCTCAGGGTCTACT
CAACTGTGGAGTATATGGCTGGACGCAGCACAAATCCACCAACTAAAGCAGGAGGCTCGGCGTGATGCA
GATACCCAGACACCATTATTATGCTCACAGAAGAGATTCTATAGCAGGGGCTTAAATTCAGTGAATCCA
CCCTGACTTTTCTGCCAGTACTTCTACCATTTT

ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT
ACAAGGATGACGACGATAAGGTTTAA

Protein Sequence: >RC236897 representing NM_001294314
Red=Cloning site Green=Tags(s)


MKHTQSGQSTSPLVIDYTCRVCQMAFVFSLLIPLLLMTPVFCLGNTSECFQNFSSQSHKCILMHSPPSAMA
ELPPSANTSVCSLYFYGIAIFLGSFVLSLLTIMVLLIRAQTLYKFKVSTGFLGSEQWAVIHIIVDQVRV
FYPVAFCCWGPVILMIKLTQPDKLHMLYVLQAL TATSQGLLNCGVYGTQHKFHQLKQEARRDA
DTQTPLLSQKRFYSRGLNSLESLTFPASTSTIF

TRTRPLEQKLISEEDLAANDILDYKDDDDKV




Restriction Sites: SgfI-MluI

Cloning Scheme:



Plasmid Map:



ACCN: NM_001294314

ORF Size: 735 bp

OTI Disclaimer: The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

OTI Annotation: This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.

Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
Reconstitution Method:	<ol style="list-style-type: none">1. Centrifuge at 5,000xg for 5min.2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.3. Close the tube and incubate for 10 minutes at room temperature.4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
RefSeq:	<u>NM_001294314.1</u> , <u>NP_001281243.1</u>
RefSeq Size:	1824 bp
RefSeq ORF:	738 bp
Locus ID:	89894
UniProt ID:	<u>Q8NCL8</u>
Cytogenetics:	12q24.12-q24.13
Protein Families:	Druggable Genome, Transmembrane
MW:	27.9 kDa