

Product datasheet for **RC236924**

TMEM40 (NM_001284406) Human Tagged ORF Clone

Product data:

Product Type: Expression Plasmids
Product Name: TMEM40 (NM_001284406) Human Tagged ORF Clone
Tag: Myc-DDK
Symbol: TMEM40
Vector: pCMV6-Entry (PS100001)
E. coli Selection: Kanamycin (25 ug/mL)
Cell Selection: Neomycin
ORF Nucleotide Sequence: >RC236924 representing NM_001284406
Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCC**CGATCGCC**

ATGGCCGAGGCCGGAGCCGGCTCCCTCAGCCTTCCGGGAGGAAAAGCCATGGAGACTTCAGCATCCTCCT
CCCAGCCTCAGGACAACAGTCAAGTCCACAGAGAAACAAGATGTAGACTATGGAGAGACAGATTTCCA
CAAGCAAGATGGGAAGGCTGGACTTTTTCCCAAGAACAATATGAGAGAAACAAGTCTTCTCCTCCTCC
TCCTTCTCCTCATCCTCCTCATCTTTCATCCTCCTCCTCCTCAGAGAGCAATGATGAGGACCAGC
AACCCAGAGCAACCGGAAAACATCGACGGAGCCTGGGGCTGGATACCCACGGGAACGGCTCACCCGG
TCCTGGGCATGGGGAGCCTGACGTTTTGAAGGATGAGCTTCAACTCTATGGAGATGCTCCTGGAGAGGTG
GTACCCTCTGGGGAATCAGGACTCCGAAGGAGAGGCTCTGACCCAGCAAGTGGAGAAGTGGAGGCCTCTC
AGTTAAGAAGACTGAATATAAAGAAAGATGATGAGTTTTTCCATTTTCGTCCTCCTGTGCTTTGCCATCGG
GGCCTTGCTGGTGTGTTACTATTACGCAGACTGGTTCATGTCTCTTGGGGTCGGCCTGCTCACCTTC
GCCTCCCTGGAACCGTTGGCATCTACTTCGGACTAGTGTACCGTATCCACAGCGTCTCCAAGGCTTCA
TCCCCCTTCCAGAAGTTTAGGCTGACAGGGTTCAGGAAGACTGAC

ACGCGTACGCGGCCGCTCGAGCAGAAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT
ACAAGGATGACGACGATAAGGTTTAA

Protein Sequence: >RC236924 representing NM_001284406
Red=Cloning site Green=Tags(s)

MAEAGAGSLSLPGGKAMETSASSSQPDNSQVHRETEDVDYGETDFHKQDGKAGLFSQEYERNKSSSSSS
SSSSSSSSSSSSSESNDEDQPRATGKHRRSLGAGYPHNGSPGPGHGPDLKDELQLYGDAPGEV
VPSGESGLRRRSDPASGEVEASQLRRLNIKKDDEFFHFVLLCFAIGALLVCYHYADWFMSLGVLLTF
ASLETVGIYFGLVYRIHSVLQGFIPLFQKFRLTGFRKTD


TRTRPLEQKLISEEDLAANDILDYKDDDDKV




[View online »](#)

Restriction Sites: SgfI-MluI

Cloning Scheme:



Plasmid Map:



ACCN: NM_001284406

ORF Size: 747 bp

OTI Disclaimer: The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

OTI Annotation: This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.

| | |
|-------------------------------|---|
| Components: | The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water). |
| Reconstitution Method: | <ol style="list-style-type: none">1. Centrifuge at 5,000xg for 5min.2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.3. Close the tube and incubate for 10 minutes at room temperature.4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C. |
| RefSeq: | <u>NM_001284406.1</u> , <u>NP_001271335.1</u> |
| RefSeq Size: | 2161 bp |
| RefSeq ORF: | 750 bp |
| Locus ID: | 55287 |
| UniProt ID: | <u>Q8WWA1</u> |
| Cytogenetics: | 3p25.2 |
| Protein Families: | Transmembrane |
| MW: | 27.3 kDa |