

## Product datasheet for **RC236950**

### **C22orf25 (TANGO2) (NM\_001283154) Human Tagged ORF Clone**

#### Product data:

Product Type:	Expression Plasmids
Product Name:	C22orf25 (TANGO2) (NM_001283154) Human Tagged ORF Clone
Tag:	Myc-DDK
Symbol:	TANGO2
Synonyms:	C22orf25; MECRCN
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)
Cell Selection:	Neomycin
ORF Nucleotide Sequence:	>RC236950 representing NM_001283154 Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC  
GCC**CGGATCGC**C

ATGTGCATCATCTTCTTTAAGTTTGATCCTCGCCCTGTTTCCAAAACGCGTACAGGCTCATCTTGGCAG  
CCAACAGGGATGAATTCTACAGCCGACCCTCCAAGTTAGCTGACTTCTGGGGGAACAACAACGAGATCCT  
CAGTGGGCTGGACATGGAGGAAGGCAAGGAAGGAGGCACATGGCTGGGCATCAGCACACGTGGCAAGCTG  
GCAGCACTCACCAACTACCTGCAGCCGAGCTGGACTGGCAGGCCGAGGGCGAGGTGAACCTGTACACC  
ACTTTCTGACCACTGACGTGGACAGCTTGTCTACCTGAAGAAGGTCTCTATGGAGGGCCATCTGTACAA  
TGGCTTCAACCTCATAGCAGCCGACCTGAGCACAGCAAAGGGAGACGTCAATTTGCTACTATGGGAACCGA  
GGGGAGCCTGATCCTATCGTTTTGACGCCAGGCACCTACGGGCTGAGCAACGCGCTGCTGGAGACTCCCT  
GGAGGAAGCTGTGCTTTGGGAAGCAGCTCTTCTGGAGGCTGTGGAACGGAGCCAGGCGCTGCCAAGGA  
TGTGCTCATCGCCAGCCTCCTGGATGTGCTCAACAATGAAGAGGCAGTGCCAGACCCGGCCATCGAGGA  
CCAGGGTGGGGAGTACGTGCAGCCCATGCTGAGCAAGTACGCGGCTGTGTGCGTGCCTGCCCTGGCTAC  
GGCACCAGAACCAACTATCATCTCTGGTAGATGCGGACGGCCACGTGACCTTAC

AG**CGGACCG**ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCC  
TGGATTACAAGGATGACGACGATAAGGTTTAA



[View online »](#)


**Protein Sequence:** >RC236950 representing NM\_001283154  
Red=Cloning site Green=Tags(s)

MCIFFKFDPRPVSKNAYRLILAAANRDEFYSRPSKLADFWGNNNEILSGLDMEEGKEGGTTLGISTRGKL  
 AALTNYLQPQLDWQARGRGLVTHFLTDDVDSL SYLKKVSMEGHLYNGFNLIAADLSTAKGDVICYYGNR  
 GEPDPIVLTPGTYGLSNALLETPWRKLCFGKQLFLEAVERSQUALPKDVLIASLLDVLNNEEAAARPGHRG  
 PGWGVRAAAHEQVRGCVRALPWL RHQNQHYHPGRRCGRPRDLH


SGPTRTRPLEQKLI SEEDLAANDILDYKDDDDKV

**Restriction Sites:** SgfI-RsrII

**Cloning Scheme:**



**Plasmid Map:**



**ACCN:** NM\_001283154  
**ORF Size:** 756 bp

<b>OTI Disclaimer:</b>	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. <a href="#">More info</a>
<b>OTI Annotation:</b>	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
<b>Components:</b>	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
<b>Reconstitution Method:</b>	<ol style="list-style-type: none"><li>1. Centrifuge at 5,000xg for 5min.</li><li>2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.</li><li>3. Close the tube and incubate for 10 minutes at room temperature.</li><li>4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.</li><li>5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.</li></ol>
<b>RefSeq:</b>	<a href="#">NM_001283154.2</a> , <a href="#">NP_001270083.1</a>
<b>RefSeq Size:</b>	2327 bp
<b>RefSeq ORF:</b>	759 bp
<b>Locus ID:</b>	128989
<b>UniProt ID:</b>	<a href="#">Q6ICL3</a>
<b>Cytogenetics:</b>	22q11.21
<b>MW:</b>	28.7 kDa
<b>Gene Summary:</b>	This gene belongs to the transport and Golgi organization family, whose members are predicted to play roles in secretory protein loading in the endoplasmic reticulum. Depletion of this gene in <i>Drosophila</i> S2 cells causes fusion of the Golgi with the ER. In mouse tissue culture cells, this protein co-localizes with a mitochondrially targeted mCherry protein and displays very low levels of co-localization with Golgi and peroxisomes. Allelic variants of this gene are associated with rhabdomyolysis, metabolic crises with encephalopathy, and cardiac arrhythmia. [provided by RefSeq, Apr 2016]