

## Product datasheet for **RC236968**

### LILRB4 (NM\_001278430) Human Tagged ORF Clone

#### Product data:

**Product Type:** Expression Plasmids  
**Product Name:** LILRB4 (NM\_001278430) Human Tagged ORF Clone  
**Tag:** Myc-DDK  
**Symbol:** LILRB4  
**Synonyms:** CD85K; ILT-3; ILT3; LIR-5; LIR5  
**Vector:** pCMV6-Entry (PS100001)  
**E. coli Selection:** Kanamycin (25 ug/mL)  
**Cell Selection:** Neomycin  
**ORF Nucleotide Sequence:** >RC236968 representing NM\_001278430  
**Red=Cloning site Blue=ORF Green=Tags(s)**

TTTTGTAATACGACTCACTATAGGGCGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC  
GCC**CGATCGCC**

ATGATCCACCTTCACGGCTCTGCTCTGCCTCGGGCTGAGTCTGGGCCCCAGGACCCACATGCAGGCAG  
GGCCCCCCCCAAACCCACCTCTGGGCTGAGCCAGGCTCTGTGATCAGCTGGGGAACTCTGTGACCAT  
CTGGTGTGAGGGACCTGGAGGCTCGGGAGTACCGTCTGGATAAAGAGGAAAGCCAGCACCCCTGGGAC  
AGACAGAACCCACTGGAGCCCAAGAACAAGGCCAGATTCTCCATCCCATCCATGACAGAGGACTATGCAG  
GGAGATACCGCTGTTACTATCGCAGCCCTGTAGGCTGGTACAGCCCAGTGACCCCTGGAGCTGGTGTGAT  
GACAGGAGCCTACAGTAAACCCACCTTTAGCCCTGCCGAGTCTCTTGTGACCTCAGGAAAGAGCGTG  
ACCCTGCTGTGTGATCAGTACGGAGCCCAATGGACTTTTCTTCTGATCAAGGAGCGGGCAGCCCATCCCC  
TACTGCATCTGAGATCAGAGCACGGAGCTCAGCAGCACCAGGCTGAATCCCCATGAGTCTGTGACCTC  
AGTGCACGGGGGACCTACAGGTGCTTCAGCTCACACGGCTTCTCCACTACCTGCTGTACACCCCACT  
GACCCCTGGAGCTCATAGTCTCAGGATCCTTGGAGGGTCCCAGGCCCTACCCACAAGGTCCGCTCAA  
CAGCTGCAGGCCCTGAGGACCAGCCCTCATGCCTACAGGGTCACTCCCCACAGTGGTGAG

**ACGCGT**ACGCGGCCGCTCGAGCAGAAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT  
ACAAGGATGACGACGATAAGGTTTAA



[View online »](#)

Protein Sequence: >RC236968 representing NM\_001278430  
 Red=Cloning site Green=Tags(s)

MIPTFTALLCLGLSLGPRTHMQAGPLPKPTLWAEPGSVISWNSVTIWCQGTLEAREYRLDKEESPAPWD  
 RQNPLEPKNKARFSIPSMTEYAGRYRCYRSPVGWSQPSDPLELVMGTGAYSKPTLSALPSPLVTSKSV  
 TLLCQSRSPMDTFLLIKERAAHPLLHLRSEHGAQQHQAEFPMSPVTSVHGGTYRCFSSHGFSHYLLSHPS  
 DPLELIVSGSLEGPRPSPTRSVSTAAGPEDQPLMPTGSVPHSGE


TRTRPLEQKLISEEDLAANDILDYKDDDDKV

Restriction Sites: SgfI-MluI

Cloning Scheme:



Plasmid Map:



ACCN: NM\_001278430  
 ORF Size: 762 bp

<b>OTI Disclaimer:</b>	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. <a href="#">More info</a>
<b>OTI Annotation:</b>	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
<b>Components:</b>	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
<b>Reconstitution Method:</b>	<ol style="list-style-type: none"><li>1. Centrifuge at 5,000xg for 5min.</li><li>2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.</li><li>3. Close the tube and incubate for 10 minutes at room temperature.</li><li>4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.</li><li>5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.</li></ol>
<b>RefSeq:</b>	<a href="#">NM_001278430.3</a> , <a href="#">NP_001265359.2</a>
<b>RefSeq Size:</b>	999 bp
<b>RefSeq ORF:</b>	765 bp
<b>Locus ID:</b>	11006
<b>UniProt ID:</b>	<a href="#">Q8NHJ6</a>
<b>Cytogenetics:</b>	19q13.42
<b>Protein Families:</b>	Transmembrane
<b>MW:</b>	28.2 kDa
<b>Gene Summary:</b>	This gene is a member of the leukocyte immunoglobulin-like receptor (LIR) family, which is found in a gene cluster at chromosomal region 19q13.4. The encoded protein belongs to the subfamily B class of LIR receptors which contain two or four extracellular immunoglobulin domains, a transmembrane domain, and two to four cytoplasmic immunoreceptor tyrosine-based inhibitory motifs (ITIMs). The receptor is expressed on immune cells where it binds to MHC class I molecules on antigen-presenting cells and transduces a negative signal that inhibits stimulation of an immune response. The receptor can also function in antigen capture and presentation. It is thought to control inflammatory responses and cytotoxicity to help focus the immune response and limit autoreactivity. Multiple transcript variants encoding different isoforms have been found for this gene. [provided by RefSeq, Jul 2008]