

## Product datasheet for **RC237006**

### Chimaerin 2 (CHN2) (NM\_001293081) Human Tagged ORF Clone

#### Product data:

**Product Type:** Expression Plasmids  
**Product Name:** Chimaerin 2 (CHN2) (NM\_001293081) Human Tagged ORF Clone  
**Tag:** Myc-DDK  
**Symbol:** CHN2  
**Synonyms:** ARHGAP3; BCH; CHN2-3; RHOGAP3  
**Vector:** pCMV6-Entry (PS100001)  
**E. coli Selection:** Kanamycin (25 ug/mL)  
**Cell Selection:** Neomycin  
**ORF Nucleotide Sequence:** >RC237006 representing NM\_001293081  
**Red=Cloning site Blue=ORF Green=Tags(s)**

TTTTGTAATACGACTCACTATAGGGCGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC  
GCC**GCGATCGCC**

ATGCGAGTCATTATCGTGGTCTCTGGCACCCAGTCCACACGTTCCGAGGCCACACTGGTGTGAATATT  
GTGCCAATTTTCATGTGGGGGCTCATCGCCCAAGGGTCCGGTGCTCAGACTGTGGATTGAACGTACACAA  
ACAGTGTCCAAGCACGTTCCCAATGACTGCCAACCTGATCTCAAGAGGATCAAGAAAGTGTACTGTTGT  
GACCTCACAACACTTGTGAAGGCTCACAACACTCAGAGACCCATGGTGGTAGACATATGCATTCCGGAAA  
TTGAAGCAAGAGGATTAATAATCGGAAGGCCTTTACAGAGTCTCTGGGTTCACTGAACACATTGAAGATGT  
CAAAATGGCATTGACAGAGATGGTAAAAGGCCGATATATCTGCCAATGTCTATCCAGACATAAACATC  
ATCACTGGAGCCCTTAACTGTATTTTCAGAGACTTACCCATCCCTGTCATCACATATGATACCTATTCCA  
AATTTATAGATGCAGCAAAAATCTCCAATGCAGATGAGAGGCTGGAAGCCGTCCATGAAGTGCTGATGCT  
GCTGCCTCCTGCCACTATGAAACCCCTCCGGTACCTAATGATCCACCTCAAAAAGTTACTATGAATGAA  
AAAGACAATTTTCATGAATGCAGAAAATCTGGGGATCGTGTGGGGCCACTCTGATGAGGCCCTGAGG  
ACAGCACCTGACCACCTGCATGATATGCGGTACCAAAAGCTGATTGTGCAGATTTTAATAGAAAACGA  
AGACGTTTTATTTC

**ACGCGTACGCGGCCGCTCGAGCAGAAA**ACTCATCTCAGAAGAGGATCTGGCAGCAATGATATCCTGGATT  
ACAAGGATGACGACGATAAGGTTTAA



[View online »](#)

Protein Sequence: >RC237006 representing NM\_001293081  
 Red=Cloning site Green=Tags(s)

MRVTYRGLWHPVHTFRGPHWCEYCANFMWGLIAQGVRCSDCGLNVHKQCSKHVPNDQCQDLKRIKKVYCC  
 DLTTLVKAHNTQRPMVVDICIREIEARGLKSEGLYRVSGFTEHIEDVKMAFDRDGEKADISANVYPDINI  
 ITGALKLYFRDLPVITYDYTSKFIIDAAKISNADERLEAVHEVLMMLPPAHYETLRYLMIHLKKVTMNE  
 KDNFMNAENLGIIVFGPTLMRPPEDSTLTLHDMRYQKLIVQILINEDVLF


TRTRPLEQKLISEEDLAANDILDYKDDDDKV

Restriction Sites: SgfI-MluI

Cloning Scheme:



Plasmid Map:



ACCN: NM\_001293081  
 ORF Size: 783 bp

<b>OTI Disclaimer:</b>	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. <a href="#">More info</a>
<b>OTI Annotation:</b>	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
<b>Components:</b>	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
<b>Reconstitution Method:</b>	<ol style="list-style-type: none"><li>1. Centrifuge at 5,000xg for 5min.</li><li>2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.</li><li>3. Close the tube and incubate for 10 minutes at room temperature.</li><li>4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.</li><li>5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.</li></ol>
<b>RefSeq:</b>	<a href="#">NM_001293081.1</a> , <a href="#">NP_001280010.1</a>
<b>RefSeq Size:</b>	2530 bp
<b>RefSeq ORF:</b>	786 bp
<b>Locus ID:</b>	1124
<b>Cytogenetics:</b>	7p14.3
<b>MW:</b>	30.6 kDa
<b>Gene Summary:</b>	This gene encodes a guanosine triphosphate (GTP)-metabolizing protein that contains a phorbol-ester/diacylglycerol (DAG)-type zinc finger, a Rho-GAP domain, and an SH2 domain. The encoded protein translocates from the cytosol to the Golgi apparatus membrane upon binding by diacylglycerol (DAG). Activity of this protein is important in cell proliferation and migration, and expression changes in this gene have been detected in cancers. A mutation in this gene has also been associated with schizophrenia in men. Alternative transcript splicing and the use of alternative promoters results in multiple transcript variants. [provided by RefSeq, May 2014]