

Product datasheet for **RC237067**

CLK2 (NM_001294339) Human Tagged ORF Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	CLK2 (NM_001294339) Human Tagged ORF Clone
Tag:	Myc-DDK
Symbol:	CLK2
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)
Cell Selection:	Neomycin
ORF Nucleotide Sequence:	>RC237067 representing NM_001294339 Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCC**CGATCGC**C

ATGTTTGACTGGTTTGACTACCATGGCCACATGTGTATCTCCTTTGAGCTTCTGGGCCTTAGCACCTTCG
ATTTCTCAAAGACAACAACTACCTGCCCTACCCATCCACCAAGTGCGCCACATGGCCTTCCAGCTGTG
CCAGGCTGTCAAGTTCTCCATGATAACAAGCTGACACATACAGACCTCAAGCCTGAAAAATTCTGTTT
GTGAATTCAGACTATGAGCTCACCTACAACCTAGAGAAGAAGCGAGATGAGCGCAGTGTGAAGAGCACAG
CTGTGCGGGTGGTAGACTTTGGCAGTGCCACCTTTGACCATGAGCACCATAGCACCATTGTCTCCACTCG
CCATTACCGAGCACCAGAAGTCATCCTTGAGTTGGGCTGGTCACAGCCTTGTGATGTGTGGAGTATAGGC
TGCATCATCTTTGAATACTATGTGGGATTACCCTCTTCCAGACCCATGACAACAGAGAGCATCTAGCCA
TGATGAAAGGATCTGGGTCCTATCCCTTCCCGGATGATCCGAAAGACAAGAAAGCAGAAATATTTTAA
CCGGGGTCGCCTGGATTGGGATGAGAACACATCAGCTGGGCGCTATGTTTCGTGAGAACTGCAAACCGCTG
CGGCGGTATCTGACCTCAGAGGCAGAGGAACACCACAGCTCTTCGATCTGATTGAAAGCATGCTAGAGT
ATGAACCAGCTAAGCGGCTGACCTTGGGTGAAGCCCTTCAGCATCTTTCTTCGCCCGCCTTCGGGCTGA
GCCGCCAACAAGTTGTGGGACTCCAGTCGGGATATCAGTCGG

ACGCGTACGCGGCCGCTCGAGCAGAAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT
ACAAGGATGACGACGATAAGGTTTAA



[View online »](#)

Protein Sequence: >RC237067 representing NM_001294339
Red=Cloning site Green=Tags(s)

MFDWFDYHGHMCISFELLGLSTFDLKDNNYLPYPIHQVRHMAFQLCQAVKFLHDNKLTHDLPENILF
 VNSDYELTYNLEKKRDRERSVKSTAVRVVDFGSATFDHEHHSTIVSTRHYRAPEVILELGWSQPCDVWSIG
 CIIFEYVYVGFTLFQTHDNREHLAMMERILGPIPSRMIRKTRKQKYFYRGRLDWDENTSAGRYVRENCKPL
 RRYLTSEAEHHQLFDLIESMLEYEPKRLTLGEALQHPFFARLRAEPPNKLWDSSRDISR


TRTRPLEQKLISEEDLAANDILDYKDDDDKV

Restriction Sites: SgfI-MluI

Cloning Scheme:



Plasmid Map:



ACCN: NM_001294339

ORF Size: 813 bp

OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
Reconstitution Method:	<ol style="list-style-type: none">1. Centrifuge at 5,000xg for 5min.2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.3. Close the tube and incubate for 10 minutes at room temperature.4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
RefSeq:	NM_001294339.2
RefSeq Size:	2115 bp
RefSeq ORF:	816 bp
Locus ID:	1196
UniProt ID:	P49760
Cytogenetics:	1q22
Protein Families:	Druggable Genome, Protein Kinase
MW:	32.8 kDa
Gene Summary:	This gene encodes a dual specificity protein kinase that phosphorylates serine/threonine and tyrosine-containing substrates. Activity of this protein regulates serine- and arginine-rich (SR) proteins of the spliceosomal complex, thereby influencing alternative transcript splicing. Chromosomal translocations have been characterized between this locus and the PAFAH1B3 (platelet-activating factor acetylhydrolase 1b, catalytic subunit 3 (29kDa)) gene on chromosome 19, resulting in the production of a fusion protein. Note that this gene is distinct from the TELO2 gene (GeneID:9894), which shares the CLK2 alias, but encodes a protein that is involved in telomere length regulation. There is a pseudogene for this gene on chromosome 7. Alternative splicing results in multiple transcript variants. [provided by RefSeq, Jun 2014]