

Product datasheet for **RC237120**

DBNDD1 (NM_001288709) Human Tagged ORF Clone

Product data:

Product Type: Expression Plasmids
Product Name: DBNDD1 (NM_001288709) Human Tagged ORF Clone
Tag: Myc-DDK
Symbol: DBNDD1
Vector: pCMV6-Entry (PS100001)
E. coli Selection: Kanamycin (25 ug/mL)
Cell Selection: Neomycin
ORF Nucleotide Sequence: >RC237120 representing NM_001288709
Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCC**CGATCGCC**

ATGACCCAGCAGTGGAGGGAAGTGAATGAGGACAGAGAAGCAGTGGCCAGGCCAGGAGCTGGCAATCGC
TGAATCTGGGCGATAGGGACATGGGGTCACTTTCATCGAAACTGCCTTGGCGCTCTCTGCAGCCCGGGC
CGCCCGCCTCTCCTTACCTCCATCCCCATCCCGCCCTTCCAGCCTCCCAAGCACACGCCCTCGCCC
CCACCCGCCCCCGGTTAATTCGCGCCATCCAGGGCGCTCTGTGAGTCTGTATTACCCAGCGCGAGTC
AGGCCCTCTGTGCCCTCGCCAGCCCGAGGGCGCCCGCACACAGTCGGCGCCCCACACGCAAGTAC
CTGGAAGAGAACAAGAGGACGCCAGCGGGCGGGCCGCAAGATCGTTAAGGAGGCTGAGGTGCCGACG
GCTGCGCTGGGCGTCCCAGCCAGGGGACAGGGGACAATGGCCACACGCCTGTGGAGGAGGAGGTGGGG
GCATCCCAGTACCAGCACCGGGCTCCTGCAGGTACGGAGAGGAGGCAGCCTCTGAGCAGCGTCTCCTC
TCTGGAGGTCCACTTCGACCTCCTGGACCTCACTGAGCTACCGACATGTCGGACCAGGAGCTGGCCGAG
GTCTTTGCTGACTCGGACGACGAGAACCTCAACACCGAGTCCCCAGCAGGTCTGCACCCGCTGCCCGGG
CCGGCTACCTGCGCTCCCTTCTGGACGAGGACAAGGGCTGAGCAGAGCCACGAGAAGCAGCCCTAGG
CGACCCCGAGCGGACGGCCACAGTCCTGGACACGTTTCTCACTGTGGAGAGGCCCCAGGAGGAC

ACGCGTACGCGGCCGCTCGAGCAGAAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT
ACAAGGATGACGACGATAAGGTTTAA



[View online »](#)


Protein Sequence: >RC237120 representing NM_001288709
 Red=Cloning site Green=Tags(s)

MTQQWREVNEDREAVAQARSWQSLNLGDRDMGVSF IETALALSAARARRPLLTSHPSPRPFQPPKHTPSP
 PPAPRVNSRHPGRSVSLSPASQAPPVPLASPEGRPHTVGAPHASTWKRTRGRQRARPPKIVKEAEVPQ
 AALGVPAQGTGDNGHTPVVEEEVGGIPVPAPGLLQVTERRQPLSSVSSLEVHFDLLDLTEL TDMSDQELAE
 VFADSDDENLNTESPAGLHPLPRAGYLRSPSWTRTRAEQSHEKQPLGDPERQATVLDFTL TVERPQED


TRTRPLEQKLISEEDLAANDILDYKDDDDKV

Restriction Sites: SgfI-MluI

Cloning Scheme:



Plasmid Map:



ACCN: NM_001288709

ORF Size: 834 bp

OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
Reconstitution Method:	<ol style="list-style-type: none">1. Centrifuge at 5,000xg for 5min.2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.3. Close the tube and incubate for 10 minutes at room temperature.4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
RefSeq:	NM_001288709.1 , NP_001275638.1
RefSeq Size:	2338 bp
RefSeq ORF:	231 bp
Locus ID:	79007
UniProt ID:	Q9H9R9
Cytogenetics:	16q24.3
MW:	30.8 kDa