

Product datasheet for RC237187

PRMT1 (NM_001207042) Human Tagged ORF Clone

Product data:

Product Type: Expression Plasmids
Product Name: PRMT1 (NM_001207042) Human Tagged ORF Clone
Tag: Myc-DDK
Symbol: PRMT1
Synonyms: ANM1; HCP1; HRMT1L2; IR1B4
Vector: pCMV6-Entry (PS100001)
E. coli Selection: Kanamycin (25 ug/mL)
Cell Selection: Neomycin
ORF Nucleotide Sequence: >RC237187 representing NM_001207042
Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCC**CGGATCGCC**

ATGGCGGCAGCCGAGGCCGGAAGTGCATCATGGAGGTGTCTGTGGCCAGGCGGAAAGCAGTGAGAAGC
CCAACGCTGAGGACATGACATCCAAAGATTACTACTTTGACTCCTACGCACACTTTGGCATCCACGAGGA
GATGCTGAAGGACGAGGTGCGCACCCCTCACTTACCGCAACTCCATGTTTCATAACCGGCACCTCTTCAAG
GACAAGGTGGTGTGCTGGACGTGCGCTCGGCACCGCATCCTCTGCATGTTTGTGCTGCAAGCCGGGGCCC
GCAAGGTATCGGGATCGAGTGTCCAGTATCTCTGATTATGCGGTGAAGATCGTCAAAGCCAACAAGTT
AGACCACGTGGTGACCATCATCAAGGGGAAGGTGGAGGAGGTGGAGCTCCAGTGGAGAAGGTGGACATC
ATCATCAGCGAGTGGATGGGCTACTGCCTTTCTACGAGTCCATGCTCAACACCGTGTCTATGCCCGGG
ACAAGTGGCTGGAGGTGGACATCTATACCGTCAAGGTGGAAGACCTGACCTTACCTCCCGTTCTGCCT
GCAAGTGAAGCGGAATGACTACGTGCACGCCCTGGTGGCTACTTCAACATCGAGTTCACACGCTGCCAC
AAGAGGACCGCTTCTCCACGACCCCGAGTCCCGTACACGCACTGGAAGCAGACGGTGTCTACATGG
AGGACTACCTGACCGTGAAGACGGGCGAGGAGATCTTCGGCACCATCGGCATGCGGCCAACGCCAAGAA
CAACCGGACCTGGACTTACCATCGACCTGGACTTCAAGGGCCAGCTGTGCGAGCTGTCTCTCCACC
GACTACCGGATGCGC

ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT
ACAAGGATGACGACGATAAGGTTTAA



[View online »](#)

OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
Reconstitution Method:	<ol style="list-style-type: none">1. Centrifuge at 5,000xg for 5min.2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.3. Close the tube and incubate for 10 minutes at room temperature.4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
RefSeq:	NM_001207042.3
RefSeq Size:	1192 bp
RefSeq ORF:	858 bp
Locus ID:	3276
UniProt ID:	Q99873
Cytogenetics:	19q13.33
MW:	33.1 kDa
Gene Summary:	This gene encodes a member of the protein arginine N-methyltransferase (PRMT) family. Post-translational modification of target proteins by PRMTs plays an important regulatory role in many biological processes, whereby PRMTs methylate arginine residues by transferring methyl groups from S-adenosyl-L-methionine to terminal guanidino nitrogen atoms. The encoded protein is a type I PRMT and is responsible for the majority of cellular arginine methylation activity. Increased expression of this gene may play a role in many types of cancer. Alternatively spliced transcript variants encoding multiple isoforms have been observed for this gene, and a pseudogene of this gene is located on the long arm of chromosome 5. [provided by RefSeq, Dec 2011]