

## Product datasheet for RC237212

### Chimaerin 2 (CHN2) (NM\_001293073) Human Tagged ORF Clone

#### Product data:

Product Type:	Expression Plasmids
Product Name:	Chimaerin 2 (CHN2) (NM_001293073) Human Tagged ORF Clone
Tag:	Myc-DDK
Symbol:	CHN2
Synonyms:	ARHGAP3; BCH; CHN2-3; RHOGAP3
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)
Cell Selection:	Neomycin
ORF Nucleotide Sequence:	>RC237212 representing NM_001293073 Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC  
GCC**GCGATCGCC**

ATGTTCTCTGAAGAACTGTGGCTGGAAAATGAGAAAAAGTGTGCTGTGGTTCGGAAGTCTAAGCAGGGCA  
GGAAACGCCAAGAAGTGTGGCCGTAGCCTTCGGGGTGAAGGTCCACACGTTCCGAGGCCACACTGGTG  
TGAATATTGTGCCAATTCATGTGGGGCTCATCGCCAAAGGGTCCGGTCTCAGACTGTGGATTGAAC  
GTACACAAACAGTGTCCAAGCAGTTCCTCAATGACTGCCAACCTGATCTCAAGAGGATCAAGAAAGTGT  
ACTGTTGTGACCTCACAACACTTGTGAAGGCTCACAACACTCAGAGACCCATGGTGGTAGACATATGCAT  
TCGGGAAATTGAAGCAAGAGGATTAATTCGGAAGGCCTTACAGAGTCTCTGGGTTCACTGAACACATT  
GAAGATGTCAAAATGGCATTGACAGAGATGGTAAAAGGCCGATATATCTGCCAATGTCTATCCAGACA  
TAAACATCACTGGAGCCCTTAACTGTATTCAGAGACTTACCCATCCCTGTCATCACATATGATAC  
CTATTCAAAATTTATAGATGCAGCAAAAATCTCCAATGCAGATGAGAGGCTGGAAGCCGTCATGAAGTG  
CTGATGCTGCTGCCTCCTGCCACTATGAAACCCTCCGGTACCTAATGATCCACCTCAAAAAGGTTACTA  
TGATGAAAAAGACAATTCATGAATGCAGAAAATCTGGGGATCGTGTTTGGGCCACTCTGATGAGGCC  
CCCTGAGGACAGCACCTGACCACCTGCATGATATGCGGTACCAAAGCTGATTGTGCAGATTTTAATA  
GAAAACGAAGACGTTTTATTC

**ACGCGT**ACGCGGCCGCTCGAGCAGAAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT  
ACAAGGATGACGACGATAAGGTTTAA



[View online »](#)

**Protein Sequence:** >RC237212 representing NM\_001293073  
 Red=Cloning site Green=Tags(s)

MFSEELWLENEKKCAVVRKSKQGRKRQELLAVAFGVKVTFRGPHWCEYCANFMWGLIAQGVRCSDCGLN  
 VHKQCSKHVPNDCQPD LKR I KKVYCCDL TTLVKAHNTQRPMMVVDICIRE IEARGLKSEGL YRVSGFTEHI  
 EDVKMAFDRDGEKADISANVYPDINIITGALKLYFRDLPIP VITYDTYSKF IDAAKISNADERLEAVHEV  
 LMLLPPAHYETLRYLMIHLK KVTMNEKDNFMNAENLGIVFGPTLMRPPEDSTLTLTHDMRYQKLIVQILI  
 ENEDVLF

TRTRPLEQKLISEEDLAANDILDYKDDDDKV


**Restriction Sites:**

Sgfl-MluI

**Cloning Scheme:**



**Plasmid Map:**



**ACCN:** NM\_001293073

**ORF Size:** 861 bp

<b>OTI Disclaimer:</b>	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. <a href="#">More info</a>
<b>OTI Annotation:</b>	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
<b>Components:</b>	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
<b>Reconstitution Method:</b>	<ol style="list-style-type: none"><li>1. Centrifuge at 5,000xg for 5min.</li><li>2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.</li><li>3. Close the tube and incubate for 10 minutes at room temperature.</li><li>4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.</li><li>5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.</li></ol>
<b>RefSeq:</b>	<a href="#">NM_001293073.1</a> , <a href="#">NP_001280002.1</a>
<b>RefSeq Size:</b>	2880 bp
<b>RefSeq ORF:</b>	864 bp
<b>Locus ID:</b>	1124
<b>UniProt ID:</b>	<a href="#">P52757</a>
<b>Cytogenetics:</b>	7p14.3
<b>MW:</b>	33.5 kDa
<b>Gene Summary:</b>	This gene encodes a guanosine triphosphate (GTP)-metabolizing protein that contains a phorbol-ester/diacylglycerol (DAG)-type zinc finger, a Rho-GAP domain, and an SH2 domain. The encoded protein translocates from the cytosol to the Golgi apparatus membrane upon binding by diacylglycerol (DAG). Activity of this protein is important in cell proliferation and migration, and expression changes in this gene have been detected in cancers. A mutation in this gene has also been associated with schizophrenia in men. Alternative transcript splicing and the use of alternative promoters results in multiple transcript variants. [provided by RefSeq, May 2014]