

Product datasheet for **RC237256**

FAM110A (NM_001289147) Human Tagged ORF Clone

Product data:

Product Type: Expression Plasmids
Product Name: FAM110A (NM_001289147) Human Tagged ORF Clone
Tag: Myc-DDK
Symbol: FAM110A
Synonyms: bA371L19.3; C20orf55; F10
Vector: pCMV6-Entry (PS100001)
E. coli Selection: Kanamycin (25 ug/mL)
Cell Selection: Neomycin
ORF Nucleotide Sequence: >RC237256 representing NM_001289147
Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCC**GCGATCGCC**

ATGCCTGTGCACACGCTGAGCCCCGGAGCCCCGTCCGCCCGCCCTACCTTGCCGCTGCGGACCAGGG
TCCCTGGCTACCTGCTACGGGGCCGGCAGATGGTGGAGCCCGAAACCGAGCGCTGTGGAGCGCTGGA
GGCCGACAAGCCAAGTACGTCAAGAGCCTGCAGTGGCCAACACCCGCCAGGAGCTGTGCAGCCCTG
CTGTCCAAACAGCCGCTCTTAGCCCTGAGACTCGCCGCACAGTCTCACGCCAGCCGCCGAGCCCTGC
CTGGCCCCCTGCCGACGGCCCCAGCTGGACCTGGACATCCTCAGCAGCCTCATCGACTTGTGTGACAGCCC
CGTGTCCCCTGCCGAGGCCAGCCGCACTCCTGGACGGGCCGAGGGAGCCGGCCGTCTCCCCAGCCACC
CCTCCGCGACCGCCGCCAGTACCTCTGCGGTCCGCCGGTGGACGTCCGCCCCCTGCCCGCTCGCCTG
CCCGGCCCTGCCATCACCCGGCCCTGCCGCCCTCCAGCCAGCCCGGCCCGGGTTTGAACGTC
CAAGTCGGACTTGAGCGAGCGCTTTCTAGGGCAGCCGTGATCTCGAGCGCTTTTTAACTTCTGCGGC
CTGGACCCGGAGGAGCGAGAGGTTGGTGTGGCCACCTGGCACGGGCCAGCTCGGATATCGTGTCCC
TGGCAGGGCCAGTCTGGGCCGGCAGCTGAAGGGGGTCTCCCGCCGAGCTCGGTGACTGTTGA
GGAGCGGCCCGGAGCGCTTCCCTATGGCGTGTGGTGGAGCGCAATGCCCGCTGATCAAGTGG
TTGATGGGCTAAGGCAGGCTCGGGAGAGCCAGCAGCTGAAGGC

ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT
ACAAGGATGACGACGATAAGGTTTAA



[View online »](#)

Protein Sequence: >RC237256 representing NM_001289147
 Red=Cloning site Green=Tags(s)

MPVHTLSPGAPSAPALPCRLRTRVPGYLLRGPADGGARKPSAVERLEADKAKYVKSLHVANTRQEPVQPL
 LSKQPLFSPETRRTVLTPSRRALPGPCRRPQLDLIDLSSLIDLCDSPVSPAESRTPGRAEGAGRPPAT
 PPRPPPSTSAVRRVDVRPLPASPAPPCSPGPAASSPARPPGLQRKSDLSERF SRAAADLERFFNFCG
 LDPEEARGLGVAHLARASSDIVSLAGPSAGPGSSEGGCSRRSSVTVEERARERVPYGVSVVERNARVIKW
 LYGLRQARESPAPEG


TRTRPLEQKLISEEDLAANDILDYKDDDDKV

Restriction Sites: Sgfl-MluI

Cloning Scheme:



Plasmid Map:



ACCN: NM_001289147

ORF Size: 885 bp

OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
Reconstitution Method:	<ol style="list-style-type: none">1. Centrifuge at 5,000xg for 5min.2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.3. Close the tube and incubate for 10 minutes at room temperature.4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
RefSeq:	NM_001289147.1 , NP_001276076.1
RefSeq Size:	1614 bp
RefSeq ORF:	888 bp
Locus ID:	83541
UniProt ID:	Q9BQ89
Cytogenetics:	20p13
MW:	31.7 kDa