

Product datasheet for **RC237310**

PPP1R7 (NM_001282413) Human Tagged ORF Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	PPP1R7 (NM_001282413) Human Tagged ORF Clone
Tag:	Myc-DDK
Symbol:	PPP1R7
Synonyms:	SDS22
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)
Cell Selection:	Neomycin
ORF Nucleotide Sequence:	>RC237310 representing NM_001282413 Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCC**CGATCGCC**

ATGGAAGAACATGAGCTGCCTGTGGACATGAAACCATCAACCTGGACAGAGATGCAGAGGATGTTGATT
TGAATCACTATCGCATAGGGAAGATTGAAGGATTTGAGGTAAGAAAGTGAAGACTCTCTGCCTCCG
CCAAAATTTAATTAATGCATTGAGAATCTGGAGGAGCTACAGAGTCTTCGAGAGCTGGATCTTTACGAC
AACCAGATCAAGAAGATTGAGAACTGGAGGCGCTAACAGAGCTGGAGATTCTAGATATTTCTTTAATC
TGCTGAGAAACATCGAAGGGTTGACAAGTTGACACGACTGAAAAAACTCTTCTGGTCAACAATAAAAT
CAGTAAAATTGAGAACTTAAGCAACTTACATCAACTACAGATGCTAGAGCTGGGATCTAACCCGATCCGG
GCAATCGAAAATATCGACACCTTAACCAACCTGGAGAGTTTGTGTTTGGGAAAAACAAAATTAATAAC
TTCAGAACCTGGATGCGCTACCAACCTGACAGTCTCAGTATGCAGAGCAACCGGCTGACCAAGATCGA
GGTCTGCAGAACCTGGTGAACCTGCGGGAGCTGTACCTAGCCACAATGGCATCGAGGTCATCGAGGGC
CTGGAGAACAATAACAACTCACGATGTTGGACATTGCATCAAAATAGAATCAAAAAGATTGAAAATATCA
GCCATCTAACAGAGCTGCAAGAGTTCTGGATGAACGACAATCTCCTTGAGAGCTGGAGCGACCTCGACGA
GCTGAAGGGAGCCAGGAGCTGGAGACAGTGTACCTGGAGCGGAACCCCTTGAGAAAGGACCCCCAGTAC
CGCGGAAGGTCATGCTCGCCCTCCCCTCCGTGCGGCAGATCGATGCCACGTTCTGCAGGTTCT

ACGCGTACGCGGCCGCTCGAGCAGAAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT
ACAAGGATGACGACGATAAGGTTTAA



[View online »](#)

Protein Sequence: >RC237310 representing NM_001282413
 Red=Cloning site Green=Tags(s)

MEEHELPMETINLDRDAEDVDLNHYRIGKIEGFEVLKVKYKTLCLRQNLKCIENLEELQSLRELDLYD
 NQIKKIENLEALTELEILDISFNLLRNIEGVDKLTRLKKLFLVNNKISKIENLSNLHQLQMLELGSNRIR
 AIENIDTLTNLESLFLGKNKITKLQNLDALTNLTVLSMQSNRLTKIEGLQNLVNLRELYLSHNGIEVIEG
 LENNNKLTMLDIASNRICKIENISHLTELQEFWMNDNLLESWSDDLDELKGARSLETVYLERNPLQKDPQY
 RRKVMALALPSVRQIDATFVRF

TRTRPLEQKLISEEDLAANDILDYKDDDDKV


Restriction Sites:

Sgfl-MluI

Cloning Scheme:



Plasmid Map:



ACCN: NM_001282413

ORF Size: 903 bp

OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
Reconstitution Method:	<ol style="list-style-type: none">1. Centrifuge at 5,000xg for 5min.2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.3. Close the tube and incubate for 10 minutes at room temperature.4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
RefSeq:	NM_001282413.1 , NP_001269342.1
RefSeq Size:	1878 bp
RefSeq ORF:	906 bp
Locus ID:	5510
UniProt ID:	Q15435
Cytogenetics:	2q37.3
Protein Families:	Druggable Genome, Phosphatase
MW:	35.6 kDa
Gene Summary:	This gene encodes a protein subunit that regulates the activity of the serine/threonine phosphatase, protein phosphatase-1. The encoded protein is required for completion of the mitotic cycle and for targeting protein phosphatase-1 to mitotic kinetochores. Alternate splicing results in multiple transcript variants. [provided by RefSeq, Sep 2013]