

Product datasheet for RC237323

Serum Response Factor (SRF) (NM_001292001) Human Tagged ORF Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	Serum Response Factor (SRF) (NM_001292001) Human Tagged ORF Clone
Tag:	Myc-DDK
Symbol:	SRF
Synonyms:	MCM1
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)
Cell Selection:	Neomycin
ORF Nucleotide Sequence:	>RC237323 representing NM_001292001 Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCC**CGGATCGCC**

ATGATCACCAGTGAGACCGCAAGGCACTGATTGAGACCTGCCTCAACTCGCCAGACTCTCCACCCCGTT
CAGACCCACAACAGACCAGAGAATGAGTGCCACTGGCTTGAAGAGACAGATCTCACCTACCAGGTGTC
GGAGTCTGACAGCAGTGGGGAGACCAAGGACACACTGAAGCCGGGTTTACAGTACCAACTGCCGGT
ACAACCTCCACCATCAAACAGCACCTAGCACCTACCACCATGCAAGTCAGCAGCGGCCCTCCTTTT
CCATCACAACACTACCTGGCACCAGTGTCTGCTAGTGTGACCCCAAGTGTGTGTCAGCAGTGCATG
TGTGCTGAAGAGTACAGGCAGCGGCCCTGTCTCCTCTGGGGGCCTTATGCAGCTGCCTACCAGCTTACC
CTCATGCCTGGTGGGCAGTGGCCAGCAGGTCCCAGTGCAGGCCATTCAAGTGCACCAGGCCACAGC
AAGCGTCTCCCTCCCGTGACAGCAGCACAGACCTCACGCAGACCTCCTCCAGCGGGACAGTGACGCTGCC
CGCCACCATCATGACGTCATCCGTGCCCAACTGTGGGTGGCCACATGATGTACCCTAGCCCGCATGCG
GTGATGTATGCCCCACCTCGGGCCTGGGTGATGGCAGCCTCACCGTGTGAATGCCTTCTCCAGGCAC
CATCCACCATGCAGGTGTCACACAGCCAGTCCAGGAGCCAGGTGGCGTCCCCAGGTGTTCTTGACAGC
ATCATCTGGGACAGTGCAGATCCCTGTTTCAGCAGTTCAGCTCCACCAGATGGCTGTGATAGGGCAGCAG
GCCGGGAGCAGCAGCAACCTCACCGAGCTACAGGTGGTGAACCTGGACACCGCCACAGCACCAAGAGTG
AA

ACGCGTACGCGGCCGCTCGAGCAGAAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT
ACAAGGATGACGACGATAAGGTTTAA



[View online »](#)

Protein Sequence: >RC237323 representing NM_001292001
 Red=Cloning site Green=Tags(s)

MITSETGKALIQTCLNSPDSPPRSDPTTDQRMASATGFEEEDLTYQVSESDSSGETKDTLKPAFTVTNLPG
 TTSTIQAPSTSTTMQVSSGSPFITNYLAPVSASVSPSAVSSANGTVLKSTGSGPVS SGGMLQLPTSFT
 LMPGGAVAQQVPVQAIQVHQAPQQASPSRDSSTDLTQTSSSGTVLTPATIMTSSVPTTVGGHMMYPSHA
 VMYAPTSGLGDGSLTVLNAFSQAPSTMQVSHSQVQEPGGVPQVFLTASSGTVQIPVSAVQLHQMAVIGQQ
 AGSSSNLTELQVNLDTAHSTKSE

TRTRPLEQKLISEEDLAANDILDYKDDDDKV


Restriction Sites:

Sgfl-MluI

Cloning Scheme:



Plasmid Map:



ACCN: NM_001292001

ORF Size: 912 bp

OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
Reconstitution Method:	<ol style="list-style-type: none">1. Centrifuge at 5,000xg for 5min.2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.3. Close the tube and incubate for 10 minutes at room temperature.4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
RefSeq:	NM_001292001.2
RefSeq Size:	3417 bp
RefSeq ORF:	915 bp
Locus ID:	6722
UniProt ID:	P11831
Cytogenetics:	6p21.1
Protein Families:	Druggable Genome, Transcription Factors
Protein Pathways:	MAPK signaling pathway
MW:	31.5 kDa
Gene Summary:	This gene encodes a ubiquitous nuclear protein that stimulates both cell proliferation and differentiation. It is a member of the MADS (MCM1, Agamous, Deficiens, and SRF) box superfamily of transcription factors. This protein binds to the serum response element (SRE) in the promoter region of target genes. This protein regulates the activity of many immediate-early genes, for example c-fos, and thereby participates in cell cycle regulation, apoptosis, cell growth, and cell differentiation. This gene is the downstream target of many pathways; for example, the mitogen-activated protein kinase pathway (MAPK) that acts through the ternary complex factors (TCFs). Two transcript variants encoding different isoforms have been found for this gene. [provided by RefSeq, May 2014]