

## Product datasheet for **RC237337**

### G protein alpha 12 (GNA12) (NM\_001282440) Human Tagged ORF Clone

#### Product data:

|                          |  |
|--------------------------|--|
| Product Type:            | Expression Plasmids  |
| Product Name:            | G protein alpha 12 (GNA12) (NM_001282440) Human Tagged ORF Clone               |
| Tag:                     | Myc-DDK  |
| Symbol:                  | GNA12  |
| Synonyms:                | gep; NNX3; RMP   |
| Vector:                  | pCMV6-Entry (PS100001)   |
| E. coli Selection:       | Kanamycin (25 ug/mL)   |
| Cell Selection:          | Neomycin   |
| ORF Nucleotide Sequence: | >RC237337 representing NM_001282440<br>Red=Cloning site Blue=ORF Green=Tags(s) |

TTTTGTAATACGACTCACTATAGGGCGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC  
GCC**CGATCGCC**

ATGGCTCCACGAGCAGGGCTCGGAGCTGGCCCTCCTCTGCTTCGGAGCCGCGCTGACTGGGTGTGTGGCG  
GATGGGAGGCCCCCTGAGTGTCTCCTCCCTTTTGTACAGCCTGGGTGAAAGGACGGCCTGGTGCA  
CTTCTGCGGGAGGGAGGACATGGGGCTGCTTCCCTCTGGCGCTCCTCCACCCGGCAGCGTGGTTTCT  
TTCCAAATGATTTATCCTCTCATCGTCTGCCGCTCTGTGCTTCATGGCATTGGGGTCTTATCCGGAAGG  
GGATTGTCTTTTTGGGGGAGTCGGTGAAGTACTTCTTGACAACCTGGACCGGATCGGCCAGCTGAA  
TTACTTTCCTAGTAAGCAAGATATCCTGCTGGCTAGGAAAGCCACCAAGGGAATTGTGGAGCATGACTTC  
GTTATTAAGAAGATCCCCTTAAGATGGTGGATGTGGGCGCCAGCGGTCCCAGCGCCAGAAGTGGTTCC  
AGTGCTTCGACGGGATCACGTCCATCCTGTTTATGGTCTCCTCCAGCGAGTACGACCAGGTCTCATGGA  
GGACAGGCGCACCAACCGGCTGGTGGAGTCCATGAACATCTTCGAGACCATCGTCAACAACAAGCTCTTC  
TTCAACGTCTCCATCATTCTCTCCTCAACAAGATGGACCTCCTGGTGGAGAAGGTGAAGACCGTGA  
TCAAGAAGCACTTCCCGACTTCAGGGGCGACCCGACAGGCTGGAGGACGTCCAGCGCTACCTGGTCCA  
GTGCTTCGACAGGAAGAGACGGAACCGCAGCAAGCCACTTCCACCCTTCACCACCGCCATCGACACC  
GAGAACGTCCGCTTCGTGTTCCATGCTGTGAAAGACACCATCCTGCAGGAGAACCTGAAGGACATCATGC  
TGCAG

**ACGCGT**ACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT  
ACAAGGATGACGACGATAAGGTTTAA



[View online »](#)

**Protein Sequence:** >RC237337 representing NM\_001282440  
Red=Cloning site Green=Tags(s)

MAPRAGLGAGPLLLRSRADWVCGWEAALSAPPPFCTAWVERTGLVHFLREGGHGACFPLALLHPGSVVS  
 FQMIYPLIVLPLCASWHWGLIRKGIIVLFWGESVKYFLDNLDRIGQLNYFPSPKQDILLARKATKGI VEHDF  
 VIKKIPFKMVDVGGQRSQRQKWFQCFDGI TSILFMVSSSEYDQVLMEDRRTNRLVESMNIFETIVNNKLF  
 FNVSIILFLNKMDLLVEKVKTVISIKKHFPDFRGGPHRLEDVQRYLVQCFDRKRRNRSKPLFHHFTTAIDT  
 ENVRFVHAVKDTILQENLKDIMLQ


TRTRPLEQKLISEEDLAANDILDYKDDDDKV

**Restriction Sites:** SgfI-MluI

**Cloning Scheme:**



**Plasmid Map:**



**ACCN:** NM\_001282440

**ORF Size:** 915 bp

|                               |   |
|-------------------------------|---|
| <b>OTI Disclaimer:</b>        | The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. <a href="#">More info</a>  |
| <b>OTI Annotation:</b>        | This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.  |
| <b>Components:</b>            | The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).  |
| <b>Reconstitution Method:</b> | <ol style="list-style-type: none"> <li>1. Centrifuge at 5,000xg for 5min.</li> <li>2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.</li> <li>3. Close the tube and incubate for 10 minutes at room temperature.</li> <li>4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.</li> <li>5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.</li> </ol> |
| <b>RefSeq:</b>                | <a href="#">NM_001282440.1</a> , <a href="#">NP_001269369.1</a>   |
| <b>RefSeq Size:</b>           | 4246 bp   |
| <b>RefSeq ORF:</b>            | 918 bp  |
| <b>Locus ID:</b>              | 2768  |
| <b>Cytogenetics:</b>          | 7p22.3-p22.2  |
| <b>Protein Families:</b>      | Druggable Genome  |
| <b>Protein Pathways:</b>      | Long-term depression, MAPK signaling pathway, Regulation of actin cytoskeleton, Vascular smooth muscle contraction  |
| <b>MW:</b>                    | 35.6 kDa  |

**Gene Summary:**

Guanine nucleotide-binding proteins (G proteins) are involved as modulators or transducers in various transmembrane signaling systems (PubMed:22609986, PubMed:15525651, PubMed:15240885, PubMed:17565996, PubMed:12515866, PubMed:16787920, PubMed:16705036, PubMed:23762476, PubMed:27084452). Activates effector molecule RhoA by binding and activating RhoGEFs (ARHGEF12/LARG) (PubMed:15240885, PubMed:12515866, PubMed:16202387). GNA12-dependent Rho signaling subsequently regulates transcription factor AP-1 (activating protein-1) (By similarity). GNA12-dependent Rho signaling also regulates protein phosphatase 2A activation causing dephosphorylation of its target proteins (PubMed:15525651, PubMed:17565996). Promotes tumor cell invasion and metastasis by activating RhoA/ROCK signaling pathway and up-regulating proinflammatory cytokine production (PubMed:23762476, PubMed:16787920, PubMed:16705036, PubMed:27084452). Inhibits CDH1-mediated cell adhesion in process independent from Rho activation (PubMed:11976333, PubMed:16787920). Together with NAPA promotes CDH5 localization to plasma membrane (PubMed:15980433). May play a role in the control of cell migration through the TOR signaling cascade (PubMed:22609986).[UniProtKB/Swiss-Prot Function]