

Product datasheet for **RC237350**

ZSCAN30 (NM_001288711) Human Tagged ORF Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	ZSCAN30 (NM_001288711) Human Tagged ORF Clone
Tag:	Myc-DDK
Symbol:	ZSCAN30
Synonyms:	ZNF-WYM; ZNF397OS; ZNF917
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)
ORF Nucleotide Sequence:	>RC237350 representing NM_001288711 Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
 GCC**CGATCGCC**

ATGGTGGCTGGCAAAGTGTGATGGCAAAGCAAGAAATTGTTGAATGTGTAGCCTCAGCAGCTATGATAT
 CGCCGGGAAAACCTTCTGGAGAAACACATTCCTCAACGGATAGCTGAAGAAGCTTTGGGAGGTCTGGACAA
 CTCTAAGAAGCAAAGGGAAATGCTGCAGGGAACAAAATCAGTCAGCTTCCTTCCAGGACAGACATTTT
 AGTCTGGCAACCTTCAACAGAAGAATCCCAACAGAACACAGTGTCTTGAATCTCATGAGAGTGAAGGAA
 GTTTCAGTATGAATTCAAATGATATTACACAACAGAGCGTTGACACTAGGAAAAGCTCTATGAGTGCTT
 TGAATGTGGGAAGGCCTTTTGCCAGAGCTCAAAGCTGATTAGACATCAGAGAAATTCATACTGGAGAGAGA
 CCTTATGCATGTAAAGAATGTGGCAAAGCCTTCAGTTTGAGCTCAGACCTTGTTAGACATCAGAGAAATC
 ATAGTGGTGAAAAACCTATGAATGTTGTGAATGTGGAAGCCTTCAGGGGCAGTTCAGAGCTCATCAG
 GCATCGGAGAAATCCACTGGAGAGAAGCCTTATGAATGTGGAGAATGTGGGAAGGCCTTCAGCCGGAGC
 TCAGCCCTTATTCAGCATAAGAAAATTCACACTGGAGATAAAAGCTATGAATGTATTGCATGTGGTAAGG
 CTTTTGGTAGGAGTCTATCCTTATTGAACATCAAAGAATTCACACTGGAGAGAAACCTTATGAATGTAA
 TGAATGTGGAAAATCCTTCAATCAGAGCTCAGCCCTCACCCAGCACCAGAGAATTCACACTGGAGAGAAG
 CCTTATGAATGTAGTGAATGTAGAAAAACATTTAGGCATAGGTCAGGCCTTATGCAGCATCAGAGAACAC
 ACACCAGAGTT

ACGCGTACGCGGCCGCTCGAGCAGAAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT
 ACAAGGATGACGACGATAAGGTTTAA


[View online »](#)


Protein Sequence: >RC237350 representing NM_001288711
Red=Cloning site Green=Tags(s)

MVAGKVLMAKQEIVECVASAAMISPGKLPGETHSQRIAEALGGLDNSKKQKGNAAGNKISQLPSQDRHF
SLATFNRRIPTEHSVLESHESEGSFSMNSNDITQQSVDTREKLYECFDCGKAFQSSKLIRHQRIHTGER
PYACEKCGKAFSLSSDLVRHQRIHSGEKPYECCECGKAFRGSSSELIRHRIHTGEKPYECGECGKAFSRS
SALI QHKKIHTGDKSYECIACGKAFGRSSILIEHQRIHTGEKPYECNECGKSFNQSSALTQHRIHTGEK
PYECSECRKTRHRSGLMQHQRTHTRV

TRTRPLEQKLISEEDLAANDILDYKDDDDKV

Restriction Sites: SgfI-MluI

Cloning Scheme:



ACCN: NM_001288711

ORF Size: 921 bp

OTI Disclaimer: The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

OTI Annotation: This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

Reconstitution Method:

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

Note: Plasmids are not sterile. For experiments where strict sterility is required, filtration with 0.22um filter is required.

RefSeq: NM_001288711.1, NP_001275640.1

RefSeq Size: 3439 bp

RefSeq ORF: 924 bp

Locus ID: 100101467


UniProt ID: Q86W11

Cytogenetics: 18q12.2

MW: 34.9 kDa

Gene Summary: May be involved in transcriptional regulation.[UniProtKB/Swiss-Prot Function]

Product images:



Circular map for RC237350