

Product datasheet for **RC237404**

TEX9 (NM_001286449) Human Tagged ORF Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	TEX9 (NM_001286449) Human Tagged ORF Clone
Tag:	Myc-DDK
Symbol:	TEX9
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)
ORF Nucleotide Sequence:	>RC237404 representing NM_001286449 Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
 GCC**CGATCGCC**

ATGAAATCATGTGATGACGAAGATGATTACAGTTTAAGAGGTCTGTTACCATCTGAAGGGATAGTTCATC
 TTCATTGAGAACTAAGCCAAAGACCAAAACATTGATCCTGTAAACAAGGTTCAAAACAAATTACACTC
 TGCAATTAAGGAAGGAAACAAATCAAGTGTTAAATTGAAATACTCTGATGTCCAACTGCCGACGAT
 GTTGCCATTCCAGAGGATTTCTCAGACTTTCCCTTGCAAAAACAATTAGCAAAATTGAAGGGCAACTGG
 AGGAAGAAGGCTTACCTGAATATATAGATGATTTTTTCTGGTGTTAGTAATGACATTGGAACAGAAGC
 ACAGATCAGATTTCTAAAGGCCAACTCCATGTTATGCAGGAGGAATTGGATAATGTTGTATGTGAATGC
 AATAAAAAGGAGGATGAAATTCAGAATTTAAAGTCTCAAGTAAAAATTTGAAGAAGATTTTATGAGAC
 AGCAGCGAACAATTAATATGCAACAGTCTCAAGTAGAAAAATACAAACTCTTTTCGAAGAAGCAAAACA
 AAAGTATGATGGATTACAGCAACAGTTGTCTTCAGTAGAAAGGAATTAGAAAAATAAAGAAGACTGCAA
 AAACAGGCTGCAAGTAGTCAAAGTGCCACAGAGGTTGCTTGAATAGAGCTCTAGAAGAAGCAGAAAAGT
 ATAACTGGAGTTAAGTAAATTAAGGCAAAATAACAAGGACATAGCAATGAAGAACACAAAAAATTGA
 AGTGTTAAATCAGAAAACAAGAAGCTAGAAAAACAAAAGGAGAATTAATGATAGGGTTCAAGAAACAG
 TTAATTAATTGATGTTTTAAAGGCAAAAGATGCATATTGAAGCTGCCAAGATGCTATCTTTCACTG
 AGGAGGAATTTATGAAAGCACTGAATGGGAAATTC

ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAATGATATCCTGGATT
 ACAAGGATGACGACGATAAGGTTTAA


[View online »](#)


Protein Sequence: >RC237404 representing NM_001286449
 Red=Cloning site Green=Tags(s)

MKSCDDEDDYSLRGLLPSEGIVHLHSETKPKTKNIDPVNKVQNKLSANKGRKTNSSVKLYSDVQTADD
 VAIPEDFSDFSLAKTISKIEGQLEEEGLPEYIDDI FSGVNDIGTEAQIRFLKAKLHVMQEELDNVCEC
 NKKEDIQNLKSQVKNFEEDFMRQRTINMQSQVEKYKTLFEEANKKYDGLQQQLSSVERELNKRRLO
 KQAASSQSATEVRLNRALEEAKEYKLELSKLRQNNKDIAEEHKKIEVLKSENKKLEKQKGELMIGFKKQ
 LKLIDVLKRQKMHEAAKMLSFTEEFMKALEWGNS

TRTRPLEQKLISEEDLAANDILDYKDDDDKV

Restriction Sites: SgfI-MluI

Cloning Scheme:



ACCN: NM_001286449

ORF Size: 948 bp

OTI Disclaimer: The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

OTI Annotation: This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

Reconstitution Method:

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

Note: Plasmids are not sterile. For experiments where strict sterility is required, filtration with 0.22um filter is required.

RefSeq: NM_001286449.1, NP_001273378.1

RefSeq Size: 1782 bp

RefSeq ORF: 951 bp


Locus ID: 374618

UniProt ID: Q8N6V9

Cytogenetics: 15q21.3

MW: 37 kDa

Product images:



Circular map for RC237404