

Protein Sequence: >RC237408 protein sequence
Red=Cloning site Green=Tags(s)

MKVLWAALLVTFLAGCQAKVEQAVETEPEPELRQQTEWQSGQRWELALGRFWDYLRWVQTLSEQVQEELL
 SSQVTQELRALMDETMKELKAYKSELEEQLTPVAEETRARLSKELQAAQARLGADMEDVCGRLVQYRGEV
 QAMLGQSTEELRVRLASHLRKLRKLLRDADDLQKRLAVYQAGAREGAERGLSAIRERLGPLVEQGRVRA
 ATVGLAGQPLQERAQAWGERLRARMEEMGSRTDRDLDEVKEQVAEVRAKLEEQAQQIRLQAEAFQARLK
 SWFEPLVEDMQRWAGLVEKVQAAVGTSAAPVPSDNH

TRTRPLEQKLISEEDLAANDILDYKDDDDDKV

Chromatograms: https://cdn.origene.com/chromatograms/mk6267_a06.zip

Restriction Sites: SgfI-MluI

Cloning Scheme:

Cloning sites used for ORF Shuttling:



* The last codon before the Stop codon of the ORF

ACCN: NM_001302689

ORF Size: 951 bp

OTI Disclaimer: The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

OTI Annotation: This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

Reconstitution Method:

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

RefSeq: [NM_001302689.2](#)

RefSeq Size: 1165 bp

RefSeq ORF: 954 bp

Locus ID: 348

UniProt ID: [P02649](#)

Cytogenetics: 19q13.32


Protein Families: Adult stem cells, Druggable Genome, Secreted Protein, Stem cell - Pluripotency

Protein Pathways: Alzheimer's disease

MW: 36.2 kDa

Gene Summary: The protein encoded by this gene is a major apoprotein of the chylomicron. It binds to a specific liver and peripheral cell receptor, and is essential for the normal catabolism of triglyceride-rich lipoprotein constituents. This gene maps to chromosome 19 in a cluster with the related apolipoprotein C1 and C2 genes. Mutations in this gene result in familial dysbetalipoproteinemia, or type III hyperlipoproteinemia (HLP III), in which increased plasma cholesterol and triglycerides are the consequence of impaired clearance of chylomicron and VLDL remnants. [provided by RefSeq, Jun 2016]

Product images:



Circular map for RC237408