

Product datasheet for **RC237421**

zinc finger protein 138 (ZNF138) (NM_001271639) Human Tagged ORF Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	zinc finger protein 138 (ZNF138) (NM_001271639) Human Tagged ORF Clone
Tag:	Myc-DDK
Symbol:	ZNF138
Synonyms:	pHZ-32
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)
Cell Selection:	Neomycin
ORF Nucleotide Sequence:	>RC237421 representing NM_001271639 Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCC**CGGATCGCC**

ATGGGACCACTGACATTTATGGATGTGGCCATAGAGTTCTCTTTGGAGGAGTGCCAGTGCCTGGACACTG
CACAGCGGAATGTATATAGGCATGTGATGTTAGAGAACTACAGAAACCTGGTTTTCTGGACCTGATTAC
CTGTCTGGAACAAGGAAAAGAGCCCTGGAATATGAAGAGACATGAGATGGTGGTAGCCAAACATTCAGCT
CTGTGTTCTCGTTTTGCCAAGACCTTTGGCTAGAGCAGAACATAAAAAGATTCTTTCCAAAAAGTGACAC
TGAGCAGATATGGAAAATATGGACATAAGAATTTACAGTTAAGAAAAGGCTGTAAAAGTGTGGATGAGTG
TAAGGGACACCAAGGAGTTATAATGGACTTAACCAATGTTTGAAAATTACCACAAGCAAAAATTTCAA
TGTAATAAATATGTAAAAGTCATGCATAAATTTCAAATTCAAAATAGACACAAGATAAGACATACTGAAA
ATAAACATTTTCAGATGTAAAGAAATGTGACAAATCACTTTGCATGCTTTCACGCCTAACTCAACATAAAAA
AATTCATACTAGAGAGAATTTCTACAAATGTGAAGAGTGTGGAAAAACCTTTAACTGGTCCACAAACCTT
TCTAAACCTAAGAAAATTCATACTGGAGAAAAACCTACAAATGTGAAGTATGTGGAAAAGCCTTTCACC
AATCCTCAATCCTTACTAAACATAAGATAATTCGTAAGGAGAAAAACCTATAAATGTGCACACTGTGG
CAAAGCCTTTAAACAGTCTCACACCTTACTAGACATAAGATAATTCATACTGAAGAGAAAACCTACAAA
TGTGAACAATGTGGCAAGGTCTTTAAGCAGTCCCAACCTTACTAAACATCAGATAATTTATACTGGAG
AGGAACCATACAAATGTGAGGAATGTGGCAAAGCTTTTAACTATCT

ACGCGTACGCGGCCGCTCGAGCAGAAACTCATCTCAGAAGAGGATCTGGCAGCAATGATATCCTGGATT
ACAAGGATGACGACGATAAGGTTTAA



[View online »](#)

OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
Reconstitution Method:	<ol style="list-style-type: none">1. Centrifuge at 5,000xg for 5min.2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.3. Close the tube and incubate for 10 minutes at room temperature.4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
RefSeq:	NM_001271639.1 , NP_001258568.1
RefSeq Size:	2620 bp
RefSeq ORF:	960 bp
Locus ID:	7697
Cytogenetics:	7q11.21
Protein Families:	Transcription Factors
MW:	37.9 kDa
Gene Summary:	May be involved in transcriptional regulation as a repressor.[UniProtKB/Swiss-Prot Function]