

Product datasheet for RC237437

MYADM (NM_001290192) Human Tagged ORF Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	MYADM (NM_001290192) Human Tagged ORF Clone
Tag:	Myc-DDK
Symbol:	MYADM
Synonyms:	SB135
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)
ORF Nucleotide Sequence:	>RC237437 ORF sequence Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCC**CGGATCGCC**

RCATGCCAGTGACGGTAACCCGCACCACCATCACAACCACCACGACGTCATCTTCGGGCCTGGGGTCCCC
CATGATCGTGGGGTCCCCTCGGGCCCTGACACAGCCCCTGGGTCTCCTTCGCCTGCTGCAGCTGGTGTCT
ACCTGCGTGGCCTTCTCGCTGGTGGCTAGCGTGGGCGCCTGGACGGGTCCATGGCAACTGGTCCATGT
TCACCTGGTGCTTCTGCTTCTCCGTGACCCTGATCATCCTCATCGTGGAGCTGTGCGGGCTCCAGGCCCC
CTTCCCCCTGTCTTGGCGCAACTCCCCATCACCTTCGCCTGCTATGCGGCCCTTCTGCTCTCGGCC
TCCATCATCTACCCACCACCTATGTCCAGTTCCTGTCCACGGCGTTTCGCGGGACCACGCCATCGCCG
CCACCTTCTTCTCCTGCATCGCGTGTGTGGCTTACGCCACCAGTGGCCTGGACCCGGGCCCGGCCGG
CGAGATCACTGGCTATATGGCCACCGTACCCGGGCTGCTGAAGGTGCTGGAGACCTTCGTTGCCTGCATC
ATCTTCGCGTTCATCAGCGACCCCAACCTGTACCAGCACCAGCCTGCCCTGGAGTGGTGCCTGGCGGTGT
ACGCCATCTGTTTCATCCTAGCGGCCATCGCCATCCTGCTGAACCTGGGGAGTGACCAACGTGCTACC
CATCCCCCTCCCCAGCTTCTGTGCGGGCTGGCCTTGTGTCTGTCTCCTCTATGCCACCGCCCTGTT
CTCTGGCCCTCTACCAGTTCGATGAGAAGTATGGCGGCCAGCCTCGGCGCTCGAGAGATGTAAGCTGCA
GCCGCAGCCATGCCTACTACGTGTGTGCCTGGGACCGCCGACTGGCTGTGGCCATCCTGACGGCCATCAA
CCTACTGGCGTATGTGGCTGACCTGGTGCACCTGCCCCACCTGGTTTTTGTCAAGGTC

ACGCGTACGCGGCCGCTCGAGCAGAAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT
ACAAGGATGACGACGATAAGGTTTAA



[View online »](#)

Reconstitution Method:

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

RefSeq: [NM_001290192.1](#), [NP_001277121.1](#)

RefSeq Size: 3100 bp

RefSeq ORF: 969 bp

Locus ID: 91663


UniProt ID: [Q96S97](#)

Cytogenetics: 19q13.42

Protein Families: Transmembrane

MW: 35.3 kDa

Product images:



Circular map for RC237437