

Product datasheet for **RC237540**

CCDC198 (NM_001283056) Human Tagged ORF Clone

Product data:

Product Type: Expression Plasmids
Product Name: CCDC198 (NM_001283056) Human Tagged ORF Clone
Tag: Myc-DDK
Symbol: CCDC198
Synonyms: C14orf105
Vector: pCMV6-Entry (PS100001)
E. coli Selection: Kanamycin (25 ug/mL)
Cell Selection: Neomycin
ORF Nucleotide Sequence: >RC237540 representing NM_001283056
 Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
 GCC**CGGATCGCC**

ATGGGCCTGAGTCACTCTAAGACTCACCTTAGGGTGATCAAAGTAGCACCTTTGCAAAACAAGAGGTAG
 AGACTCCCTCGGCTGGCCGTGTGGACTTTGCATTCAATCAGAATTTGGAAGAAAAGACTTCATATCACT
 GGCAAGACTGCAGGACCAGAATAAAGCCTTGAAGGGCAGCTGCCACCTTTACAAGAAAAGTGGTATGGA
 AGATATTTCTACAGCATCCAGAGACATGTATTTGACATCCCACTGGAACACAGAGAAACAAGTATTATTA
 AAAGGCATCCACCCCAAAGACTTCAAAGCTTGAACCCATTGACTTGCCACGAGTAATTACTTCAGGAAG
 ACTCCTGAGCCAGCGAGAAGCCAGGACAATGCACAAAGCAAAGGTACTAGAAAAGAAAATGCAAACCTCCA
 ATGTATACTTCTGAGAACAGACAATATTTGCATAAGATGCAAGTGCTGGAAATGATCCGTAAGACAAG
 AGGCCAAATGGAGTTAAAGAAAAGTCTTCATGGAGAGGCAAGAATTAATAAGCAAAGTCCAAGGGACCA
 TAAAGCCAAGAAAACCTTCAAAGCACCCCAAGGAATGATGACCATGACCTTCTAACCATGTTGCCTGAT
 GAAATCTTGAACAGAGGTCCCGACAGAGTTTCGCTCCTGTTACCCAGGCTGGAGTGCAATGGCATGATCT
 CGGCTCACGCAAACTTCGCCTCCAAGTCAAATGATTCTCCTGCCTCAGCCTCCGAGTAGCTGGGAT
 TACAGGAAATTCAAAGAATACAGAATTTTTGAAACATCAAGCAGTGAATAACTACTGTCCCTGGAAAATT
 GGCAAAATGGAACATGGCTTCATGAACAAGAGGCCAGGGACAGCTTCTCTGGGACAGTTCAGCTCTG
 ACTCAGATGAGCAGGGGAAAGATGAGAAGAAGCCACGAGCACTGGTGAGGACCAGGACAGAGAGAATCCC
 ACTTTTCGATGAGTTTTTTGATCAAGAA

ACGGTACGGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAATGATATCCTGGATT
 ACAAGGATGACGACGATAAGGTTTAA



[View online »](#)

Protein Sequence: >RC237540 representing NM_001283056
Red=Cloning site Green=Tags(s)


MGLSHSKTHLRVIKVAPLQNKEVETPSAGRVDFAFNQNLEEKTSYSLARLQDQNKALEGQLPPLQENWYG
 RYSTASRDMYFDIPLHRETSIIKRHPQRLQKLEPIDLPRVITSGRLLSQREARTMHKAKVLEKKMQTP
 MYTSENRYLHKMQVLEMIRKRQEAQMELKKS LHGEARINKQSPRDHKAKKTLQSTPRNDDHLLTMLPD
 EILNRGPDRVSLLLPRLECNMISAHRKLRPLSSNDSPASASRVAGITGNSKNTEFLKHQAVNNYCPWKI
 GKMETWLHEQEAQGQLLWSSSSDSDEQKDEKKPRALVRTRTERIPLFDEFDQE

TRTRPLEQKLISEEDLAANDILDYKDDDDKV


Restriction Sites:

Sgfl-MluI

Cloning Scheme:



Plasmid Map:



ACCN: NM_001283056

ORF Size: 1008 bp

OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
Reconstitution Method:	<ol style="list-style-type: none">1. Centrifuge at 5,000xg for 5min.2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.3. Close the tube and incubate for 10 minutes at room temperature.4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
RefSeq:	NM_001283056.2
RefSeq Size:	3214 bp
RefSeq ORF:	1011 bp
Locus ID:	55195
UniProt ID:	Q9NVL8
Cytogenetics:	14q22.3
MW:	39.4 kDa