


Protein Sequence: >RC237583 representing NM_001286536
 Red=Cloning site Green=Tags(s)

MAAMLKCVMSGSQVKVFGKAVQALSRISEFWLDPSKKGLALRCVNSSRSAYGCVLFSPVFFQHYQWSAL
 VKMSENELDTTLHLKCKLGMKSILPIFRCLNSLERNIEKCRIFTRSDKCKVVIQFFYRHDLNAVHSEMF
 VGSDEFDFQIGMDTEITFCFKELKGILTFSEATHAPISIIYDFPGKPLALSIDDMLEAVNFILATLADE
 QSRASSPQSLCLSQKRKRSDLIEKKAGKNVTGQALECISKKAAPRRLYPKETLTNISALENCGSPAMKRV
 DGDVSEVSESSVSNTTEVPGSLCLRKCAAGKNLMEVMPNISVLSDIEQYPPLPRQPEFGVE

TRTRPLEQKLISEEDLAANDILDYKDDDDKV

Restriction Sites: SgfI-MluI

Cloning Scheme:



ACCN: NM_001286536

ORF Size: 1023 bp

OTI Disclaimer: The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

OTI Annotation: This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

Reconstitution Method:

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

Note: Plasmids are not sterile. For experiments where strict sterility is required, filtration with 0.22um filter is required.

RefSeq: NM_001286536.2

RefSeq Size: 1582 bp

RefSeq ORF: 1026 bp


Locus ID: 144715

UniProt ID: Q6WBX8

Cytogenetics: 12q24.11

MW: 38.5 kDa

Product images:



Circular map for RC237583