

Product datasheet for **RC237585**

RAD52 (NM_001297421) Human Tagged ORF Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	RAD52 (NM_001297421) Human Tagged ORF Clone
Tag:	Myc-DDK
Symbol:	RAD52
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)
Cell Selection:	Neomycin
ORF Nucleotide Sequence:	>RC237585 representing NM_001297421 Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCC**GCGATCGCC**

ATGTTTGGTTACAATGGCTGGGCACACTCCATCACGCAGCAGAATGTGGATTTTGTGGACCTCAACAATG
GCAAGTTCTACGTGGGAGTCTGTGCATTTGTGAGGTCACAGTGAAGGATGGTTCATATCATGAAGATGT
TGTTATGGTGTAGTGAGGGCCTCAAGTCCAAGGCTTTATCTTTGGAGAAGCAAGGAAGGAGGGCGGTG
ACAGACGGGCTGAAGCGAGCCCTCAGGAGTTTTGGGAATGCACTTGAAAAGTGTATTCTGGACAAAGACT
ACCTGAGATCACTAAATAAGCTTCCACGCCAGTTGCCTCTTGAAGTGGATTTAACTAAAGCGAAGAGACA
AGATCTTGAACCGTCTGTGGAGGAGCAAGATACAACAGCTGCCGACCGAACATGGCCCTGGGACACCCA
CAGCTGCAGCAGGTGACCTCCCCTTCCAGACCCAGCCATGCTGTGATACCGGGGACCAGGACTGCAGCT
CCCGAAGCCTGAGCTCATCCGCCGTGGAGAGCGAGGCCACGCACCAGCGGAAGCTCCGGCAGAAGCAGCT
GCAGCAGCAGTCCGGGAGCGGATGGAGAAGCAGCAGGTTGAGTCTCCACGCCGTGAGTGTGAGAAGAGT
GAGGCAGCGCCTCCGGCCCTCCTGTGACGCACAGCACTCCTGTAAGTGTCTCAGAACCCTCCTGGAGA
AAGACTTCTTGCAGGAGTGACTCAAGAATTAATCAAGACTCTTGAAGACAACCTTGAAAAGTGGGCTGT
GACTCCCGATGCAGGGGATGGTGTGGTCAAGCCCTCGTCTAGAGCAGACCCAGCCAGACCTCTGACACA
TTAGCCTTGAACAACAGATGGTGAACCAAGGACTCCACACAGCGTTTGCCACCAGAAACCAAG
CAAAATCTGGATCTTGGGACCTCCAAACTTATAGCGCTGACCAACGCACAACAGGAAAAGTGGGAATCTCA
TAGGAAGACCCAGGACATGAAGAAAAGGAAATATGATCCATCT

ACGCGTACGCGGCCGCTCGAGCAGAAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT
ACAAGGATGACGACGATAAGGTTTAA



[View online »](#)

Protein Sequence: >RC237585 representing NM_001297421
 Red=Cloning site Green=Tags(s)

MFGYNGWAHSITQQNVDFVDLNNKGFYVGVC AFVRVQLKDGSYHEDVGYGVSEGLKSKALSLEKARKEAV
 TDGLKRALRSFGNALGNCILDKDYLRSLNKLPRQLPLEVDLTKAKRQDLEPSVEEARYNSCRPNMALGHP
 QLQQVTSRPSHAVIPADQDCSSRSLSSSAVESEATHQRKLRQKQLQQQFRERMEKQQVVRVSTPSAEKS
 EAAPPAPVTHSTPVTVSEPLLEKDFLAGVTQELIKTLEDNSEKWAVTPDAGDGVVVKPSSRADPAQTSDT
 LALNNQMVTVNRPVSVCHQKPKQAKSGSWDLQTY SADQRTTGNWESHRKSDMKKRKYDPS


TRTRPLEQKLISEEDLAANDILDYKDDDDKV

Restriction Sites: SgfI-MluI

Cloning Scheme:



Plasmid Map:



ACCN: NM_001297421

ORF Size:	1023 bp
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
Reconstitution Method:	<ol style="list-style-type: none">1. Centrifuge at 5,000xg for 5min.2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.3. Close the tube and incubate for 10 minutes at room temperature.4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
RefSeq:	NM_001297421.1 , NP_001284350.1
RefSeq Size:	2847 bp
RefSeq ORF:	1026 bp
Locus ID:	5893
Cytogenetics:	12p13.33
Protein Families:	Druggable Genome
Protein Pathways:	Homologous recombination
MW:	38.4 kDa
Gene Summary:	The protein encoded by this gene shares similarity with <i>Saccharomyces cerevisiae</i> Rad52, a protein important for DNA double-strand break repair and homologous recombination. This gene product was shown to bind single-stranded DNA ends, and mediate the DNA-DNA interaction necessary for the annealing of complementary DNA strands. It was also found to interact with DNA recombination protein RAD51, which suggested its role in RAD51 related DNA recombination and repair. A pseudogene of this gene is present on chromosome 2. Alternative splicing results in multiple transcript variants. Additional alternatively spliced transcript variants of this gene have been described, but their full-length nature is not known. [provided by RefSeq, Jul 2014]