

Product datasheet for **RC237661**

C7orf26 (NM_001303039) Human Tagged ORF Clone

Product data:

Product Type: Expression Plasmids
Product Name: C7orf26 (NM_001303039) Human Tagged ORF Clone
Tag: Myc-DDK
Symbol: C7orf26
Vector: pCMV6-Entry (PS100001)
E. coli Selection: Kanamycin (25 ug/mL)
Cell Selection: Neomycin
ORF Nucleotide Sequence: >RC237661 representing NM_001303039
Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCC**GCGATCGCC**

ATGAGCGACATCCGCCACTCGCTGCTGCGCCGCGATGCGCTGAGCGCCGCCAAGGAGGTGTGTACCACC
TGGACATCTACTCAGCAGCCAGCTGCAGAGCGCGCCGCTGCCATCGTGGACAAGGGCCCCGTGGAGCT
GCTGGAGGAGTTCGTGTTCCAGGTGCCAAGGAGCGCAGCGCGCAGCCCAAGGAGCAAACCAAGGACTCT
GTTCCGGCAGATTATTTTTTCATCCCTTTTCAGCCCTCAAGGGAACAAAGCCGATGACAGCCGGATGAGCT
TGTTGGGAAAACCTGGTCTCCATGGCGGTGGCTGTGTGTCGAATCCCGGTGTTGGAGTGTGCTGCCTCCTG
GCTTCAGCGGACGCCCGTGGTTACTGTGTGAGGTTAGCCAAGGCCCTTGTAGATGACTACTGCTGTTTTG
GTGCCGGGATCCATTCAGACGCTGAAGCAGATATTCAGTGCCAGCCGAGATTCTGCTGCCAGTTTCATCA
CCTCCGTTACCGCGCTCTATGACCTGTCATCAGATGACCTCATTCCACCTATGGACTTGCTTGAATGAT
TGTCACCTGGATTTTTGAGGACCAAGCGTCTGCAAGTGCTCATGACGCTGCAGCTGCACCTGACCGAG
AAGAATCTGTATGGCGCCTGGGGCTGATCCTCTCGACCACATGGTCCCCTGGTAGAGGAGATCAACA
GGTTGGCGGATGAACGAACCCCTCAACGCCCTCCAGGAGATTGAGCTCTCGCTGGACCGGCTGGCGCA
GGCTCTGCAGGTGGCCATGGCCTCAGGAGCTCTGCTGTGCACGAGAGATGACCTGAGAACCTTGTGCTCC
AGGCTGCCCATATAACCTCCTCCAGCTGGTGATCTCGGGTCCCCTGTCAGCAGTGCCTCACGCCGCGC
TCCCCCGGGGTTCTACCCACATCCACACGCCCGCTGGGCTACGGGGCTGTCCCGGCCACCCCGC
CGCCACCCCGCCTGCCACGCATCCCGGCCACACCTTCATCTCCGGCGTGACCTTTCCCTCAGGCC
ATCCGC

ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCTGGATT
ACAAGGATGACGACGATAAGGTTTAA



[View online »](#)


Protein Sequence: >RC237661 representing NM_001303039
 Red=Cloning site Green=Tags(s)

MSDIRHSLLRDALSAAKEVLYHLDIYFSSQLQSAPLPIVDKGPVELLEEFVFPVKERSAQPKEQTKDS
 VRQIIFSSLFSPQGNKADDSRMSLLGKLVSMVAVCRIPLVECAASWLQRTPVVYCVRLAKALVDDYCCCL
 VPGSIQTLKQIFASPRFCCQFITSVTALYDLSSDDLIPPMDLLEMIVTWIFEDPSVLQVLMTLQLHLTE
 KNLYGRLGLILFDHMVPLVEEINRLADELNPLNASQEIELSLDRLAQAALQVAMASGALLCTRDDLRTLCS
 RLPHNLLQLVISGPVQQSPHAALPPGFYPHIHTPPLGYGAVPAHPAHPALPTHGHTFISGVTFFFRP
 IR


TRTRPLEQKLISEEDLAANDILDYKDDDDKV

Restriction Sites: SgfI-MluI

Cloning Scheme:



Plasmid Map:



ACCN: NM_001303039

ORF Size:	1056 bp
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
Reconstitution Method:	<ol style="list-style-type: none">1. Centrifuge at 5,000xg for 5min.2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.3. Close the tube and incubate for 10 minutes at room temperature.4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
RefSeq:	NM_001303039.2
RefSeq Size:	1904 bp
RefSeq ORF:	1059 bp
Locus ID:	79034
UniProt ID:	Q96N11
Cytogenetics:	7p22.1
MW:	39.5 kDa