

## Product datasheet for **RC237810**

### ARFIP1 (NM\_001287431) Human Tagged ORF Clone

#### Product data:

**Product Type:** Expression Plasmids  
**Product Name:** ARFIP1 (NM\_001287431) Human Tagged ORF Clone  
**Tag:** Myc-DDK  
**Symbol:** ARFIP1  
**Synonyms:** HSU52521  
**Vector:** pCMV6-Entry (PS100001)  
**E. coli Selection:** Kanamycin (25 ug/mL)  
**Cell Selection:** Neomycin  
**ORF Nucleotide Sequence:** >RC237810 representing NM\_001287431  
Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC  
GCC**GCGATCGCC**

ATGGCTCAAGAATCTCCAAAAATTCAGCAGCAGAAATCCAGTGACTAGTAATGGAGAAGTTGATGACT  
CTCGTGAACATAGCTTTAATAGGGATTTGAAGCATTCAATACCATCTGGACTTGGTCTCTCAGAAACCCA  
AATTACATCTCATGGCTTTGACAATACCAAAGAGGGTGTATTGAAGCAGGAGCATTTCAGGTTCCCA  
GCACCGCCACTGCCATCTGTTATGTCTCTAGCAGGGTTGCAGCTAGTCGACTGGCTCAGCAAGGAAGTG  
ATTTAATTGTTCTGCAGGTGGCCAGAGAACACAGACAAAAAGTGGACCAGTTATTCTAGCAGATGAAAT  
TAAAAATCCTGCAATGGAAAAGTTAGAAGTTGTTAGAAAATGGAGTCTAAACACCTATAAGTGTACTCGA  
CAGATTATCTCTGAGAAGCTAGGCCGTGGCTCAAGAACTGTGGACCTTGAAGCTCAGATTGATA  
TATTAAGGGATAACAAGAAAAATATGAAAATTTTTAAAAGTGGCTCAAACATTGTGACCCAGCTTTT  
CCAGATGGTACATACCAAAGGCAACTTGGAGATGCATTTGCTGACCTGAGTTTGAAGTCACTAGAAGT  
CATGAAGAATTTGGCTATAATGCCGATACCCAGAACTGCTGGCTAAAAATGGAGAGACTCTTCTGGGG  
CCATTAATTTTTTCATTGCTAGTGTGAACACTTTGGTGAATAAAACCATTGAAGATACATTAATGACTGT  
GAAACAGTATGAAAGTGCCAGGATTGAATATGATGCATATCGCACTGATTTGGAAGAAGTGAATCTTGG  
CCACGTGACGCAAACTCTGCCAAAGATTGAGCAGTCACAGCATCTCTCCAAGCACATAAGGAAAAAT  
ATGATAAAATGCGCAATGATGTTTCTGTCAAATTTGAAATTTCTAGAAGAAAAAAGGTTAAAGTATTGCA  
CAATCAGCTGGTCTTTTCCACAATGCCATTGCCGCTTACTTTGCTGGGAATCAGAAGCAGCTTGAACAG  
ACACTTAAACAGTTCCATATCAAATTGAAAACCCCTGGAGTGGATGCCCATCTTGGCTTGAAGAACAG

**ACGCGT**ACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCTGGATT  
ACAAGGATGACGACGATAAGGTTTAA



[View online »](#)

**Protein Sequence:** >RC237810 representing NM\_001287431  
 Red=Cloning site Green=Tags(s)

MAQESPKNSAAEIPVTSNGEVDDSDREHSFNRLKHSPLPSGLGLSETQITSHGFDNTKEGVIEAGAFQGGSP  
 APPLPSVMSPSRVAASRLAQQGSDLIVPAGGQRTQTKSGPVILADEIKNPAMEKLELVRKWSLNTYKCTR  
 QIISEKLRGRSRTVDLELEAQIDILRDNKKKYENILKLAQTLSTQLFQMVHTQRQLGDAFADLSLKSLEL  
 HEEFGYNADTQKLLAKNGETLLGAINFFIASVNTLVNKTIEDTLMTVKQYESARIEYDAYRTDLEELNLG  
 PRDANTLPKIEQSQHLFQAHKKEYDKMRNDVSVKLFLEENKVKVLHNQLVLFHNAIAAYFAGNQKQLEQ  
 TLKQFHlKlKTPGVDAPSWLEEQ

TRTRPLEQKLISEEDLAANDILDYKDDDDKV


**Restriction Sites:**

SgfI-MluI

**Cloning Scheme:**



**Plasmid Map:**



**ACCN:** NM\_001287431

<b>ORF Size:</b>	1119 bp
<b>OTI Disclaimer:</b>	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. <a href="#">More info</a>
<b>OTI Annotation:</b>	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
<b>Components:</b>	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
<b>Reconstitution Method:</b>	<ol style="list-style-type: none"><li>1. Centrifuge at 5,000xg for 5min.</li><li>2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.</li><li>3. Close the tube and incubate for 10 minutes at room temperature.</li><li>4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.</li><li>5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.</li></ol>
<b>RefSeq:</b>	<a href="#">NM_001287431.1</a> , <a href="#">NP_001274360.1</a>
<b>RefSeq Size:</b>	2992 bp
<b>RefSeq ORF:</b>	1122 bp
<b>Locus ID:</b>	27236
<b>UniProt ID:</b>	<a href="#">P53367</a>
<b>Cytogenetics:</b>	4q31.3
<b>MW:</b>	42.2 kDa
<b>Gene Summary:</b>	Putative target protein of ADP-ribosylation factor.[UniProtKB/Swiss-Prot Function]