

## Product datasheet for **RC237831**

### ANGEL2 (NM\_001300758) Human Tagged ORF Clone

#### Product data:

**Product Type:** Expression Plasmids  
**Product Name:** ANGEL2 (NM\_001300758) Human Tagged ORF Clone  
**Tag:** Myc-DDK  
**Symbol:** ANGEL2  
**Synonyms:** Ccr4d; KIAA0759L  
**Vector:** pCMV6-Entry (PS100001)  
**E. coli Selection:** Kanamycin (25 ug/mL)  
**Cell Selection:** Neomycin  
**ORF Nucleotide Sequence:** >RC237831 representing NM\_001300758  
Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC  
GCC**GCGATCGCC**

ATGTCCTATAATACTTTCACAAGATTTACTGGAAGATAACTCTCACCTTTATAGACATTGCCGGCGGC  
CAGTATTACACTGGAGTTTTAGGTTTCCCAATATCTGAAAGAAATTAACATTTTGATGCAGACGCTACT  
TTGTTTGAAGAAGTTCAAGAAGATCATTATGGAGCAGAGATCAGCCAAGTTTGAATCACTGGGTTAT  
CACTGTGAATATAAGATGCGGACAGGAAGGAACTGATGGCTGTGCTATTTGCTTCAAACATTCCAAAT  
TTTCACTCTTGTCAGTGAACCCAGTGAATTCTCCGCCCTGATATTTCTCTGTTGGACAGAGACAATGT  
TGGATTAGTTTTACTCTTACAGCCAAAATTCATATGCTGCCTGCCCTGCAATCTGCGTAGCAAAATACG  
CATCTGTTGTATAATCCAAGGCGAGGTGATTAAGCTGACGCAATTGGCAATGCTACTGGCAGAGATTT  
CCAGTGTGCCACCAGAAAGATGGCAGCTTCTGCCCTATTGTTATGTGTGGTGACTTTAATTCTGTTCC  
TGGTCTCCACTATATAGTTTCATAAAGGAAGGAAAATGAATTATGAAGGACTTCCCATAGGAAAGGTA  
TCTGGCCAGGAACAGTCTTACGGGGACAAAGAATTTATCTATTCCAATTTGGCCCCAAACCTAGGTA  
TCTCACAGAAGTGTGTATGAGGTACAGCAGGTACCAAAAGTAGAAAAGACAGACAGTGTCTGACACA  
AACACAGCTGAAGCAAACAGAGGTCTAGTGACAGCTGAAAAATTTGCTTCAAATTTACAGCACCATTTT  
AGTTTGCATCTGTTTATTCACATTACTTTCTGACACTGGAATTCCAGAAGTGACCACCTGTCATTCCC  
GAAGTGCCATAACTGTGGATTATATTTTCTACTCTGCAGAAAAGGAAGATGTTGCTGGGCACCCAGGAGC  
TGAAGTTGCTTTGGTTGGTGGCTTCAAACCTTCTAGCTAGACTGTCACTTCTTACAGAACAAGACTTATGG  
ACTGTTAATGGACTTCAAACGAAAATAACTCTTCAGATCATCTGCCTTTATTGGCAAAGTTCAGACTTG  
AGCTC

**ACGCGT**ACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT  
ACAAGGATGACGACGATAAGGTTTAA



[View online »](#)

**Protein Sequence:** >RC237831 representing NM\_001300758  
 Red=Cloning site Green=Tags(s)

MSYNILSQDLLEDNSHL YRHCRRPVLHWSFRFPNLIK EIKHFDADVLC LQEVQEDHYGAEIRPSLES LGY  
 HCEYKMRTRGRKPDGCAICFKHSKFSLLSVNPVEFFRPDISLLDRDNVGLVLLLPKIPYAACPAICVANT  
 HLLYNPRRGDIKLTQLAMLLAEISSVAHQKDGSCPIVMCGDFNSVPGSPLYSFIKEGKLNIEGLPIGKV  
 SGQEQSSRGQRILSIPIWPPNLGISQNCVYEVQVPKVEKTDSDLTQTQLKQTEVLVTAEKLSNLQHMF  
 SLSSVYSHYFPDTGIPEVTTCHSRSAITVDYIFYSAEKEDVAGHPGAEVALVGGLKLLARLSLLTEQDLW  
 TVNGLPNNSSDHLPLLAKFRLEL

TRTRPLEQKLISEEDLAANDILDYKDDDDKV


**Restriction Sites:**

Sgfl-MluI

**Cloning Scheme:**



**Plasmid Map:**



**ACCN:** NM\_001300758

<b>ORF Size:</b>	1125 bp
<b>OTI Disclaimer:</b>	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. <a href="#">More info</a>
<b>OTI Annotation:</b>	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
<b>Components:</b>	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
<b>Reconstitution Method:</b>	<ol style="list-style-type: none"><li>1. Centrifuge at 5,000xg for 5min.</li><li>2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.</li><li>3. Close the tube and incubate for 10 minutes at room temperature.</li><li>4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.</li><li>5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.</li></ol>
<b>RefSeq:</b>	<a href="#">NM_001300758.1</a> , <a href="#">NP_001287687.1</a>
<b>RefSeq Size:</b>	4391 bp
<b>RefSeq ORF:</b>	1128 bp
<b>Locus ID:</b>	90806
<b>UniProt ID:</b>	<a href="#">Q5VTE6</a>
<b>Cytogenetics:</b>	1q32.3
<b>MW:</b>	42.6 kDa