

Product datasheet for **RC237840**

TRMT1 (TRMU) (NM_001282785) Human Tagged ORF Clone

Product data:

Product Type: Expression Plasmids
Product Name: TRMT1 (TRMU) (NM_001282785) Human Tagged ORF Clone
Tag: Myc-DDK
Symbol: TRMU
Synonyms: LCAL3; MTO2; MTU1; TRMT; TRMT1
Vector: pCMV6-Entry (PS100001)
E. coli Selection: Kanamycin (25 ug/mL)
Cell Selection: Neomycin
ORF Nucleotide Sequence: >RC237840 representing NM_001282785
Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCC**CGATCGCC**

ATGCAGGCCTTGCAGCAGTCTGTGCGCCCTGTCGGCGGCGTGGACAGCGCCGTGGCCGCGTCTGCTGC
TGAGGCGGAGAGGTTACCAGGTGACAGGGGTGTTATGAAGAAGTGGACTCACTGGATGAACATGGGGT
CTGTACTGCCGACAAAGACTGTGAAGATGCTTACAGAGTTTGCCAGATCTTAGACATCCCTTTCCATCAA
GTGTCCTACGTAAGGAGTATTGGAATGATGTGTTCAAGTACTTTTGAATGAGTATGAAAAAGGAAGGA
CTCCCAATCCTGACATAGTTTGCAACAAGCACATCAAATTTAGTTGCTTTTTTCATTATGCTGTGGATAA
TCTTGGGGCAGATGCCATTGCCACAGGTCACTATGCAAGAAGTCCCTGGGAAGTGAAGAAGTCTTTGAG
CAGAAGCACGTTAAGAAGCCGAAGGGCTTTTCAGAAATCGGTTTGAAGTTAGAAATGCGGTAAAACCTCC
TCCAGGCAGTGACAGCTTAAAGACCAGACCTTCTTTCAGCCAGGTTTCCAGGATGCCCTGAGGAG
AACCATCTCCCTCTGGGGGATTAACGAAAGAGTTTGTAAAGAAAATCGCTGCTGAGAATAGACTTCAT
CATGTGCTTCAAGAAGAGAGCATGGCATGTGTTTCATCGGAAGAGGAATTTGAACATTTCTCTTC
TTCAGTATCTGCAGCCTCGACCTGGTCACTTTATTTCCATAGAAGACAATAAGGTTCTGGGAACACATAA
AGGTTGGTTCTGTATACCTTGGGCCAGAGAGCAACATAGGTGGCCTGAGAGAGCCTGGTACGTGGTG
GAGAAGGACAGCGTCAAGGGTGACGTGTTTGTGGCCCCCGACAGACCACCCAGCCCTGTACAGGGACC
TGCTGAGGACCAGCCGCTGCACTGGATTGCGGAGGAGCCTCCCGCAGCACTGGTCCGGGACAAGATGAT
GGAGTGCCACTTCCGATTCGCCACCAGATGGCACTAGTTTGTGTTCTACAAGGGGACGAGTGCCT
GGGCAGCGGAAGATCCTGCGGCTGGGGCCGTCTGCCTACACGCTCCAGAAGGGCCAGCGCAGAGCTGGG
ATGGCCAC

ACGCGTACGCGGCCGCTCGAGCAGAAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT
ACAAGGATGACGACGATAAGGTTTAA



[View online »](#)

Protein Sequence: >RC237840 representing NM_001282785
 Red=Cloning site Green=Tags(s)

MQALRHVVICALSGGVDSAVAALLRRRGYQVTGVFMKNWDSLDEHGVCTADKDCEDAYRVCQILDIPFHQ
 VSYVKEYWNVDFSDFLNEYEKGRTPNPDIVCNKHIKFSFFFHYAVDNLGADAIATGHYARTSLEDEEVFE
 QKHVKKPEGLFRNRFVNRNAVKLLQAADSFKDQTFLLSQVSDALRRTIFPLGGLKEFVKKIAAENRLH
 HVLQKKE SMG MCF IGKRNFEHFLQLYLQPRPGHF I SIEDNKVLGTHKGWFLYTLGQRANIGGLREPWYVV
 EKDSVKGDVVFVAPRTDHPALYRDLLRTRSRVHWIAEPPAALVRDKMMECHFRRHQMALVCCVLQGGRRV
 GQREDPAAGAVCLHAPEGPAQSWDGH

TRTRPLEQKLISEEDLAANDILDYKDDDDKV


Restriction Sites:

SgfI-MluI

Cloning Scheme:



Plasmid Map:



ACCN: NM_001282785

ORF Size:	1128 bp
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
Reconstitution Method:	<ol style="list-style-type: none">1. Centrifuge at 5,000xg for 5min.2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.3. Close the tube and incubate for 10 minutes at room temperature.4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
RefSeq:	NM_001282785.2
RefSeq Size:	1891 bp
RefSeq ORF:	1131 bp
Locus ID:	55687
UniProt ID:	O75648
Cytogenetics:	22q13.31
MW:	43.4 kDa
Gene Summary:	This nuclear gene encodes a mitochondrial tRNA-modifying enzyme. The encoded protein catalyzes the 2-thiolation of uridine on the wobble positions of tRNA(Lys), tRNA(Glu), and tRNA(Gln), resulting in the formation of 5-taurinomethyl-2-thiouridine moieties. Mutations in this gene may cause transient infantile liver failure. Polymorphisms in this gene may also influence the severity of deafness caused by mitochondrial 12S ribosomal RNA mutations. Alternative splicing results in multiple transcript variants. [provided by RefSeq, Sep 2013]