

Product datasheet for RC237947

ETNK2 (NM_001297760) Human Tagged ORF Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	ETNK2 (NM_001297760) Human Tagged ORF Clone
Tag:	Myc-DDK
Symbol:	ETNK2
Synonyms:	EKI2; HMFT1716
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)
Cell Selection:	Neomycin
ORF Nucleotide Sequence:	>RC237947 representing NM_001297760 Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCC**CGGATCGCC**

ATGGCTGTGCCCCCTTCGGCCCCTCAGCCGCGCGCTCCTTTACCTGAGGAGGCACACGCCTTGCCCGC
AGTGCTCATGGGCATGGAGGAGAAGGCGCGGCCAGCCAGCTGCCGGGAGCCCGGGCCCCCGAG
GGCCGCCCGCGTGCCTACTTCGGCATTCCGTGGACCCGGACGACATCCTTCCCGGGCCCTGCGCCTC
ATCCAGGAGCTGCGGCCGATTGGAAACCCGAGCAAGTTCGGACCAAGCGCTTACGGATGGCATACCA
ACAAGCTGGTGGCCTGCTATGTGGAGGAGGACATGCAGGACTGCGTGTGGTCCGGGTGTATGGGGAGCG
GACGGAGCTGCTGGTGGACCGGGAGAATGAGGTCAGAACTTCCAGCTGCTGCGAGCACACAGCTGTGCC
CCCAAACCTACTGCACCTTCCAGAAATGGGCTGTGCTATGAGTACATGCAGGGTGTGGCCCTGGAGCCTG
AGCACATCCGTGAGCCCCGGCTTTTCAGGTTAATCGCCTAGAAATGGCAAAGATTCATACTATCCACGC
CAACGGCAGCCTGCCAAGCCATCCTCTGGCACAAGATGCACAATTATTCACGCTTGTGAAGAACGAG
ATCAACCCAGCCTTTCTGCAGATGTCCCTAAGGTAGAGGTGTTGGAACGGGAGCTGGCCTGGCTGAAGG
AGCATCTGTCCCAGCTGGAGTCCCCTGTGGTGTTCACAAATGACCTGCTGCAAGAATATCATCTA
TGACAGCATCAAAGGTCACGTGCGGTTCAATTGACTATGAATATGCTGGCTACAACCTACCAAGCTTTTGAC
ATTGGCAACCATTCAATGAGTTTGCAGGCGTGAATGAGGTGGATTACTGCCTGTACCCGCGCGGGAGA
CCCAGCTGCAGTGGCTGCACTACTACCTGCAGGCACAAAAGGGGATGGCCGTGACCCCGAGGAGGTGCA
AAGGCTCTACGTGCAAGTCAACAAGTTTCCCTGGGTCTAGCTGTGTCTTCCACAATGACTGCATCC
CTCCAGTGTGTAGAGTCGGAACAGGCATGGGGAGATTGCCAGGCTGACCCTCTCTGGTCTGTTCCAG
CGCTCTCACTTCTCTGGGCTCTCTGGGCCCTCATCCAGAACCAGTACTCCACCATCGACTT

ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT
ACAAGGATGACGACGATAAGGTTTAA



[View online »](#)


Protein Sequence: >RC237947 representing NM_001297760
 Red=Cloning site Green=Tags(s)

MAVPPSAPQPRA SFHLRRHTPCPQCSWGMEEKAAASASCREPPGPPRAAAVAYFGISVDPDDILPGALRL
 IQELRPHWKPEQVRTKRF TDGITNKL VACYVEEDMQDCVLRVYGERTELLVDRENEVRNFQLLRAHS
 CA PKLYCTFQNGLCYEYMQGVALEPEHIREPRLFRLIALEMAKIHTIHANGSLPKPILWHKMHNFTLVKNE
 INPSLSADVPKVEVLERELAWLKEHLSQLESPVVFCCHNDLLCKNIYDSIKGHVRFIDYEYAGYNYQAFD
 IGNIHFNEFAGVNEVDYCLYPARETQLQWLHYYLQAQKGMVTPREVQRLVYQVNFALGPSCVSSMTAS
 LQCCR VGNRHGEIARLTLSGLFPGVSLLLGSLGPHPEPVLHHRLL


TRTRPLEQKLISEEDLAANDILDYKDDDDKV

Restriction Sites: SgfI-MluI

Cloning Scheme:



Plasmid Map:



ACCN: NM_001297760

ORF Size:	1182 bp
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
Reconstitution Method:	<ol style="list-style-type: none">1. Centrifuge at 5,000xg for 5min.2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.3. Close the tube and incubate for 10 minutes at room temperature.4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
RefSeq:	NM_001297760.2
RefSeq Size:	2744 bp
RefSeq ORF:	1185 bp
Locus ID:	55224
UniProt ID:	Q9NVF9
Cytogenetics:	1q32.1
Protein Families:	Druggable Genome
Protein Pathways:	Glycerophospholipid metabolism, Metabolic pathways
MW:	45.3 kDa
Gene Summary:	The protein encoded by this gene is a member of choline/ethanolamine kinase family which catalyzes the first step of phosphatidylethanolamine (PtdEtn) biosynthesis via the cytidine diphosphate (CDP) ethanolamine pathway. Alternative splicing results in multiple transcript variants. [provided by RefSeq, Jul 2014]