

## Product datasheet for RC237968

### ZNF302 (NM\_001289186) Human Tagged ORF Clone

#### Product data:

Product Type:	Expression Plasmids
Product Name:	ZNF302 (NM_001289186) Human Tagged ORF Clone
Tag:	Myc-DDK
Symbol:	ZNF302
Synonyms:	HSD16; MST154; MSTP154; ZNF135L; ZNF140L; ZNF327
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)
ORF Nucleotide Sequence:	>RC237968 ORF sequence Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC  
GCCCGCATCGCC

ATGTCCTCAGGTGACATTTAGTGATGTGGCTATAGACTTCTCTCATGAAGAGTGGGCATGCCTAGATTCTG  
CTCAGAGGGACTTATACAAGGATGTGATGGTCCAGAATTATGAGAACCTGGTCTCTGTAGGTCTTTCCGT  
AACTAAGCCATATGTGATCATGTTATTGGAGGATGGAAAAGAGCCCTGGATGATGGAGAAAAACTGTCA  
AAAGATTGGGAATCAAGATGGGAAAACAAGGAATTATCAACAAAAGAGGATTTTATGATGAAGATTCAC  
CCCAACAGTAACAATGGAAAAGTTGTAACAAAGTTATGAATTTTCAAATTCATAAAGAATTTGGA  
ATATACAGAATGCGACACATTTAGAAGCACCTTTCATTCAAAGTCTACTCTTTCTGAACCACAAAACAAT  
TCTGCTGAAGGAATTCACACAAATATGATATATTAAGAAGAATTTATCAAAAAGTCAGTTATAAAAA  
GTGAGAGAATAAATGGTGGAAAAGAACTTTTAAATTTCTAATAAAAAGTGGGCAGCCTTCAACCAGAGCAA  
ATCTCTTACCCTTCCCCAGACTTGTAAATAGAGAGAAAACTATACATGCAGTGAATGTGGAAAAGCCTTT  
GGCAAACAGTCAATCCTCAGTCGCCACTGGAGAATTCATACAGGAGAGAAGCCCTATGAATGTCGTGAAT  
GTGGGAAGACTTTTAGCCATGGTTCATCCCTTACACGACATCAGATAAGCCATAGTGGAGAGAAACCTTA  
CAATGCATTGAATGTGGGAAGGCCTTTAGCCATGGCTCATCACTTACTAACCATCAGAGCACTCACAGC  
GGAGAGAAACCGTATGAATGTGAATGTGAAAGTCTTTTAGTGTGTCCCTTCTCATTGAGCATC  
TAAGAATTCATACGCAAGAAAAACGCTATGAGTGTGATATATGTGAAAAGCCTTCAATTCATAGTTGCTC  
TCTCATTACCATCAGAAAAGCCATACTGGAGAGAAGCCTTATGAATGTAGAGAATGTGGAAAAGCTTTC  
TGCTGTAGCTCACACCTTACTCAACATCAAGAATTCACAGTATGAAGAAAAATATGAATGCAACAAAT  
GTCTCAAGGTCTTTAGTAGCTTCTCATTCTTGTTC AACATCAGAGTATTCATACTGAAGAAAAACCGTT  
TGAAGTT

ACGCGTACGCGGCCGCTCGAGCAGAAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT  
ACAAGGATGACGACGATAAGGTTTAA



[View online >](#)

**Protein Sequence:** >RC237968 protein sequence  
Red=Cloning site Green=Tags(s)

MSQVTFSDVAIDFSHEEWACLDSAQRDL YKDV MVQNYENLVSVGLSVTKPYVIMLLEDGKEPMMMEKLS  
 KDWESRWENKELSTKKDIYDEDSPQPVTMEKVVKQSYEF SNSNKNLEYTECDTFRSTFHSTLSEPONN  
 SAEGNSHKYDILKKNLSKKSIVIKSERINGGKLLNSNKSAAFNQSKSL TLPQTCNREKIYTCSECGKAF  
 GKQSILSRHWRIHTGEKPYECRECGKTF SHGSSL TRHQI SHSGEKPYKCI ECGKAF SHGSSL TNHQSTHT  
 GEKPYECMNCGKSF SRVSLLIQHLRIHTQEKRYE CRICGKAF IHSSSLIHHQKSHTGEKPYECRECGKAF  
 CCSSHLTQHRIHSMKKKYECNKCLKVFSSFSFLVQHQSITHEEKPFEV

TRTRPLEQKLISEEDLAANDILDYKDDDDKV

**Chromatograms:** [https://cdn.origene.com/chromatograms/mk6536\\_a07.zip](https://cdn.origene.com/chromatograms/mk6536_a07.zip)

**Restriction Sites:** SgfI-MluI

**Cloning Scheme:**



**ACCN:** NM\_001289186

**ORF Size:** 1197 bp

**OTI Disclaimer:** The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

**OTI Annotation:** This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.

**Components:** The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

**Reconstitution Method:**

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

**RefSeq:** [NM\\_001289186.1](#), [NP\\_001276115.1](#)

**RefSeq Size:** 2586 bp

**RefSeq ORF:** 1200 bp

**Locus ID:** 55900

**UniProt ID:** [Q9NR11](#)


**Cytogenetics:** 19q13.11

**Protein Families:** Transcription Factors

**MW:** 45.9 kDa

**Gene Summary:** This gene encodes a member of the zinc-finger protein family. The encoded protein contains seven C2H2-type zinc fingers and a KRAB domain, but its function has yet to be determined. Alternatively spliced transcript variants have been described. [provided by RefSeq, Mar 2014]

### Product images:



Circular map for RC237968