

Product datasheet for **RC238128**

ATE1 (NM_001288735) Human Tagged ORF Clone

Product data:

Product Type: Expression Plasmids
Product Name: ATE1 (NM_001288735) Human Tagged ORF Clone
Tag: Myc-DDK
Symbol: ATE1
Vector: pCMV6-Entry (PS100001)
E. coli Selection: Kanamycin (25 ug/mL)
Cell Selection: Neomycin
ORF Nucleotide Sequence: >RC238128 representing NM_001288735
Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCC**CGATCGCC**

ATGTTGAAATTTCTAGCTAAAGGGGAGGTTCCAAAGGAAGTTGTGAGGATGAGCCCATGGATTCCACAA
TGGATGATGCTGTTGCGGGTGACTTTGCATTGATAAATAAACTGGATATACAGTGTGATCTTAAACACT
CAGTGATGACATCAAAGAGAGTTTAGAGAGTGAAGGAAAAAATCAAAGAAAGAAGAACCTCAGGAATTA
CTTCAGTCACAAGATTTGCTAGGAGAGAAGTTGGGCTCTGGTGAACCGTCACATTCAGTTAAAGTTCACA
CAGTTCCTAAGCCAGGCAAAGGGGCTGATTTGAGTAAGCCTCCATGTCGAAAAGCAAAGGAAATCCGGAA
AGAAAGGAAAAGGTTAAACTAATGCAGCAGAACCAGCTGGAGAACTTGAGGGTTCCAGGCTCAAGGT
CACCCACCATCTTTGTTCCACCAAAGGCTAAATCCAACCAGCCAAAATCACTCGAAGATTTAATTTTTG
AGTCTTTACCAGAGAATGCATCACACAAGTTAGAGGTGAGGTTAGTACCTGTCTCCTTTGAGGACCCAGA
GTTCAAGTCGTCTTTAGCCAGTCCTTTCTTTGTATGTCAAGTATCAAGTGGCCATACACCAGGATCCA
CCCGATGAATGTGGGAAGACTGAGTTCACAAGATTCCTTTGCAGTTCACCCTTGAGGCAGAGACTCCCC
CTAATGGGCCAGATTGTGGCTATGGCTCCTTTACCAGCAGTACTGGCTTGACGGAAGATCATTGCTGT
GGGGTGATTGACATCTCCAACTGTGTATCATCTGTGATTTGTAACGATCCTGATTATTGTTTT
TTGCTTTGGGCGTCTACTCTGCACTACGAGAAATGCTTTTACTAGGCAGCTTCATGAGAAAATCTCTC
AACTCAGCTATTATTATATGGGTTTCTACATTCATTATGTCCCAAGATGAAATATAAGGGTCAGTATAG
ACCTTCTGATTTGCTGTGCCCTGAGACATATGTTGGGTACCCATTGAGCAATGCCTGCCTTCACTTGAA
AACTCCAAGTACTGCCGTTTCAACCAGGACCCAGAAGCAGTGGATGAGGATCGCAGTACGGAACTGACC
GATTGCAGGTGTTTCAAGAGAGCCATCATGCCTTACGGTGTATAAGAAACAGCAGAAAAGACCCAAAG
TGAGGAGGCTGCTGTTCTGCAGTACGCCAGCCTGGTGGGCAGAAAGTGTCCGAGCGGATGCTGCTGTTC
AGAAAC

ACGCGTACGCGGCCGCTCGAGCAGAAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT
ACAAGGATGACGACGATAAGGTTTAA



[View online »](#)

Protein Sequence: >RC238128 representing NM_001288735
Red=Cloning site Green=Tags(s)


MLKFLAKGEVPGKSGCEDEPMDSTMDAVAGDFALINKLDIQCDLKTLSDDIKESLESEGKNSKKEEPQEL
 LQSQDFVGEKLGSGEPSSHVHTVPKPGKADLSKPPCRKAKEIRKERKRLKLMQONPAGELEGFQAQG
 HPPSLFPPKAKSNQPKSLEDLIFEESLPENASHKLEVRLVPVSFEDPEFKSSFQSFSLYVKYQVAIHQDP
 PDECGKTEFTRFLCSSLPEAETPPNGPDCGYGSFHQQYWLDGKIIAVGVIDILPNCVSSVYLYYDPDYSF
 LSLGVYSALREIAFTRQLHEKTSLSYYMIFYIHSCPMMKYKQYRPSDLLCPETYVWVPIEQCLPSLE
 NSKYCRFNQDPEAVDEDRSTEPDRLQVFKRAIMPYGVYKQKQDPSEEAALVQYASLVGQKCSERMLLFRN

TRTRPLEQKLISEEDLAANDILDYKDDDDKV


Restriction Sites:

SgfI-MluI

Cloning Scheme:



Plasmid Map:



ACCN:	NM_001288735
ORF Size:	1266 bp
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
Reconstitution Method:	<ol style="list-style-type: none">1. Centrifuge at 5,000xg for 5min.2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.3. Close the tube and incubate for 10 minutes at room temperature.4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
RefSeq:	NM_001288735.1 , NP_001275664.1
RefSeq Size:	4803 bp
RefSeq ORF:	1269 bp
Locus ID:	11101
UniProt ID:	O95260
Cytogenetics:	10q26.13
MW:	48.3 kDa
Gene Summary:	This gene encodes an arginyltransferase, an enzyme that is involved in posttranslational conjugation of arginine to N-terminal aspartate or glutamate residues. Conjugation of arginine to the N-terminal aspartate or glutamate targets proteins for ubiquitin-dependent degradation. Alternative splicing results in multiple transcript variants. [provided by RefSeq, Dec 2013]