

Product datasheet for **RC238149**

Aminoadipate aminotransferase (AADAT) (NM_001286683) Human Tagged ORF Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	Aminoadipate aminotransferase (AADAT) (NM_001286683) Human Tagged ORF Clone
Tag:	Myc-DDK
Symbol:	Aminoadipate aminotransferase
Synonyms:	KAT2; KATII; KYAT2
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)



[View online »](#)

ORF Nucleotide Sequence:

>RC238149 representing NM_001286683
 Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
 GCC**CGCATCGCC**

ATGAATTACGCACGGTTCATCACGGCAGCGAGCGCAGCCAGAAACCCCTTCTCCCATCCGGACCATGACTG
 ACATATTGAGCAGAGACCAAAATCGATGATCTCCTTGGCTGGTGGCTTACCAAATCCAAACATGTTTCC
 TTTAAGACTGCCGTAATCACTGTAGAAAATGGAAAGACCATCCAATTTGGAGAAGAGATGATGAAGAGA
 GCACCTCAGTATTCTCCGAGTGTGGAATTCAGAGCTTTTGTCTGGCTAAAACAGTTACAAAATAAAAT
 TGCATAATCCTCCTACCATCCATTACCCACCCAGTCAAGGACAAATGGATCTATGTGTCACATCTGGCAG
 CCAACAAGGTCTTTGTAAGGTGTTTGAATGATCATTAACTCTGGAGATAATGTCCTCCTAGATGAACCT
 GCTTATTCAGGAAGTCTCAAAGTCTGCACCCACTGGGCTGCAACATTATTAATGTTGCCAGTGTGAAA
 GTGGGATTGTTCCAGATCCCTAAGAGACATACTTTCCAGATGGAACCAGAAGATGCAAGAATCCCA
 GAAAAACACCCCAAAATTTCTTTATACTGTTCCAATGGCAACAACCCCTACTGAAAACCTATTAACCACT
 GAACGCAAAAAGGAAATCTATGAGCTTGAAGAAAATATGATTTCTCATAATAGAAGATGATCCTTACT
 ATTTTCTCCAGTTTAAACAGTTCAGGGTACCAACATTTCTTTCCATGGATGTTGATGGACGTGTCATCAG
 AGCTGACTCTTTTTCAAAAATCATTTCCTCTGGGTTGAGAATAGGATTTTAACTGGTCCAAAACCCCTTA
 ATAGAGAGAGTTATTTTACACATACAAGTTTCAACATTGCACCCAGCACTTTTAAACAGCTCATGATAT
 CACAGCTTCTACACGAATGGGGAGAAGAAGTTTCATGGCTCATGTAGACAGGGTTATTGATTTCTATAG
 TAACCAGAAGGATGCAACTACTGGCAGCTGCAGACAAGTGGTTAACTGGTTTGGCAGAATGGCATGTTCT
 GCTGCTGGAATGTTTTTATGGATTAAGTTAAAGGCATTAATGATGTAAGAAGTGAAGAAAAGG
 CCGTTAAGATGGGGTATTAATGCTCCCTGGAAATGCTTTCTACGTCGATAGCTCAGCTCCTAGCCCTTA
 CTGAGAGCATCCTTCTCTCAGCTTCTCCAGAACAGATGGATGTGGCCTTCCAGGTATTAGCACAACTT
 ATAAAAGAAATCTTTA

ACGCGTACGCGGCCGCTCGAGCAGAAACTCATCTCAGAAGAGGATCTGGCAGCAATGATATCCTGGATT
 ACAAGGATGACGACGATAAGGTTTAA

Protein Sequence:

>RC238149 representing NM_001286683
 Red=Cloning site Green=Tags(s)

MNYARFITAASAARNPSPIRTMTDILSRGPKSMISLAGGLPNPNMFPKTAVITVENGKTIQFGEEMMKR
 ALQYSPSAGIPELLSWLKQLQIKLHNPTIHYPPSQGQMDL CVTSGSQGLCKVFEMIINPGDNVLLDEP
 AYSGLQSLHPLGCNIINVASDESGIVPDSL RDILSRWKPEDAKNPQKNTPKFLYTPNGNPTGNSLTS
 ERKKEIYELARKYDFLI IEDDPYYFLQFNKFRVPTFLSMDVDGRVIRADSF SKIISSGLRIGFLTGPKPL
 IERVILHIQVSTLHPSTFNQLMISQLLHEWGEEGFMAHVDRVIDFYSNQKDAILAAADKWL TGLAEWHVP
 AAGMFLWIKVKGINVKELIEEKAVKMGVLMPLGNAFYVDSSAPSPYLRASFSSASPEQMDVAFQVLAQL
 IKESL

TRTRPLEQKLISEEDLAANDILDYKDDDDKV

Chromatograms:

https://cdn.origene.com/chromatograms/mk6276_b06.zip

Restriction Sites:

SgfI-MluI

Cloning Scheme:


ACCN: NM_001286683

ORF Size: 1275 bp

OTI Disclaimer: The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

OTI Annotation: This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

Reconstitution Method:

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

RefSeq: [NM_001286683.1](#), [NP_001273612.1](#)

RefSeq Size: 2193 bp

RefSeq ORF: 1278 bp

Locus ID: 51166

UniProt ID: [Q8N5Z0](#)


Cytogenetics: 4q33

Protein Pathways: Lysine biosynthesis, Lysine degradation, Metabolic pathways, Tryptophan metabolism

MW: 47.4 kDa

Gene Summary: This gene encodes a protein that is highly similar to mouse and rat kynurenine aminotransferase II. The rat protein is a homodimer with two transaminase activities. One activity is the transamination of alpha-amino adipic acid, a final step in the saccharopine pathway which is the major pathway for L-lysine catabolism. The other activity involves the transamination of kynurenine to produce kynurenine acid, the precursor of kynurenic acid which has neuroprotective properties. Several transcript variants encoding two different isoforms have been found for this gene. [provided by RefSeq, Nov 2013]

Product images:



Circular map for RC238149