

Product datasheet for **RC238158**

TTC39A (NM_001297662) Human Tagged ORF Clone

Product data:

Product Type: Expression Plasmids
Product Name: TTC39A (NM_001297662) Human Tagged ORF Clone
Tag: Myc-DDK
Symbol: TTC39A
Synonyms: C1orf34; DEME-6
Vector: pCMV6-Entry (PS100001)
E. coli Selection: Kanamycin (25 ug/mL)
Cell Selection: Neomycin
ORF Nucleotide Sequence: >RC238158 representing NM_001297662
Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGCCGGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCC**CGGATCGCC**

ATGACCTCGGCTGGCGGGCCCCAGGAGCCCTGCCCGGGGACTCCTGAGAGCAGCCTCCATGAGGCC
TGGACCAGTGCATGACCGCCCTGGACCTCTTCTCACCAACCAGTTCTCAGAAGCACTCAGTACCTCAA
GCCAGAACCAAGGAAAGCATGTACCACTCACTGACATATGCCACCATCCTGGAGATGCAGGCCATGATG
ACCTTTGACCTCAGGACATCCTGCTTGCCGGCAACATGATGAAGGAGGCACAGATGCTGTGTGAGAGGC
ACCGGAGGAAGTCTTCTGTAACAGATTCTTCAGCAGCCTGGTGAACCGCCCCACGCTGGGCCAATTAC
TGAAGAGGAAATCCACGCTGAGGTCTGCTATGCAGAGTGCCTGCTGCAGCGAGCAGCCCTGACCTTCTG
CAGGGTTCCTCACACGGAGGGGAGTCAAGGAGCCTTGCATGATCCCTCTCACGCCTGCAGCTGCC
CACCTGGGCCAGGCCGTCAGCATCTTTTCTCCTGCAGGACGAGAACATGGTGAAGCTTCAAAAGGCGG
CATCAAAGTTCGAAACAGCTACCAGACCTACAAGGAGCTGGACAGCCTTGTTCAGTCCTCACAATACTGC
AAGGGTGAGAACCACCCGCACTTTGAAGGAGGAGTGAAGCTTGGTGTAGGGGCCCTCAACCTGACACTGT
CCATGCTTCTACTAGGATCCTGAGGCTGTTGGAGTTTGTGGGTTTTTCAGGAAACAAGGACTATGGGCT
GCTGCAGCTGGAGGAGGAGCGTCAGGGCACAGCTTCCGCTCTGTGCTCTGTGTCATGCTCCTGCTGTGC
TACCACACCTTCTCACCTTCGTGCTCGGTAAGGAGTCAACATCGAGGAGGCCGAGAAGCTTGA
AGCCCTACCTGAACCGGTACCCTAAGGGTGCCATCTTCTGTTCTTTGCAGGGAGGATTGAAGTCATTAA
AGGCAACATTGATGCAGCCATCCGGCTTTTCGAGGAGTGTGTGAGGCCAGCAGCACTGGAAGCAGTTC
CACCACATGTGCTACTGGGAGCTGATGTGGTCTTACCTACAAGGGCCAGTGAAGATGCTCCTACTTCT
ACGCCGACCTGCTCAGCAAGGAGAACTGCTGGTCCAAGGCCACCTACATTTACATGAAGGCCGCTACCT
CAGCATGTTTGGGAAGGAGGACCACAAGCCGTTTCGGGACGACGAAGTGAATTATTTTCGGTGGTCCCA
CATCGTGTCCCACTGTG

ACGCGTACGCGGGCCGCTCGAGCAGAAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT
ACAAGGATGACGACGATAAGGTTTAA



Protein Sequence: >RC238158 representing NM_001297662
Red=Cloning site Green=Tags(s)


MTSAGGAPGALPAGTPPESSLHEALDQCMTALDLFLTNQFSEALSYLKPRTKESMYHSLTYATILEMQAMM
 TFDPDQDILLAGNMMKEAQMQLCQRHRRKSSVTDSSSLVNRPTLGQFTEEEIHAEVCYAECLLQRAALTFL
 QGSSHGGAVRPRALHDP SHACSCPPGGRQHLFLLQDENMVSFIKGGIKVRNSYQTYKELDSL VQSSQYC
 KGENHPHFEGGVKLVGVAFNLTL SMLPTRLRLLEFVGFSGNKDYGLLQLEEGASGHSFRS V LCVMLLLC
 YHTFLTFVLGTGNV NIEEA EKLLKPYLNRYPKGAIFFLFFAGRIEVIKGNIDAAIRRFEECCEAQQHWKQF
 HHMCYWELMWCFTYKQWKMSYFYADLLSKENCWSKATYIYMKAAYLSMFGKEDHKPF GDDEVELFRWVP
 HRAPT V

TRTRPLEQKLISEEDLAANDILDYKDDDDKV


Restriction Sites:

SgfI-MluI

Cloning Scheme:



Plasmid Map:



ACCN: NM_001297662

ORF Size:	1278 bp
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
Reconstitution Method:	<ol style="list-style-type: none">1. Centrifuge at 5,000xg for 5min.2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.3. Close the tube and incubate for 10 minutes at room temperature.4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
RefSeq:	NM_001297662.2
RefSeq Size:	1892 bp
RefSeq ORF:	1281 bp
Locus ID:	22996
Cytogenetics:	1p32.3
MW:	48.8 kDa