

Product datasheet for **RC238181**

Aminoadipate aminotransferase (AADAT) (NM_001286682) Human Tagged ORF Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	Aminoadipate aminotransferase (AADAT) (NM_001286682) Human Tagged ORF Clone
Tag:	Myc-DDK
Symbol:	AADAT
Synonyms:	KAT2; KATII; KYAT2
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)
Cell Selection:	Neomycin
ORF Nucleotide Sequence:	>RC238181 representing NM_001286682 Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCC**CGCATCGCC**

ATGAATTACGCACGGTTCATCACGGCAGCGAGCGCAGCCAGAAACCCCTTCTCCATCCGGACCATGAGTG
AGAAACGGGCTGACATATTGAGCAGAGGACCAAAATCGATGATCTCCTGGCTGGTGGCTTACCAATCC
AAACATGTTTCTTTAAGACTGCGTAATCACTGTAGAAAATGGAAGACCATCCAATTTGGAGAAGAG
ATGATGAAGAGAGCACTTCAGTATTCTCCGAGTCTGGAATCCAGAGCTTTTGTCTGGCTAAAACAGT
TACAAATAAAATGCATAATCCTCCTACCATCCATTACCCACCCAGTCAAGGACAATGGATCTATGTGT
CACATCTGGCAGCCAACAAGGTCTTTGTAAGGTGTTTGAATGATCATTAACTCTGGAGATAATGTCCTC
CTAGATGAACCTGCTTATTCAGGAACCTCTCAAAGTCTGCACCCACTGGGCTGCAACATTATTAATGTTG
CCAGTGATGAAAGTGGGATTGTTCCAGATTCCTAAGAGACATACTTTCCAGATGGAAACCAGAAGATGC
AAAGAATCCCCAGAAAAACCCCCAAATTTCTTTATACTGTTCCAAATGGCAACAACCCCTACTGGAAC
TCATTAACCAGTGAACGCAAAAAGGAAATCTATGAGCTTGAAGAAAATATGATTTCCCTATAAGAA
ATGATCCTTACTATTTTCTCCAGTTTAAACAAGTTCAGGGTACCAACATTTCTTTCCATGGATGTTGATGG
ACGTGTCATCAGAGCTGACTCTTTTCAAAAATCATTTCCTCTGGGTTGAGAATAGGATTTTAACTGTT
CCAAAACCTTAATAGAGAGAGTTATTTTACACATACAAGTTTCAACATTGCACCCAGCACTTTTAAAC
AGCTCATGATATCACAGCTTCTACACGAATGGGAGAAGAAGTTTTCATGGCTCATGTAGACAGGGTTAT
TGATTTCTATAGTAACCAGAAGGATGCAATACTGGCAGCTGCAGACAAGTGGTTAACTGGTTTGGCAGAA
TGGCATGTTCTGCTGCTGGAATGTTTTTATGGATTAAGTTAAAGGCATTAATGATGTAAGAAGTGA
TTGAAGAAAAGGCCGTTAAGATGGGGTATTAATGCTCCCTGGAATGCTTTCTACGTCGATAGCTCAGC
TCCTAGCCCTTACTTGAGAGCATCCTTCTTTCAGCTTCTCCAGAACAGATGGATGTGCCTTCCAGGTA
TTAGCACAACTATAAAGAATCTTTA

ACGCGTACGCGGCCGCTCGAGCAGAAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT
ACAAGGATGACGACGATAAGGTTTAA



Protein Sequence: >RC238181 representing NM_001286682
Red=Cloning site Green=Tags(s)


MNYARFITAASAARNPSPIRTMSEKRADILSRGPKSMISLAGGLPNPNMFPFKTAVITVENGKTIQFGEE
 MMKRALQYSPSAGIPELLSWLKQLQIKLHNPPTIHYPPSQGMDCVTSQSQGLCKVFEMIINPGDNVL
 LDEPAYSGTLQSLHPLGCNIINVASDESGIVPDSLRLDLSRWKPEDAKNPQKNTPKFLYTPNGNPTGN
 SLTSEKKEIYELARKYDFLIEDDPYYFLQFNKFRVPTFLSMDVDGRVIRADSFSKIISSGLRIGFLTG
 PKPLIERVILHIQVSTLHPSTFNQLMISQLLHEWGEEGFMHVDRVIDFYSNQKDAILAAADKWLTLGLAE
 WHVPAAGMFLWIKVKGINDVKELIEEKAVKMGVLMPLPGNAFYVDSAPSPLYLRASFSSASPEQMDVAFQV
 LAQLIKESL

TRTRPLEQKLISEEDLANDILDYKDDDDKV


Restriction Sites:

Sgfl-MluI

Cloning Scheme:



Plasmid Map:



ACCN: NM_001286682

ORF Size:	1287 bp
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
Reconstitution Method:	<ol style="list-style-type: none">1. Centrifuge at 5,000xg for 5min.2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.3. Close the tube and incubate for 10 minutes at room temperature.4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
RefSeq:	NM_001286682.1 , NP_001273611.1
RefSeq Size:	2120 bp
RefSeq ORF:	1290 bp
Locus ID:	51166
UniProt ID:	Q8N5Z0
Cytogenetics:	4q33
Protein Pathways:	Lysine biosynthesis, Lysine degradation, Metabolic pathways, Tryptophan metabolism
MW:	48.3 kDa
Gene Summary:	This gene encodes a protein that is highly similar to mouse and rat kynurenine aminotransferase II. The rat protein is a homodimer with two transaminase activities. One activity is the transamination of alpha-amino adipic acid, a final step in the saccharopine pathway which is the major pathway for L-lysine catabolism. The other activity involves the transamination of kynurenine to produce kynurenine acid, the precursor of kynurenic acid which has neuroprotective properties. Several transcript variants encoding two different isoforms have been found for this gene. [provided by RefSeq, Nov 2013]