

Product datasheet for RC238208

GALK2 (NM_001289031) Human Tagged ORF Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	GALK2 (NM_001289031) Human Tagged ORF Clone
Tag:	Myc-DDK
Symbol:	GALK2
Synonyms:	GK2
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)
Cell Selection:	Neomycin
ORF Nucleotide Sequence:	>RC238208 representing NM_001289031 Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCC**GCGATCGC**C

ATGTTTAACTCCAAGTTTGGATCTATCCCAAGTTTTATGTTTCGAGCACCAGGAAGAGTCAACATAATAG
GAGAGCATATAGATTATTGTGGATATTCTGTTCTTCTATGGCTGTAGAACAAGATGTGCTAATAGCTGT
AGAACCTGTGAAAACGTACGCTCTCCAAGTGGCCAATACAAATCCCTTGTATCCGGACTTCAGTACTAGT
GCTAATAACATCCAGATTGATAAAACCAAGCCTTTGTGGCACAACATTTCTTATGTGGACTTAAAGGAA
TTCAGGAACACTTTGGTCTTAGTAACCTGACTGGAATGAACTGCCTGGTAGATGGAAATATCCCACCAAG
TTCTGGCCTCTCCAGCTCCAGTGCTTTGGTCTGTTGTGCTGGCTTGGTGCAGCTCACAGTGTGGGAAGG
AATCTATCCAAGTGGAAGTGCAGAAATCTGTGCCAAGAGTGAGCGTTACATTGGCACTGAAGGAGGAG
GCATGGACCAGTCTATATCATTCTTGCAGAAGAAGGAAGTGCCTGATAGAATTTAGTCTCTGAG
GGCAACCGATGTAAAACCTCCAAGTGGAGCAGTGTGTTGTGATTGCCAACAGTTGTGTGGAGATGAATAAG
GCAGCAACTTCCCATTTCAATATCAGGGTGTGGAGTGTGGCTGGCTGCGAAGCTCCTGGCTAAATACA
AAAGCTTGCAATGGGACAAAGTACTGAGGCTGGAGGAGTGCAGGCTAAACTAGGGATTAGTCTAGAAGA
AATGCTGTTGGTACAGAAGATGCCCTTATCCTGAACCTATAACCCTGAGGAGATCTGCAGGTGTCTG
GGAATTAGCCTGGAGGAAGTCCGAACCCAAATCCTGAGTCCAAACACTCAAGATGTGCTCATCTTCAAAC
TCTATCAGCGGGCAAAGCATGTGTACAGCGAGGCTGCGCGAGTGTCCAGTTTAAAGAAGATATGTGAAGA
AGCACCTGAAAACATGGTCCAGCTGCTGGGAGAGTTGATGAACCAGAGCCACATGAGCTGCCGGGACATG
TATGAGTGCAGCTGCCCGAGCTGGATCAGCTGGTGGACATCTGTGCGAAGTTTGGGGCTCAAGGGTCCAC
GACTTACTGGAGCAGGATGGGGAGGCTGCACAGTATCAATGGTACCTGCGGACAAGCTGCCAGCTTTCT
AGCAAAATGTGCACAAAGCTTATTACCAGAGGAGTGTGGAAGCTTAGCACCGGAGAAGCAAAAGTTTGT
GCTACCAAACCTGGAGGTGGGGCTTTGGTTTTGCTTGAGGCC

AG**CGGACCG**ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCC
TGGATTACAAGGATGACGACGATAAGGTTTAA



Protein Sequence: >RC238208 representing NM_001289031
Red=Cloning site Green=Tags(s)


MFNSKFGSIPKFYVRAPGRVNIIGEHIIDYCGYSVLPMAVEQDVLIAVEPVKTYALQLANTNPLYPDFSTS
 ANNIQIDKTKPLWHNYFLCGLKGIQEHLGNLTGMNCLVDGNIPPSSGLSSSSALVCCAGLVTLTVLGR
 NLSKVELAEICAKSERYIGTEGGMDQISFLAEEGTAKLIEFSPLRATDVKLP SGAVFVIANSVEMNK
 AATSHFNIRVMCRLAAKLLAKYKSLQWDKVLRL EEVQAKLGISLEEMLLVTEDALHPEPYNPEICRCL
 GISLEELRTQILSPNTQDVLIFKLYQRAKHVYSEAAARVLQFKICEEAPENMVQLLGELMNQSHMSCRDM
 YECSCPELDQLVDICRFGAQGSRLTGAGWGGCTVSMVPADKLP SFLANVHKAYQRSDGSLAPEKQSLF
 ATKPGGALVLLLEA

SGPTRTRPLEQKLI SEEDLAANDILDYKDDDDKV


Restriction Sites:

SgfI-RsrII

Cloning Scheme:



Plasmid Map:



ACCN: NM_001289031

ORF Size:	1302 bp
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
Reconstitution Method:	<ol style="list-style-type: none">1. Centrifuge at 5,000xg for 5min.2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.3. Close the tube and incubate for 10 minutes at room temperature.4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
RefSeq:	NM_001289031.1 , NP_001275960.1
RefSeq Size:	3313 bp
RefSeq ORF:	1305 bp
Locus ID:	2585
Cytogenetics:	15q21.1-q21.2
Protein Families:	Druggable Genome
Protein Pathways:	Amino sugar and nucleotide sugar metabolism, Galactose metabolism, Metabolic pathways
MW:	48.1 kDa
Gene Summary:	This gene encodes a highly efficient N-acetylgalactosamine (GalNAc) kinase, which has galactokinase activity when galactose is present at high concentrations. The encoded protein is a member of the GHMP kinase family. Alternative splicing results in multiple transcript variants. [provided by RefSeq, Jun 2017]