

Product datasheet for **RC238276**

ZNF718 (NM_001289930) Human Tagged ORF Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	ZNF718 (NM_001289930) Human Tagged ORF Clone
Tag:	Myc-DDK
Symbol:	ZNF718
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)
Cell Selection:	Neomycin
ORF Nucleotide Sequence:	>RC238276 representing NM_001289930 Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCC**CGATCGCC**

ATGTTGGAGAACTACAGAAACCTGGTCTCCCTGGGTGTTAGTATCTCTAACCCAGACCTGGTCACCAGTC
TGGAGCAAAGAAAAGGCCCTACAATTTGAAGATACATGAAACAGCAGCCAGACCCCCAGCTGTGTGTTCT
CATTTCACCCAAAACCTTTGGACAGTGCAGGGCATAGAAGATTCATTCCACAACTTATACAAAAGGA
CATGAGAAACGTGGACATGAGAATTTAAGAAAACTTGTAAAAGTATAAATGAGTGAAGGTGCAGAAAG
GTGGTTATAATAGAATTAACCAATGCTTATTAACCTACCCAGAAAAAACAATTCAATCTAATATATGTGT
CAAAGTTTTTCATAAATTTCAAATTCAAAAGATAAGATAAGATAACTGGAGATAAAACCTTTAAAG
TGTAAAGAATGTGGCAAATCATTTCACGTGCTCTCACGCCTAACTCAACACAAAAGAATTCATACTGGAG
AGAACCCTACACATGTGAAGAATGTGGCAAAGCCTTTAATTGGTCCTCAATTCTACTAAACATAAGAG
AATTCATGCCAGAGAGAAATTCACAAGTGTGAAGAATGTGGTAAAGGCTTTACTCGTCCTCACACCTT
ACTAAACATAAGAGAATTCATACTGGAGAGAAACCTACATATGTAAAAATGTGGTAAAGCTTTTAACC
AATCCTCAACCCTTAATTTACATAAGAGAATTCATTCTGCACAAAATACTACAAATGTGAAGAATGTGG
TAAAGCCTTTAAGTGGTCCTCATCCCTAATGAACATAAGAGAATTCATGCTGGAGAGAAACCTTCTCA
TGCGAAGAATGTGGCAATGTCTTTACCACATCCTCAGACTTTGCTAAACATAAGAGAATTCATACAGGAG
AGAAACCTTACAAATGTGAAGAATGTGGAAAATCCTTTAATAGGTCCACAACCTTTACGACACATAAGAG
AATCCATACTGGAGAGAAACCTACACATGTGAAGAATGTGGAAAAGCCTTTAATTGGTCCTCAACCCTT
AATGTACACAAGAGAATTCACCTGGAAAAATCCCTACAAATGTGAAGATTGTGGCAAAGCCTTTAAAG
TGTTTGCAAACCTGCATAATCATAAGAAAATTCATACTGGAGAGAAACCTACATATGTAACAATGTGG
CAAAGCCTTTAAACAGTCCTCACACTTGAATAAACATAAGAAAATTCACACTGTAGATAAACCTACAAA
TGTAAGAATGCGGGAAGCTTTAAGCAGTACTCCAACCTTCTCAACATAAGAGAATTCATACTGGAG
GAAAATTT

ACGCGTACGCGGCCGCTCGAGCAGAAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT
ACAAGGATGACGACGATAAGGTTTAA



[View online »](#)

Protein Sequence: >RC238276 representing NM_001289930
 Red=Cloning site Green=Tags(s)


MLENYRNLVSLGVSISNPDLVTSLEQRKEPYNLKIHETAARPPAVCSHFTQNLWTVQGIEDSFHKLIPKG
 HEKRGHENLRKTCKSINECKVQKGGYNRINQCLLTQKKTIQSNICVKVFHKFSNSNKDKIRYTGDKTFK
 CKECGKSFHVL SRLTQHKRIHTGENPYTCEECKAFNWSSILTKHKRIHAREKFYKCEECKGFTRSSHL
 TKHKRIHTGEKPYICEKCGKAFNQSSSTLNLHKRIHSAQKYKCEECKGAFKWSSSLNEHKRIHAGEKPFSS
 CEECGNVFTTSSDFAKHKRIHTGEKPYKCEECKGKSFNRSTTLTTHKRIHTGEKPYTCEECKAFNWSSSTL
 NVHKRIHSGKNPYKCEDCGKAFKVFANLHNHKKIHTGEKPYICKQCGKAFKQSSHLNKHKKIHTVDKPYK
 CKECGKAFKQYSNLPQHKRTHHTGGKF

TRTRPLEQKLISEEDLAANDILDYKDDDDKV


Restriction Sites:

SgfI-MluI

Cloning Scheme:



Plasmid Map:



ACCN:	NM_001289930
ORF Size:	1338 bp
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
Reconstitution Method:	<ol style="list-style-type: none">1. Centrifuge at 5,000xg for 5min.2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.3. Close the tube and incubate for 10 minutes at room temperature.4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
RefSeq:	NM_001289930.1 , NP_001276859.1
RefSeq Size:	4292 bp
RefSeq ORF:	1341 bp
Locus ID:	255403
UniProt ID:	Q3SXZ3
Cytogenetics:	4p16.3
MW:	52 kDa
Gene Summary:	May be involved in transcriptional regulation.[UniProtKB/Swiss-Prot Function]