

Product datasheet for **RC238362**

TBC1D19 (NM_001292054) Human Tagged ORF Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	TBC1D19 (NM_001292054) Human Tagged ORF Clone
Tag:	Myc-DDK
Symbol:	TBC1D19
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)



[View online »](#)

ORF Nucleotide Sequence:

>RC238362 representing NM_001292054
 Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
 GCC**CGCATCGC**

ATGTTGCAGGAGGAGTCGGACCTCTCTCATTATTGCCAGATAGTCCAAAAGCTCAAGGGCTCCAATT
 TGTA CTCTCAGCTGGAACGGCAGGCCTGGGGAAGTTGGGAAAAAGAATTTGAAGAGTTTAAATAGTAT
 GTGCACTGAACTGAGTATCCCACTGGCAGCAAAGAGGCCAGTTGGAGAACAGAAAGAACTCTTAATAAA
 TGAATGAAATGGGAAGTATGAACAGATTTAAGCCTTTTCAGACCTGTTATGCACCTAAGGATTTTC
 TTGAGGTATTAATTAATCTTCGCAACCCAAATTATGAAAACGGTGATTCTCTAGTTTCAGGACTCATTT
 GGGTTTAATTCAAGTCCCACTGAAAGTAAAGACATCCCTGAATTGAAAGAATGCTTTGTGGAAGTTGGC
 TTAATATAGGACAAGTGGGTATAGATGATTCTACACAAGTGCCTCCTGAACTTTTGAAAATGAGCATG
 TACGTATTGGGCAAAAAGTTCTAGCAGAACAGATAGTGTCTGCTCAACAGTACATCAGACAAGGAAG
 TCCACGGCACTGAGAGCTGAATTGTGGGCTCTCATTTTGAATATTTCCAGCCAACTGAGGATGTCTTG
 TATTATGAGCAGCTTAAGACCAATGTGATACAACATGACCTTTTGGTGGACAGTCTAATCTATAAAGATG
 TGAAGTTGACAGCAAGCAATGATGATTATTATTTGTATTTGAAGATTATTTATATCAGGTATTACTTTG
 TTTTCCCGGGATACATCTGTGTTGAGTCACTTTGCATTCAACAGTGCCTCGCCACCAAAATCATACATA
 AGAGGAAAAGTAGGACTGGAAGAATATGCTGTCTTTACCCACCAAAATGGTGTTATCCCTTTTCATGGAT
 TTTCAATGATGTTGCACCACTTTGTTTTCTATACCATGAACCTTCCAAATTGTATCAGATATCCGTGA
 GATGATGTGCGTTTTTTCTTCAGACTCCATTCATCTCTCTCATCCTTCTGGTATTGTGCACTCTGT
 CTGCTGTTTGAAGTCTTCTTCAAACCTATCTTCCCAACTCTTTATCATCTACGAGAAATTGGGGCTC
 AACCACTTCGCATATCATTTAAGTGGATGGTTCGAGCTTCTCTGGATACTTAGCTACAGATCAGCTCTT
 GCTTTTATGGGATAGAATCCTAGGATACAACCTCTGGAATTCTTGCTGTGCTGGCAGCTGCCGTGTTT
 GCTTTCGAGCAGTGAACCTGATGGAGGTGACATCACTGGCTGCAGCTGAAGCAGTTCTGCTGACCTTT
 CTACTTTAAAAGTTATGCCTCTTCTTCAGATTTTTCTGTTTGCTACTGTCACC

ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT
 ACAAGGATGACGACGATAAGGTTTAA

Protein Sequence:

>RC238362 representing NM_001292054
 Red=Cloning site Green=Tags(s)


MLQEESDLSLIIAQIVQKLKGSNLYSQLERQAWGSWEKRILKSLNSMCTELSIPLARKRPVGEQKELLNK
 WNEMGTDEPDLSLFRPVYAPKDFLEVLINLRNPYENGDSLFRTHLGLIQVPLKVVDIPELKECFVELG
 LNIGQLGIDDSTQVPPELFENEHVRIGQKVLAEQDSAAAQYIRQGSPTALRAELWALILNISSQPEDVL
 YYEQLKTNVIQHDLLVDSLIIYKDVKLTA SDDYFVFEDYLYQVLLCFSRDTSVLSHFAFNSASPPKSYI
 RGKLGLEEYAVFYPPNGVIPFHGFSMYVAPLCFLYHEPSKLYQIFREMYVRRFFRLHSISSHPSGIVSLC
 LLFETLLQTYLPQLFYHLREIGAQLRISFKWMVRAFSGYLATDQLLLLWDRILGYNSLEILAVLAAAVF
 AFRVNLMEVTSLAAAEAVLADLSTLKVMPLLQIFL FATVT

TRTRPLEQKLISEEDLAANDILDYKDDDDKV

Restriction Sites:

SgfI-MluI


Cloning Scheme:



ACCN:	NM_001292054
ORF Size:	1383 bp
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
Reconstitution Method:	1. Centrifuge at 5,000xg for 5min. 2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA. 3. Close the tube and incubate for 10 minutes at room temperature. 4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom. 5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
Note:	Plasmids are not sterile. For experiments where strict sterility is required, filtration with 0.22um filter is required.
RefSeq:	NM_001292054.1 , NP_001278983.1
RefSeq Size:	2922 bp
RefSeq ORF:	1386 bp
Locus ID:	55296

UniProt ID: Q8N5T2
Cytogenetics: 4p15.2
MW: 53.2 kDa
Gene Summary: May act as a GTPase-activating protein for Rab family protein(s).[UniProtKB/Swiss-Prot Function]

Product images:



Circular map for RC238362