

Product datasheet for **RC238410**

LMBR1L (NM_001300751) Human Tagged ORF Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	LMBR1L (NM_001300751) Human Tagged ORF Clone
Tag:	Myc-DDK
Symbol:	LMBR1L
Synonyms:	LIMR
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)



[View online »](#)

ORF Nucleotide Sequence:

>RC238410 representing NM_001300751
 Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
 GCC**CGATCGCC**

ATGGAAGCACCTGACTACGAAGTCTATCCGTGCGAGAACAGCTATTCCACGAGAGGATCCGCGAGTGTA
 TTATATCAACACTTCTGTTTGAACACTGTACATCCTCTGCCACATCTTCTGACCCGCTTCAAGAAGCC
 TGCTGAGTTCACCACAGTGGATGATGAAGATGCCACCGTCAACAAGATTGCGCTCGAGCTGTGCACCTTT
 ACCCTGGCAATTGCCCTGGGTGCTGTCTGCTCCTGCCCTTCTCCATCATCAGCAATGAGGTGCTGCTCT
 CCCTGCCCTCGGAATACTACATCCAGTGGCTCAACGGCTCCCTCATCCATGGCCTCTGGAACCTTGTTTT
 TCTCTTCTCAACCTGTCCCTCATCTTCTCATGCCCTTTGCATATTTCTTCACTGAGTCTGAGGGCTTT
 GCTGGCTCCAGAAAGGTGCTCTGGGCCGGTCTATGAGACAGTGGTATGTTGATGCTCCTCACTCTGC
 TGGTGCTAGGTATGGTGTGGGTGGCATCAGCCATTGTGGACAAGAACAAGGCCAACAGAGAGTCACTCTA
 TGACTTTTGGGAGTACTATCTCCCTACCTCTACTCATGCATCTCCTTCTTGGGTTCTGCTGCTCCTG
 GTGTGTACTCCACTGGGTCTCGCCCGCATGTTCTCCGTCAGTGGGAAGCTGCTAGTCAAGCCCCGGCTGC
 TGAAGACCTGGAGGAGCAGCTGTACTGCTCAGCCTTTGAGGAGGCAGCCCTGACCCGCAGGATCTGTAA
 TCCTACTTCTGCTGGCTGCCTTTAGACATGGAGCTGCTACACAGACAGGTCCTGGCTCTGCAGACACAG
 AGGGTCTGCTGGAGAAGAGGCGGAAGGCTTACGCTGGCAACGGAACCTGGGCTACCCCTGGCTATGC
 TGTGCTTGTGCTGCTGACGGGCTGTCTGTGCTCATTGTGGCCATCCACATCCTGGAGCTGCTCATCGA
 TGAGGCTGCCATGCCCGAGGCATGCAGTTGTACTCATCTTTTACCTAATGGTGTCTCAGTTGTGGGC
 TTCTATAGCTCTCACTCTTCCGGAGCCTGCGGCCAGATGGCAGCACTGCCATGACGCAGATAATTG
 GGAAGTGTGTCTGTCTCTGGTCTTAAGTCAAGTCAAGTCTCTGTCTTCTCTGAACCTGGGGCTCACTGC
 CTTTGACCTGCTGGGTGACTTTGGACGCTTCAACTGGCTGGGCAATTTCTACATTGTGTTCTCTACAAC
 GCAGCCTTTGCAGGCCTCACCACACTCTGTCTGGTGAAGACCTTCACTGCAGCTGTGCGGCAGAGCTGA
 TCCGGGCTTTGGGCTGGACAGACTGCCGCTGCCGCTCTCCGTTTCCCCAGGCATCTAGGAAGACCCA
 GCACCAAG

ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCTGGATT
 ACAAGGATGACGACGATAAGGTTTAA

Protein Sequence:


>RC238410 representing NM_001300751
 Red=Cloning site Green=Tags(s)

MEAPDYEVLVSVREQLFHERIRECIISTLLFATLYILCHIFLTRFKKPAEFTTVDDDEDATVNKIALELCTF
 TLAIALGAVLLLPFSIISNEVLLSLPRNYIQLNGSLIHGLWNLVFLFSNLSLIFLMPFAYFFTESEGF
 AGSRKGVLRVYETVVMLMLLTLLVLGMVWVASAIVDKNKANRESLYDFWEYYPYLYSCISFLGLVLLLL
 VCTPLGLARMFVTKLLVKPRLLLEDLEEQLYCSAFEEAALTRRICNPTSCWLPLDMELLHRQVLALQTQ
 RVLLEKRRKASAWQRNLGYPLAMLCCLVLTGLSVLIVAIHILELLIDEAAMPGRMQVVLIFYLMVSSVVG
 FYSSPLFRSLRPRWHDAMTQIIIGNCVCLLVSSALPVFSRTLGLTRFDLLGDFGRFNWLGNFYIVFLYN
 AAFAGLTTLCLVKFTAAVRAELIRAFGLDRLPLPVSGFPQASRKTQHQ

TRTRPLE**QKL**ISEEDLAANDILDYKDDDDKV

Restriction Sites:

SgfI-MluI

Cloning Scheme:


ACCN: NM_001300751

ORF Size: 1407 bp

OTI Disclaimer: The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

OTI Annotation: This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

Reconstitution Method:

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

Note: Plasmids are not sterile. For experiments where strict sterility is required, filtration with 0.22um filter is required.

RefSeq: [NM_001300751.1](#), [NP_001287680.1](#)


RefSeq Size: 2301 bp

RefSeq ORF: 1410 bp

Locus ID: 55716

UniProt ID: [Q6UX01](#)
 Cytogenetics: 12q13.12
 Protein Families: Druggable Genome, Transmembrane
 MW: 53.7 kDa
 Gene Summary: Probable LCN1 receptor. May mediate LCN1 endocytosis.[UniProtKB/Swiss-Prot Function]

Product images:



Circular map for RC238410