

Product datasheet for **RC238512**

SLC38A1 (NM_001278388) Human Tagged ORF Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	SLC38A1 (NM_001278388) Human Tagged ORF Clone
Tag:	Myc-DDK
Symbol:	SLC38A1
Synonyms:	ATA1; NAT2; SAT1; SNAT1
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)



[View online »](#)

ORF Nucleotide Sequence:

>RC238512 ORF sequence
 Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
 GCC**GCGATCGCC**

ATGATGCATTTCAAAGTGGACTCGAATTAAGTGGTTCGAAAACATGACAGTGCCCGAGGATGATAACA
 TTAGCAATGACTCCAATGATTTACCGAAGTAGAAAATGGTCAGATAAATAGCAAGTTTATTCTGATCG
 TGAAGTAGAAGAAGTCTCACAACAGCCATTTGGAAAAAAGAAGTGTGATGAGTATATCCAGGTACA
 ACCTCCTTAGGCATGTCTGTTTTAACCTAAGCAACGCCATTATGGGCAGTGGGATTTGGGACTCGCT
 TTGCCCTGGCAAACACTGGAATCCTACTTTTTCTGGTACTTTTGACTTCAGTGACATTGTGTCTATATA
 TTCAATAAACCTCCTATTGATCTGTTCAAAGAAACAGGCTGCATGGTGTATGAAAAGCTGGGGAAACA
 GTCTTTGGCACCACAGGAAGTTCGTAATCTTTGGAGCCACCTCTACAGAACTGGAGCAATGCTGA
 GCTACCTTTCATCGTAAAAATGAACTACCTCTGCCATAAAGTTTCTAATGGGAAAGGAAGACATT
 TTCAGCCTGGTACGTGGATGGCCGCTTCTGGTGGTATAGTTACCTTTGGCATAATTCTCCCTCTGTGT
 CTCTTGAAGAAGTTAGGGTATCTTGGCTATACTAGTGGATTTTCCTTGAGCTGTATGGTTTTTTTCTAA
 TTGTGGTTATTTACAAGAAATTTCAAATCCCTGCATTGTTCCAGAGCTAAATCAACAATAAGTGCTAA
 TTCAACAAATGCTGACACGTGTACGCCAAAATATGTTACCTTCAATTCAAAGACCGTGTATGCTTTACCC
 ACCATTGCATTTGCATTTGTTGGCACCCGTCAGTCTGCCAATTTACAGTGAGCTTAAAGACCGATCAC
 AGAAAAAATGCAGATGGTTTCAAACATCTCCTTTTTCGCCATGTTTGTATGTACTTCTTGACTGCCAT
 TTTTGGCTACTTGACATTCTATGACAACGTGCAGTCCGACCTCCTCACAAATATCAGAGTAAAGATGAC
 ATTCTCATCCTGACAGTGGCTGGCTGTCATTGTTGCTGTGATCCTCACAGTCCGGTGTATTTTTTCA
 CGGTTCTGTTTCACTTTTGAAGTGGTAAAGAAAACAAAGTTTAAATTTATGTCGTCATACCGTGGTTAC
 CTGCATACTCTTGGTTGTTATCAACTTGTGGTATCTTCATACCTCCATGAAGGATATTTTTGGAGTC
 GTAGGAGTTACATCTGCTAACATGCTTATTTTCTTCTTCTTCTTTATTTAAAAATCACAGACC
 AGGATGGAGATAAAGGAAGTCAAAGAAATTTGGGCTGCCCTTTTCTGGGCTGGGGGTGTTGTTCTCCTT
 GGTGAGCATTCCCTTGGTCATCTATGACTGGGCTGCTCATCGAGTAGTGACGAAGGCCAC

ACGCGTACGCGGCCGCTCGAGCAGAAACTCATCTCAGAAGAGGATCTGGCAGCAATGATATCCTGGATT
 ACAAGGATGACGACGATAAGGTTTAA

Protein Sequence:

>RC238512 protein sequence
 Red=Cloning site Green=Tags(s)

MMHFKSGLELTELQNMTPEDDNI SNDSNDFTEVENGQINSKFISDRESRRSLTNSHLEKKKCD EYIPGT
 TSLGMSVFNLSNAIMGSGILGLAFALANTGILLFLVLLTSVTLLSIYSINLLLICSKETGCMVYEKLG EQ
 VFGTTGKFVIFGATSLQNTGAMLSYLFIVKNELPSAIKFLMGKEETFSAWYVDGRVLLVIVTFGIILPLC
 LLKNLGYLGYTSGFSLSCMVFFLIVVIYKQFQIPCIPELNSTISANSTNADTCTPKYVTFNSKTVYALP
 TIAFAFVCHPSVLP IYSELKDRSQKMQMVSNI SFFAMFVYFLTAIFGYLTFYDNVQSDLLHKYQSKDD
 ILILTVRLAVIVAVILTVPVLFVTVRSSLFELAKKTKFNL CRHTVVTCILLVVINLLVIFIPSMKDFIGV
 VGVTSANMLIFILPSSLYLKITDQDGDKGTQRIWAALFLGLGVLFSLVSIPLVIYDWACSSSSDEGH

TRTRPLEQKLI SEEDLAANDILDYKDDDDKV

Chromatograms:

https://cdn.origene.com/chromatograms/mk6724_c05.zip

Restriction Sites:

Sgfl-MluI

Cloning Scheme:


ACCN: NM_001278388

ORF Size: 1461 bp

OTI Disclaimer: The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

OTI Annotation: This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

Reconstitution Method:

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

RefSeq: [NM_001278388.2](#)

RefSeq Size: 8385 bp

RefSeq ORF: 1464 bp

Locus ID: 81539

UniProt ID: [Q9H2H9](#)


Cytogenetics: 12q13.11

Protein Families: Transmembrane

MW: 54 kDa

Gene Summary: Amino acid transporters play essential roles in the uptake of nutrients, production of energy, chemical metabolism, detoxification, and neurotransmitter cycling. SLC38A1 is an important transporter of glutamine, an intermediate in the detoxification of ammonia and the production of urea. Glutamine serves as a precursor for the synaptic transmitter, glutamate (Gu et al., 2001 [PubMed 11325958]).[supplied by OMIM, Mar 2008]

Product images:



Circular map for RC238512