

Product datasheet for **RC238618**

VRK2 (NM_001288837) Human Tagged ORF Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	VRK2 (NM_001288837) Human Tagged ORF Clone
Tag:	Myc-DDK
Symbol:	VRK2
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)



[View online »](#)

ORF Nucleotide Sequence:

>RC238618 representing NM_001288837
 Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
 GCC**CGCATCGCC**

ATGCCACCAAAAAGAAATGAAAAATACAACTTCTATTCATTCCAGAAGCAAGTTCTGGATGATA
 TGGAAAGCAATCAGTGGTACTGGCAAGAAGATTGGCTCTGGAGATTGGATTGATATTTAGCTTT
 CCCACAAATAAACAGAGAAAGATGCAAGACATGTAGTAAAAGTGAATATCAAGAAAATGGCCCGTTA
 TTTTCAGAACTTAAATTTATCAGAGAGTTGCAAAAAAGACTGTATCAAAAAGTGGATAGAACGCAAAAC
 AACTTGATTATTTAGGAATTCCTCTGTTTTATGGATCTGGTCTGACTGAATTCAGGGAAGAAGTTACAG
 ATTTATGGTAATGGAAAGACTAGGAATAGATTTACAGAAGATCTCAGGCCAGAATGGTACCTTTAAAAAG
 TCAACTGCCTGCAATTAGGTATCCGAATGTTGGATGTACTGGAATATATACATGAAAATGAATATGTTCC
 ATGGTGTGTAAGAGCAGCAATCTACTTTGGGTTACAAAAATCCAGACCAGGTTTATCTTGCAGATTA
 TGGACTTCTACAGATATTGTCCAATGGGAACCACAAACAGTATCAGGAAAATCTAGAAAAGGCCAT
 AATGGGACAATAGAGTTTACCAGCTTGGATGCCCAAGGGAGTAGCCTTGTCCAGACGAAGTGACGTTG
 AGATCCTCGGCTACTGCATGCTGCGGTGTTGTGTGGGAACTTCCCTGGGAACAGAACCTGAAGGACCC
 TGTGGCTGTGCAGACTGCTAAAACAAATCTGTTGGACGAGCTCCCCAGTCAGTGCTTAAATGGGCTCCA
 TCTGGAAGCAGTTGCTGTGAAATAGCCCAATTTTGGTATGTGCTCATAGTTTAGCATATGATGAAAAGC
 CAACTATCAAGCCCTCAAGAAAATTTGAACCCTCATGGAATACCTTTAGGACCACTGGACTTTTCCAC
 AAAAGGACAGAGTATAAATGTCCATACTCCAAACAGTCAAAAAGTTGATTCACAAAAGGCTGCAACAAAG
 CAAGTCAACAAGGCACACAATAGGTTAATCGAAAAAAAAGTCCACAGTGAGAGAAGCGCTGAGTCTGTG
 CAACATGGAAAGTGCAGAAAGAGGAGAAACTGATTGGATTGATGAACAATGAAGCAGCTCAGGAAAGCAC
 AAGGAGAAGACAGAAAATATCAAGAGTCTCAAGAACCCTTGAATGAAGTAAACAGTTTCCACAAAAAATC
 AGCTATACACAATTTCCAAACTCATTATGAGCCTCATCAAGATTTTACCAGTCCAGATATATTCAAGA
 AGTCAAGATCTCCATCTTGGTATAAATACACTTCCACAGTCAAGCAGGGGATCACAGACTTAGAAAGTTC
 AACTGGACTTTGGCCTACAATTTCCAGTTTACTCTTAGTGAAGAGACAAACGCAGATGTTTATTATTAT
 CGCATCATACCTGCTCTTTGATGTTAGTATTTCTTGCTTTATTTTTTCTC

ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCTGGATT
 ACAAGGATGACGACGATAAGGTTTAA

Protein Sequence:

>RC238618 representing NM_001288837
 Red=Cloning site Green=Tags(s)

MPPKRNEKYKLPPIPFPEGKVLDDMEGNQWVLGKKIGSGGFGLIYLAFPTNKPEKDARHVVKVEYQENGL
 FSELKFYQRVAKKDCIKKIWERKQLDYLGIPLFYGSLTEFKGRSYRFMVMERLGIDLQKISGQNGTFKK
 STVLQLGIRMLDVLEYIHENEYVHGDVKAANLLGYKNPDQVYLADYGLSYRYCPNGNHKQYQENPRKGH
 NGTIEFTSLDAHKGVALSRRSDVEILGYCMLRWLGGKLPWEQNLKDPVAVQTAKTNLLDELQSVLKWAP
 SGSSCCEIAQFLVCAHSLAYDEKPNYQALKILNPHGIPGLDFSTKGQSINVHTPNSQKVDSSQKAATK
 QVNKAHNRLIEKKVHSERSAESCATWKVQKEEKLIGLMNNEAAQESTRRRQKYQESQEPLNEVNSFPQKI
 SYTQFPNSFYEPHQDFTSPDIFKKSRSWYKYTSTVSTGITDLESSTGLWPTISQFTLSEETNADVYYY
 RIIIPVLLMLVFLALFFL

TRTRPLEQKLISEEDLAANDILDYKDDDDKV

Chromatograms:

https://cdn.origene.com/chromatograms/mk6232_a07.zip

Restriction Sites:

Sgfl-MluI

Cloning Scheme:


ACCN: NM_001288837

ORF Size: 1524 bp

OTI Disclaimer: The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

OTI Annotation: This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

Reconstitution Method:

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

RefSeq: [NM_001288837.1](#), [NP_001275766.1](#)

RefSeq Size: 2499 bp

RefSeq ORF: 1527 bp

Locus ID: 7444

UniProt ID: [Q86Y07](#)


Cytogenetics: 2p16.1

Protein Families: Druggable Genome, Protein Kinase, Transmembrane

MW: 58.1 kDa

Gene Summary: This gene encodes a member of the vaccinia-related kinase (VRK) family of serine/threonine protein kinases. The encoded protein acts as an effector of signaling pathways that regulate apoptosis and tumor cell growth. Variants in this gene have been associated with schizophrenia. Alternative splicing results in multiple transcript variants that differ in their subcellular localization and biological activity. [provided by RefSeq, Jan 2014]

Product images:



Circular map for RC238618