

Product datasheet for **RC238632**

PRAME (NM_001291715) Human Tagged ORF Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	PRAME (NM_001291715) Human Tagged ORF Clone
Tag:	Myc-DDK
Symbol:	PRAME
Synonyms:	CT130; MAPE; OIP-4; OIP4
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)
Cell Selection:	Neomycin



[View online »](#)

ORF Nucleotide Sequence:

>RC238632 representing NM_001291715
 Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
 GCC**CGGATCGCC**

ATGGAACGAAGGCGTTTGTGGGTTCCATTCAGAGCCGATACATCAGCATGAGTGTGTGGACAAGCCAC
 GGAGACTTGTGGAGCTGGCAGGGCAGAGCCTGCTGAAGGATGAGGCCCTGGCCATTGCCGCCCTGGAGTT
 GCTGCCAGGGAGCTCTCCCGCCACTTTCATGGCAGCCTTTGACGGGAGACACAGCCAGACCCTGAAG
 GCAATGGTGCAGGCCTGGCCCTCACCTGCCTCCCTCTGGGAGTCTGATGAAGGGACAACATCTTCACC
 TGGAGACCTTCAAAGCTGTGCTTGTGACTTGTGTCTCCTTGCCAGGAGGTTGCCCCAGGAGGTG
 GAAACTTCAAGTCTGGATTTACGGAAGAACTCTCATCAGGACTTCTGGACTGTATGGTCTGAAACAGG
 GCCAGTCTGTACTATTTCCAGAGCCAGAAGCAGCTCAGCCCATGACAAAGAAGCGAAAAGTAGATGGTT
 TGAGCACAGAGGCAGAGCAGCCCTTTCATCCAGTAGAGGTCTCGTAGACCTGTTCTCAAGGAAGGTGC
 CTGTGATGAATTGTTCTCTACCTCATTGAGAAAGTGAAGCGAAAGAAAATGTACTACGCCTGTGCTGT
 AAGAAGCTGAAGATTTTTCGAATGCCCATGCAGGATATCAAGATGATCCTGAAAAATGGTGCAGCTGGACT
 CTATTGAAGATTTGGAAGTGACTTGTACCTGGAAGCTACCCACCTTGGCGAAATTTCTCCTTACCTGGG
 CCAGATGATTAATCTGCGTAGACTCCTCCTCTCCACATCCATGCATCTTCTACATTTCCCGGAGAAG
 GAAGAGCAGTATATCGCCAGTTCACCTCTCAGTTCCTCAGTCTGCAGTGCCTGCAGGCTCTATGTGG
 ACTCTTTATTTTCTTAGAGGCCCGCTGGATCAGTTGCTCAGGCACGTGATGAACCCCTTGAAACCCCT
 CTCAATAACTAAGTCCCGCTTCGGAAGGGATGTGATGCATCTGTCCAGAGTCCCAGCGTCACTCAG
 CTAAAGTCTCTGAGTCTAAGTGGGTCTGCTGACCGATGTAAGTCCCAGCCCTCAAAGCTCTGCTGG
 AGAGAGCCTTCCACCCTCCAGGACCTGGTCTTTGATGAGTGTGGGATCACGGATGATCAGCTCCTTGC
 CCTCCTGCCTTCCCTGAGCCACTGCTCCAGCTTACGACCTTAAGCTTCTACGGGAATTCATCTCCATA
 TCTGCCCTGCAGAGTCTCTGCAGCACCTCATCGGGCTGAGCAATCTGACCCACGTGCTGTATCCTGTCC
 CCCTGGAGAGTTATGAGGACATCCATGGTACCCTCCACCTGGAGAGGCTTGCCTATCTGCATGCCAGGCT
 CAGGGAGTTGCTGTGAGTTGGGGCGGCCAGCATGGTCTGGCTTAGTGCCAACCCCTGCTCCTACTGT
 GGGGACAGAACCTTCTATGACCCGGAGCCATCCTGTGCCCTGTTTCATGCCTAAT

ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT
 ACAAGGATGACGACGATAAGGTTTAA

Protein Sequence:

>RC238632 representing NM_001291715
 Red=Cloning site Green=Tags(s)


MERRRLWGSIQSRYISMSVWTSPPRRLVELAGQSLLKDEALIAALELLPRELFPPLFMAAFDGRHSQTLK
 AMVQAWPFTCLPLGVLKMGQHLHLETFKAVLDGLDVLQAQEVPRRWKLVLDLRKNSHQDFWTVWVSGNR
 ASLYSFPPEAAQPMTKRKRKVDGLSTEAQPFIPVEVLVDLFLKEGACDELFSYLIEKVKRKNVLRCC
 KKLKIFAMPMDIKMILKMQVLDSDIEDLEVTCTWKLPTLAKFSPYLGQMINLRRLLLSHIHASSYISPEK
 EEQYIAQFTSQFLSLQCLQALYVDSLFFLRGRLDQLLRHVMNPLETLSITNCRLESGDVMHLSQSPVSQ
 LSVLSLQVMTDVSPEPLQALLERASATLQDLVFDECGITDDQLLALLPSLSHCSQLTTLQSFYGNISISI
 SALQSLQHLIGLSNLTHVLYPVPLESYEDIHGTLHLERLAYLHARLRELLCELGRPSMWL SANPCPHC
 GDRTFYDPEPILCPCFMPN

TRTRPLEQKLISEEDLAANDILDYKDDDDKV

Restriction Sites:

SgfI-MluI

Plasmid Map:



ACCN: NM_001291715

ORF Size: 1527 bp

OTI Disclaimer: The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

OTI Annotation: This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

Reconstitution Method:

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

RefSeq: [NM_001291715.2](#)

RefSeq Size: 2218 bp

RefSeq ORF: 1530 bp

Locus ID: 23532

UniProt ID: [P78395](#)

Cytogenetics: 22q11.22

MW: 58.3 kDa

Gene Summary: This gene encodes an antigen that is preferentially expressed in human melanomas and that is recognized by cytolytic T lymphocytes. It is not expressed in normal tissues, except testis. The encoded protein acts as a repressor of retinoic acid receptor, and likely confers a growth advantage to cancer cells via this function. Alternative splicing results in multiple transcript variants. [provided by RefSeq, Apr 2014]