

Product datasheet for **RC238729**

ANKRD42 (NM_001300975) Human Tagged ORF Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	ANKRD42 (NM_001300975) Human Tagged ORF Clone
Tag:	Myc-DDK
Symbol:	ANKRD42
Synonyms:	PPP1R79; SARP
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)
Cell Selection:	Neomycin



[View online »](#)

ORF Nucleotide Sequence:

>RC238729 representing NM_001300975
 Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
 GCC**CGGATCGCC**

ATGCCCGGGTGGCCAATTCAGGCCCTCCACTTCTCTAGGAGACTGCAAACCCCTGTTCCAGGAAGA
 AGGTGCATTTTGGCAGCATACATGATGCAGTACGAGCTGGAGATGTAAAGCAGCTTTCAGAAATAGTGGT
 ACGTGGAGCCAGCATTAAATGAACCTTGATGTTCTCCATAAGTTTACCCCTTTACATTGGGCAGCACATTCT
 GGAAGTTTGGAGTGTCTTCATTGGCTGCTCTGGCATGGAGCTGATATCACACACGTAACAACGAGAGGTT
 GGACAGCATCTCACATAGCTGCAATCAGGGGTGAGGATGCTTGTGTACAGGCTCTATAATGAATGGAGC
 AAATCTGACAGCCAGGATGACCGGGATGCACCTCTTACATCTTGCTGCAACTCATGGACATTCTTTC
 ACTTTACAATAATGCTCCGAAGTGGAGTGGATCCAGTGTGACTGATAAGAGAGAATGGAGACCTGTGC
 ATATGCAGCTTTTCATGGGCGGCTTGGCTGCTTGCAACTTCTTGTTAAATGGGGTTGTAGCATAGAAGA
 TGTGGACTACAATGGAACCTTCCAGTTCACTTAGCAGCCATGGAAGGCCACCTTACTGTTTCAAATTC
 CTAGTCAGTAGAATGAGCAGTGCAGCAGGATTTTAAAAGCTTTCAATGATAATGGAGAAAATGACTGG
 ATTTGGCCAGAGGTTCTTCAAGCAGAACATTTTACAGTTTATCCAGGGGCTGAGTATGAAGGAAAAGA
 CCTAGAGGATCAGGAACTTTAGCATTTCCAGGTCATGTGGCTGCCTTTAAGGGTGATTTGGGGATGCTT
 AAGAAATAGTGGAGATGGAGTAATCAATATTAATGAGCGTGTGATAATGGATCAACTCCTATGCATA
 AAGCTGCTGGACAAGGCCACATAGAGTGTGGCAGTGGTAAATTAATGAGGAGCAGACAGTAATATTAC
 CAACAAAGCAGGGGAGAGACCCAGTGTGTGGCAAAGAGGTTTGGCCATTTGGCAGCAGTGAAGCTGTTA
 GAGGAGCTACAGAAATATGATATAGATGACGAAAATGAAATTGATGAAAATGATGTGAAAATTTTATAA
 GACATGGTGTGGAGGAAGCAGTGTGCAAGGATGATTTATGTCTGAGTGAATGAAAACAGATGC
 CAGAATGAGAGCTTACAAGAAAATGTAGAATTGAGACACTCCTGGAAAATGCGGAGAGCAACTATAAA
 CACTTGGGAGGCATAACAGAAGAAGATTTAAAGCAGAAGAAGAAGACAGCTTGAAGTCTGAAAAGACATCA
 AAGAAGTGCAGGGCCAGCTGGAGTATGAACGACTACGTAGAGAAAAATTAGAATGTCAGCTTGTGAATA
 TCGAGCAGAAGTTGATCAACTCAGGAAACACTGGAAAAAATCAAGTCCCAAACCTTTGTGGCTATGGTT
 GGTGCTCTTTTAAATACATTTATTTTCTCAAGAAGTATATAACAAGAGTGGCCAAGAGTACAAGCTTTGG
 GCTGGGCTCTTTGGACTTGAATCCTGGCTATAGCTTTTT

ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT
 ACAAGGATGACGACGATAAGGTTTAA

Protein Sequence:

>RC238729 representing NM_001300975
 Red=Cloning site Green=Tags(s)

MPGVANSGPSTSSRETANPCSRKKVHFGSIHDAVRAGDVKQLSEIVVRGASINELDVLHKFTPLHWAHSG
 GSLECLHWLLWHGADITHVTTRGWTASHIAAIRGQDACVQALIMNGANLTAQDDRGTPLHLAATHGHSF
 TLQIMLRSGVDPSTDKREWRPVHYAAFHGRLGCLQLLVKWCSEIEDVDYNGNLPVHLAAMEGHLHCFKF
 LVSRMSSATQVLKAFNDNGENVLDLAQRFFKQNILQFIQGAIEYEGKDLEDQETLAFPGHVAAFKGLGML
 KKLVEDGVININERADNGSTPMHKAAGQGHI ECLQWLKMGADSNITNKAGERPSDVAKRFAHLAAVKLL
 EELQKYDIDDENEIDENDVKYFIRHGVGSTDAKDDLCLSDLDKTDARMRAYKKIVELRHLLIEAESNYK
 HLGGITTEEDLKQKKEQLESEKTIKELQGQLEYERLRREKLEQLDEYRAEVDQLRETLEKIQVPNFVAMV
 GVLFNTEIFLKKYIQEWPRVQALGWGSLDLNPGYSLF

TRTRPLEQKLISEEDLAANDILDYKDDDDKV


Restriction Sites:

Sgfl-MluI

Cloning Scheme:



Plasmid Map:



ACCN:	NM_001300975
ORF Size:	1581 bp
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
Reconstitution Method:	<ol style="list-style-type: none">1. Centrifuge at 5,000xg for 5min.2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.3. Close the tube and incubate for 10 minutes at room temperature.4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
RefSeq:	NM_001300975.2
RefSeq Size:	3283 bp
RefSeq ORF:	1584 bp
Locus ID:	338699
UniProt ID:	Q8N9B4
Cytogenetics:	11q14.1
MW:	59.6 kDa