

## Product datasheet for **RC238842**

### **C17orf80 (NM\_001288771) Human Tagged ORF Clone**

#### **Product data:**

Product Type:	Expression Plasmids
Product Name:	C17orf80 (NM_001288771) Human Tagged ORF Clone
Tag:	Myc-DDK
Symbol:	C17orf80
Synonyms:	HLC-8; MIG3; SPEP1
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)



[View online »](#)

**ORF Nucleotide Sequence:**

>RC238842 representing NM\_001288771  
 Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC  
 GCC**CGGATCGCC**

ATGAGTGATAATCCACCCAGAATGGAAGTGTGTCCTTACTGTAAGAAGCCATTTAAACGATTAATAATCCC  
 ACTTGCCATACTGTAAGATGATAGGACCAACCATACCTACTGATCAAAAAGTTTATCAGTCCAAGCCAGC  
 TACACTCCCACGTGCTAAAAAGATGAAAGGACCAATCAAAGATTTAATTAAGCTAAAGGAAAAGAGTTA  
 GAGACAGAGAATGAAGAAAGAAATTCTAAGTTGGTGGTGGACAAACCAGAACAGACAGTGAAGACCTTTC  
 CACTGCCAGCTGTTGGTTTGGAAAGAGCAGCTACTACAAAGGCAGATAAAGACATCAAGAATCCAATCCA  
 ACCATCCTTCAAATGTTAAAAATACTAAACCAATGACTACTTCCAAGAAGAAACCAAGGCTCAGTTT  
 TACGCATCAGAGAAAACCTCTCTAAAAGAGAACTTGCCAAAGATTTGCCTAAATCAGGAGAAAAGTCGAT  
 GTAATCCTTCAAGCTGGAGCGTCTTTACTGGTTGGCTCAATAGAACCTTCTTTGTCAAATCAAGATAG  
 AAAATATTCTCAACTTACCTAATGATGTACAACTACCTCTGGTGATCTCAAATGGACAAAATTGAT  
 CCCCAAAGACAGGAACCTCTAGTAAAATTAAGTAGATGTGCTACTGGTGATTGTCATATTTCTCAAAGA  
 ATGTCAGTGATGGGGTAAAAGGTAAGAACATTATTAAGCAATGAGAGAGATTCCAAAGGCAGGGATCA  
 CCTCTCAGGAGTCCCTACTGATGTTACAGTTACTGAGACTCCAGAAAAGAACACAGAATCCCTCATTTTA  
 AGCCTTAAAATGAGCTCATTAGGTAAAATCCAAGTCATGGAGAAAACAAGGAAAGGACTTACCCTGGGAG  
 TAGAGACGTGTGGGAGCAAAGGAAATGCAGAGAAAAGTATGTCTGCAACAGAAAAGCAGGAACGGACTGT  
 CATGAGCCATGGCTGTGAGAACTTCAACACCAGGGATTCAAGTACAGGAAAGGAGTCTCAAGGGGAAAGA  
 CCACATTTAAGTTTGTTCATTCCGAGGGAGACGACTTACCAGTTTCAATCTGTATCGCAGTCAAGTAGTC  
 AAAGCTTGCCTCTAGCTACAACATTTCTCAAGAAAAGAAAAGCAGAAGCCAGAAATCATAATTGTGT  
 CCTGATGTAAGGCATTAATGGAGAGTCCCGAGGGACAGTTATCTCTGGAGCCCAAATCTGATAGTCAG  
 TTCCAAGCATCACACTGGGTGCCAGAGCCCTTATGTTGAGCCAGCGTCACACTCCTCAGAGCCCT  
 TCACCAATCATGCTGCAGCTGCTGGCAGGAAGACTCTTCGACGCTGCATGGGGCTGGAGTGGTTCCAGA  
 GCTCTATCCTGGTTACCTTGGACTAGGGGTGTTGCCAGGGAAGCCTCAGTGTGGAATGCAATGACCCAG  
 AAGCCACAACCTTATCAGTCCCAGGGGAAAGACTCTCACAAGGCTGGATCAGGTGCAACACCACCATAA  
 GGAAGAGTGGATTCGGTGGCATCAGATGCTCTCACAGGATACTTCGTCCTGTGTGTAGCTGGAGTTT  
 CAGACGCTGAAGCTGCAGCGGTGGCGTAAG

**ACGCGT**ACGCGGCCGCTCGAGCAGAAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT  
 ACAAGGATGACGACGATAAGGTTTAA

**Protein Sequence:**

>RC238842 representing NM\_001288771  
 Red=Cloning site Green=Tags(s)

MSDNPPRMEVCPYCKKPFKRLKSHLPYCKMIGPTIPTDQKVYQSKPATLPRAKMKKGPIDKLIKAKGKEL  
 ETENEERNKSLVVDKPEQTVKTFPLPAVGLERAATTKADKDIKNPIQPSFKMLKNTKPMTTFQEETKAQF  
 YASEKTSKRELAKDLKSGESRCNPSEAGASLLVGSIEPSLSNQDRKYSSTLPNDVQTTSGDLKLDKID  
 PQRQELLVKLLDVPTGDCHISPKNVSDGVKRVRTLLSNERDSKGRDHLGVPDVTVTETPEKNTESLIL  
 SLKMSSSLGKIQVMEKQEKGLTLGVETCGSKGNAEKSMSATEKQERTVMSHGCENFNTRDSVTGKESQGER  
 PHLSLFIPRETTYQFHVSQSSQSLASLATTFLQEKKAQAQNHNCVDPVKALMESPEGQLSLEPKSDSQ  
 FQASHTGCQSPLCSAQRHTPQSPFTNHAAAAGRKTLRSCMGLEWFPPELYPGYLGVLPGKPCWNAQTQ  
 KPQLISPQGERLSQGWIRCNNTTIRKSGFGGITMLFTGYFVLCSSWSFRRLKLRWRK

**TR**TRPLEQKLISEEDLAANDILDYKDDDDKV

**Restriction Sites:**

Sgfl-MluI

**Cloning Scheme:**


**ACCN:** NM\_001288771

**ORF Size:** 1641 bp

**OTI Disclaimer:** The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

**OTI Annotation:** This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.

**Components:** The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

- Reconstitution Method:**
1. Centrifuge at 5,000xg for 5min.
  2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
  3. Close the tube and incubate for 10 minutes at room temperature.
  4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
  5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

**RefSeq:** [NM\\_001288771.2](#)

**RefSeq Size:** 2168 bp

**RefSeq ORF:** 1644 bp


**Locus ID:** 55028

**Cytogenetics:** 17q25.1

Protein Families: Transmembrane

MW: 61.3 kDa

**Product images:**



Circular map for RC238842