

Product datasheet for **RC238886**

SLC39A12 (NM_001282734) Human Tagged ORF Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	SLC39A12 (NM_001282734) Human Tagged ORF Clone
Tag:	Myc-DDK
Symbol:	SLC39A12
Synonyms:	bA570F3.1; LZT-Hs8; ZIP-12
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)
Cell Selection:	Neomycin



[View online »](#)

ORF Nucleotide Sequence:

>RC238886 representing NM_001282734
 Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
 GCC**CGCATCGCC**

ATGAGTAATAAAGAGTATAAATTTTACCTACACAGCCTACTGAGCCTCAGGCAGGATGAAGATTCCTCTT
 TCCTTTCACAGAATGAGACAGAAGATATCTTGGCTTTCACCAGGCAGTACTTTGACACTTCTCAAAGCCA
 GTGTATGGAAACAAAACGCTGCAGAAAAAATCTGGAATAGTGAGCAGTGAAGGTGCTAATGAAAGTACG
 CTTCTCAGTTGGCAGCCATGATCATTACTTTGTCCCTCCAGGGTGTCTGGGACAAGGAACTTGC
 CTTCCCCAGACTACTTTCACAGAATATTTTCAGTTTCTGAATCGTACGAATACCCTCCGCCTATCAGA
 ACTAGACCAACTCCTCAACTCTCTGGACCAGAAGTACTTGTATCAAAAATGAGAAAATCCATCAATTT
 CAAAGGAAACAAAACAACATAATAACCCATGATCAGGACTATTCTAATTTCTCTTCATCCATGGAAAAAG
 AGTCTGAGGATGGTCCAGTTTCTGGGATCAGACCTGCTTCTCTGCTAGGCAGCTGGTGGAGATATTTCT
 ACAGAAGGGCCTCTCACTATTTCTAAGGAGGACTTTAAGCAAAATGAGTCCAGGGATCATCCAGCAGCTC
 CTAGCTGCTCCTGCCACTTACCCAAGGACCAACAAGCAAAGCTGCCACCTACCCTCTGGAGAAAATACG
 GCTACAGCACGGTGGCTGTCACCCTTCTCACACTGGGCTCCATGCTGGGGACAGCGCTGGTCTTTTCCA
 TAGCTGTGAGGAGAACTACAGGCTTATCTTACAGCTGTTTGTGGGCTTGGCCGTGGGACACTGTCTGGG
 GACGCTGTCTCCACCTTATCCCTCAGGTTCTTGGTTTACATAAGCAGGAAGCCCCAGAATTTGGGCATT
 TCCATGAAAGCAAAGGTCAATTTGGAAACTGATGGGATTAATTGGAGGCATCCATGGATTTTCTTGAT
 AGAAAAATGTTTTATTCTTGTATCACCAAAATGACAAGCAGGGCCTGTATTGGTTAATGGGCACGTG
 GGTCAATCCCACCATCTGCACTCAACTCTGAATTAAGTGACCAGGCAGGCAGAGGCAAAATCTGCTTCAA
 CTATCCAGTTGAAAAGCCAGAAGATTCACAGGCAGCTGAAATGCCTATAGGCAGTATGACAGCTCCAA
 CAGAAAAATGTAAGCCATTAGCTTGTAGCAATCATGATTCTGGTTGGGGACAGCCTGCATAATTTGCA
 GATGGCCTAGCCATAGGAGCAGCCTTCTCATCATCCGAGTCAGGAGTGACCACTACGATTGCTATCT
 TGTGTCATGAAATCCCACATGAAATGGGAGACTTTGCCGTGCTCTAAGCTCTGGACTTTCTATGAAGAC
 TGCCATCTGATGAATTTATAAGCTCCCTAAGCTGCTTTCATGGGATTATACATTGGCCTTTCCGTGTCA
 GCTGATCCATGTGTTCAAGACTGGATCTTACAGTCACTGCTGGGATGTTCTTATATTTATCCTTGGTTG
 AAATGCTTCTGAAATGACTCATGTTCAAACACAACGACCCTGGATGATGTTTCTCTGCAAAACTTTGG
 ATTGATCTAGGTTGGCTTCTCTCTGCTCTTGGCTATATATGAGCAAAATATTTAAATA

ACGCGTACGCGGCCGCTCGAGCAGAAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT
 ACAAGGATGACGACGATAAGGTTTAA

Protein Sequence:

>RC238886 representing NM_001282734
 Red=Cloning site Green=Tags(s)


MSNKEYK FYLHSLLSLRQDESSFLSQNETEDILAFTRQYFDTSQSQCMTKTLQKKSGIVSSEGANESE
 LPQLAAMIITLSLQGVCLGQGNLPSDYFTEYIFSSLNRTNLTLSLQDQLLNTLWTRSTCIKNEKIHQF
 QRKQNNIIITHDQDYSNFSMSMEKESEDGPVSWDQTCFSARQLVEIFLQKGLSLISKEDFKQMSPGIIQQL
 LSCSCHLPKDQQAQLPPTTLEKYGYSTVAVTLLTGLSMLGTALVLFHSCSEENYRLILQLFVGLAVGTL
 SG DALLHLIPQVLGLHKQEAPEFGHFHESKGIWKL MGLIGGIHGFLLIEKCFILLVSPNDKQGLSLVNGHV
 GHSHHLALNSELSDQAGRKKSASTIQLKSPEDSQAAMPISMTASNRKCKAISLLAIMILVGDLSLHNFA
 DGLAIGAASFSSSESGVTTTIAILCHEIPHEMGDFAVLLSSGLSMKTAIILMNFISLTAFMGLYIGLSVS
 ADPCVQDWIFTVTAGMFLYLSLVEMLPEMTHVQTQRPWMMFLLQNFGLILGWLSLLLLAIYEQNIKI

TRTRPLEQKLISEEDLAANDILDYKDDDDKV


Restriction Sites:

Sgfl-MluI

Cloning Scheme:



Plasmid Map:



ACCN:	NM_001282734
ORF Size:	1671 bp
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
Reconstitution Method:	<ol style="list-style-type: none">1. Centrifuge at 5,000xg for 5min.2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.3. Close the tube and incubate for 10 minutes at room temperature.4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
RefSeq:	NM_001282734.1 , NP_001269663.1
RefSeq Size:	2461 bp
RefSeq ORF:	1674 bp
Locus ID:	221074
UniProt ID:	Q504Y0
Cytogenetics:	10p12.33
Protein Families:	Transmembrane
MW:	62.1 kDa
Gene Summary:	Zinc is an essential cofactor for hundreds of enzymes. It is involved in protein, nucleic acid, carbohydrate, and lipid metabolism, as well as in the control of gene transcription, growth, development, and differentiation. SLC39A12 belongs to a subfamily of proteins that show structural characteristics of zinc transporters (Taylor and Nicholson, 2003 [PubMed 12659941]).[supplied by OMIM, Aug 2008]