

Product datasheet for **RC238921**

CLPTM1 (NM_001282176) Human Tagged ORF Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	CLPTM1 (NM_001282176) Human Tagged ORF Clone
Tag:	Myc-DDK
Symbol:	CLPTM1
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)
Cell Selection:	Neomycin



[View online »](#)

ORF Nucleotide Sequence:

>RC238921 representing NM_001282176
 Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
 GCC**CGATCGCC**

ATGAACCTGCATGTGTACATCTCAGAGCAGGACTTTACAGACTTCAACGCCACGTCGGCACTTTCT
 GGAACAGCACGATCTTGTGTATGGCGACTGGACTAGCGGCGAGAAGCTCAGACGGCTGCTACGAGCACTT
 TGCTGAGCTCGATATCCACAGAGCGTCCAGCAGAACGGCTCCATCTACATCCACGTTTACTTACCAAG
 AGTGGCTTCCACCCAGACCCCGGCAGAAAGGCCCTGTACCGCCGGCTTGCACAGTCCACATGTCCCGGA
 TGATCAACAAATACAAGCGCAGACGATTTAGAAAACCAAGAAGCTGCTGACAGGAGAGACAGAAGCGGA
 CCCAGAAATGATCAAGAGGGCTGAGGACTATGGGCTGTGGAGGTGATCTCCATTGGCACCCCAACATC
 ACCATCAACATCGTGGACGACCACGCCGTGGTGAAGGGCAGTGTCCCCCTCCCCTGGATCAATATG
 TGAAGTTCGACCGGTGAGCGGTGACTACTATCCCATCATCTACTTCAATGACTACTGGAACCTGCAGAA
 GGACTACTACCCATCAACGAGAGCCTGGCCAGCCTGCCGCTCCGCGTCTCCTTCTGCCACTCTCGCTT
 TGGCGCTGGCAGCTCTATGCTGCCAGAGCACCAAGTCGCCCTGGAACTTCTGGGTGATGAGTTGACG
 AGCAGTCAGATGAGGAGCAGGACTCGGTGAAGGTGGCCCTGCTGGAGACCAACCCCTACCTGCTGGCGCT
 CACCATCATCGTGTCTATCGTTCACAGTGTCTTCGAGTTCCTGGCCTTCAAGAATGATATCCAGTTCTGG
 AACAGCCGGCAGTCCCTGGAGGGCTGTCCGTGCGCTCCGTCTTCTTCGGCGTTTTCCAGTCATTCTGG
 TCCTCCTACATCCTGGACAACGAGACCAACTTCGTGGTCCAGGTCAGCGTCTTATTGGGGTCTCAT
 CGACCTCTGGAAGATACCAAGGTCATGGACGTCGGCTGGACCGAGAGCACAGGGTGGCAGGAATCTTC
 CCCCCTATCCTTCAAGGACAAGTCCACGTATATCGAGTCTCGACCAAGTGTATGATGATATGGCAT
 TCCGGTACCTGTCTGGATCCTCTCCCGCTCCTGGGCTGCTATGCCGCTACAGTCTTCTGTACCTGGA
 GCACAAGGGCTGGTACTCCTGGGTGCTCAGCATGCTCTACGGCTTCTGCTGACCTTCGGCTTATCACC
 ATGACGCCCCAGCTTTCATCAACTACAAGCTCAAGTCTGTGGCCACCTTCCCTGGCGCATGCTCACCT
 ACAAGGCCCTCAACACATTATCGACGACCTGTTGCGCTTGTGATCAAGATGCCGTTATGTACCGGAT
 CGGCTGCTGCGGACGATGTGGTTTTCTTCTACCTACCAACGGTGGATCTACCGCTCGACCCC
 ACCCGAGTCAACGAGTTTGGCATGAGTGGAGAAGACCCACAGCTGCCGCCCCGTGGCCGAGGTTCCCA
 CAGCAGCAGGGGCCCTCAGCCACACCTGCACCCACCAGCACCACCAGGGAGGAGGCTCCAC
 GTCCCTGCCACCAAGCCACCCAGGGGCCAGCTCTGCCAGCGAGCCCAAGAACCCCTCAAAGCCA
 GCAGAGGACAAGAAAAGGAT

AG**CGGACCG**ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCC
 TGGATTACAAGGATGACGACGATAAGGTTTAA

Protein Sequence:

>RC238921 representing NM_001282176
 Red=Cloning site Green=Tags(s)


MNLHVYISEHEHFTDFNATSALFWEQHDLVYGDWTSGENSDGCEYHFAELDIPQSVQQNGSIYIHVYFTK
 SGFHPDPRQKALYRRLATVHMSRMINKYKRRRFQKTKNLLTGTEADPEMIKRAEDYGPVEVISHWHPNI
 TINIVDDHTPWVKGSPVPLDQYVKFDVAVSGDYPIIYFNDYWNLQKDYYPINESLASLPLRVSFCLSL
 WRWQLYAAQSTKSPWNFLGDELVEQSDEEQDSVKVALLETNPYLLALTIIVSIVHSVFEFLAFKNDIQFW
 NSRQSLEGLSVRSVFFGVFQSFVLLYILDNETNFVQVSVFIVGLIDLWKITKMDVRLDREHRVAGIF
 PRLSFKDKSTYIESSTKYDDMAFRYLSWILFPLLGCYAVYSLLYLEHKGWYSWVLSMLYGFLLTFGFIT
 MTPQLFINYKLLKSAHLRMLTYKALNTFIDDLFAFVIKMPVMYRIGCLRDDVFFIYLYQRWIYRVPD
 TRVNEFGMSGEDPTAAAPVAEVPAAAGALTPTPAPTTTTATREEASTSLPTKPTQGASSASEPQEAPPKP
 AEDKKKD

SGP**TRTRPLEQKLI**SEEDLAANDILDYKDDDDKV


Restriction Sites:

Sgfl-RsrII

Cloning Scheme:



Plasmid Map:



ACCN:

NM_001282176

ORF Size:	1701 bp
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
Reconstitution Method:	<ol style="list-style-type: none">1. Centrifuge at 5,000xg for 5min.2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.3. Close the tube and incubate for 10 minutes at room temperature.4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
RefSeq:	NM_001282176.1 , NP_001269105.1
RefSeq Size:	2522 bp
RefSeq ORF:	1704 bp
Locus ID:	1209
UniProt ID:	O96005
Cytogenetics:	19q13.32
Protein Families:	Transmembrane
MW:	66 kDa
Gene Summary:	May play a role in T-cell development.[UniProtKB/Swiss-Prot Function]