

## Product datasheet for **RC238980**

### SCFD1 (NM\_001283032) Human Tagged ORF Clone

#### Product data:

Product Type:	Expression Plasmids
Product Name:	SCFD1 (NM_001283032) Human Tagged ORF Clone
Tag:	Myc-DDK
Symbol:	SCFD1
Synonyms:	C14orf163; RA410; SLY1; SLY1P; STXBP1L2
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)
Cell Selection:	Neomycin



[View online »](#)

**ORF Nucleotide Sequence:**

>RC238980 representing NM\_001283032  
 Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC  
 GCC**CGCATCGCC**

ATGTGCCTCATATTAACACAGCACAGGAGAACCAGTATGGAAGGCTTTTACACTCTGATCGAGATCCTA  
 TTCAGATGTTCCCTGCAGTATACTTTGTAATGCCAACTGAAGAAAATATTGACAGAATGTGCCAGGATCT  
 TCGAAATCAACTATATGAATCATATTATTTAAATTTTATTCTGCTATTTCAAGAAGTAAACTGGAAGAT  
 ATTGCAAATGCAGCGTTAGCAGCTAGTGCAGTAACACAAGTAGCCAAGGTTTTTGACCAATATCTCAATT  
 TTATTACTTTGGAAGATGATATGTTTGTATTATGTAATCAAAAATAAGGAGCTTGTTCATATCGTGCCAT  
 TAACAGGCCAGATATCACAGACACGAAATGGAACTGTTATGGACACTATAGTTGACAGCCTCTTCTGC  
 TTTTTTGTACTCTGGTGCTGTTCTATAATCAGATGTTCAAGAGGAACAGCAGCAGAAATGGTAGCAG  
 TGAAACTAGACAAGAACTTCGAGAAAACTAAGAGATGCAAGAAACAGTCTTTTACAGGTGATACACT  
 TGGAGCTGGCCAATTCAGCTCCAGAGGCCCTTATTAGTCCTTGTGACAGAAACATAGATTTGGCAACT  
 CCTTTACATCATACTTGGACATATCAAGCATTGGTGCACGATGTACTGGATTTCCATTTAAACAGGGTTA  
 ATTTGGAAGAATCTTCAGGAGTGGAAAACCTCCAGCTGGTGCTAGACCAAGAGAAAAACAAGAAGTC  
 TTATGATTTAACTCCGGTTGATAAATTTTGGCAAAAACATAAAGGAAGTCCATTTCCAGAAGTTGCAGAA  
 TCAGTTACAGAAAGTGAATCTTACAGAGCACAGGAAGATGAGGTCAAACGACTTAAAGCATTATGAG  
 GACTAGAAGGGGAAGATGAAGGAGCCATAAGTATGCTTTCTGACAATACCGCTAAGCTAACATCAGCTGT  
 TAGTCTTTGCCAGAACTCCTTGAGAAAAAAGACTTATTGATCTCCATACAAAATGTTGCCACTGCTGTT  
 TTAGAACATATAAAGGCAAGAAAATTGGATGTATATTTGAAATATGAAGAAAAATAATGAGCAAACTA  
 CTCTGGATAAACTCTTCTAGATATAATATCAGACCCTGATGCAGGAACCTCAGAAGATAAAATGAGGTT  
 GTTTCTTATCTATTATATAAGCACACAGCAAGCACCTTCTGAGGCTGATTTGGAGCAATATAAAAAAGCT  
 TTAAGTATGACAGGATGCAACCTTAATCCTTTACAATATATCAAAACAGTGGAAAGGCTTTTACCAAGTGG  
 CCTCAGCTCCGGCCAGCTATGGCAGCACTACCACTAAACCAATGGGTCTTTTATCACGAGTCATGAATAC  
 AGGATCACAGTTTGTGATGGAAGGAGTGAAGAACCCTGGTTTTGAAACAGCAAAATCTACCTGTTACTCGT  
 ATTTTGGACAATCTTATGGAGATGAAGTCAAACCCGAACTGATGACTATAGATATTTTATGCCAAAA  
 TGCTGCGGGCAATGACAGCTCAGTTCCAGAAAATAAAATCCATTCCAAGAGGCCATTGTTTTTGTGGT  
 GGGAGGAGCAACTACATTGAATATCAGAACTTGTGACTACATAAAGGGGAAACAAGGCAAAACACATT  
 TTATATGGCTGCAGTGAAGCTTTTAAATGCTACACAGTTCATAAAACAGTTGTCACAACCTGGACAAAAG

**ACGCGT**ACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT  
 ACAAGGATGACGACGATAAGGTTTAA

**Protein Sequence:**

>RC238980 representing NM\_001283032  
 Red=Cloning site Green=Tags(s)


MCLILKTAQENQYGRLLHSDRDPIDVPAVYFVMPTEENIDRMCQDLRNQLYESYYLNFISAISSKLED  
 IANAALAASAVTQYAKVFDQYLNFTLEDDMFVLCNQNKELVSYRAINRPDITDTEMETVMDTIVDSLFC  
 FFVTLGAVPIIIRSRGTAEMVAVKLDKLLRENLRDARNSLFTGDTLGAQSFQRPPLLVLVDRNIDLAT  
 PLHHTWYQALVHDVLDVFLNLRVNLEESSGVENSPAGARPKRKNKSYDLTPVDKFWQKHKGSPFPEVAE  
 SVQQELESYRAQEDEVKRLKSI MGLEGEDEGAISMLSDNTAKLTSAVSSLPELLEKKRLIDLHTNVATAV  
 LEHIKARKLDVYFEYEEKIMSKTTLDKSLLDIISDPDAGTPEDKMRLFLIYYISTQQAPSEADLEQYKKA  
 LTDAGCNLNPLQYIKQWKAFTKMASAPASYGSTTTKPMGLLSRVMTGSQFVMEGVKNLVLKQQNLPVTR  
 ILDNLMEMKSNPETDDYRYFDPKMLRGNDSSVPRNKNPFQEAIVFVGGGNYIEYQNLVDYIKGKQKHI  
 LYGCSELFNATQFIKQLSQLGQK

**TRTRPLEQKLI**SEEDLAANDILDYKDDDDKV


**Restriction Sites:**

SgfI-MluI

Cloning Scheme:



Plasmid Map:



ACCN: NM\_001283032  
 ORF Size: 1749 bp

<b>OTI Disclaimer:</b>	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. <a href="#">More info</a>
<b>OTI Annotation:</b>	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
<b>Components:</b>	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
<b>Reconstitution Method:</b>	<ol style="list-style-type: none"><li>1. Centrifuge at 5,000xg for 5min.</li><li>2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.</li><li>3. Close the tube and incubate for 10 minutes at room temperature.</li><li>4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.</li><li>5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.</li></ol>
<b>RefSeq:</b>	<a href="#">NM_001283032.1</a> , <a href="#">NP_001269961.1</a>
<b>RefSeq Size:</b>	2144 bp
<b>RefSeq ORF:</b>	1752 bp
<b>Locus ID:</b>	23256
<b>UniProt ID:</b>	<a href="#">Q8WVM8</a>
<b>Cytogenetics:</b>	14q12
<b>MW:</b>	66.5 kDa
<b>Gene Summary:</b>	Plays a role in SNARE-pin assembly and Golgi-to-ER retrograde transport via its interaction with COG4. Involved in vesicular transport between the endoplasmic reticulum and the Golgi (By similarity).[UniProtKB/Swiss-Prot Function]