

## Product datasheet for **RC238984**

### MMP2 (NM\_001302510) Human Tagged ORF Clone

#### Product data:

Product Type:	Expression Plasmids
Product Name:	MMP2 (NM_001302510) Human Tagged ORF Clone
Tag:	Myc-DDK
Symbol:	MMP2
Synonyms:	CLG4; CLG4A; MMP-2; MMP-II; MONA; TBE-1
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)
Cell Selection:	Neomycin



[View online »](#)

**ORF Nucleotide Sequence:**

>RC238984 representing NM\_001302510  
 Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC  
 GCC**CGGATCGCC**

ATGCAGAAGTCTTTGGACTGCCCCAGACAGGTGATCTTGACCAGAATACCATCGAGACCATGCGGAAGC  
 CACGCTGCGGCAACCCAGATGTGGCCAAC TACAAC TTTCCCTCGCAAGCCAAGTGGGACAAGAACCA  
 GATCACATA CAGGATCATTGGCTACACACCTGATCTGGACCCAGAGACAGTGGATGATGCCTTTGCTCGT  
 GCCTTCCAAGTCTGGAGCGATGTGACCCCACTGCGGTTTTCTCGAATCCATGATGGAGAGGCAGACATCA  
 TGATCAACTTTGGCCGCTGGGAGCATGGCGATGGATACCCCTTTGACGGTAAGGACGGACTCCTGGCTCA  
 TGCCTTCGCCCCAGGCACTGGTGTGGGGGAGACTCCCATTTTGTGACGATGAGCTATGGACCTTGGGA  
 GAAGGCCAAGTGGTCCGTGTGAAGTATGGGAACGCCGATGGGGAGTACTGCAAGTCCCCTTCTGTTC  
 ATGGCAAGGAGTACAACAGCTGCACTGATACCGGCCGACGCGATGGCTTCTCTGGTCTCCACCACCTA  
 CAAC TTTGAGAAGGATGGCAAGTACGGCTTCTGTCCCCATGAAGCCCTGTTACCATGGGCGGCAACGCT  
 GAAGGACAGCCCTGCAAGTTCCATTCCGCTTCCAGGGCACATCCATGACAGCTGCACCACTGAGGGCC  
 GCACGGATGGCTACCGCTGGTGGCCACCACTGAGGACTACGACCCGACAAAGAATATGGCTTCTGCC  
 TGAGACCGCCATGTTCCACTGTTGGTGGAACTCAGAAGGTGCCCCCTGTGTCTTCCCCTTCACTTCTCTG  
 GGCAACAAATATGAGAGCTGCACCAGCGCCGGCCGAGTACGCGAAAGATGTGGTGTGCGACCACAGCCA  
 ACTACGATGATGACCGCAAGTGGGGCTTCTGCCCTGACCAAGGGTACAGCCTGTTCTCTGTCGAGCCCA  
 CGAGTTTGGCCACGCCATGGGGCTGGAGCACTCCCAAGACCCTGGGGCCCTGATGGCACCATTTACACC  
 TACACCAAGAACTTCCGTCTGTCCAGGATGACATCAAGGGCATTGAGGAGCTCTATGGGCTCTCCCTG  
 ACATTGACCTTGGCACCGGCCCAACCCACGCTGGGCCCTGTCACTCCTGAGATCTGCAACACAGGACT  
 TGTATTTGATGGCATCGCTCAGATCCGTGGTGGATCTTCTTCTCAAGGACCGGTTCAATTTGGCGGACT  
 GTGACGCCACGTGACAAGCCCATGGGGCCCTGTGGTGGCCACATTCTGGCCTGAGCTCCCGGAAAAGA  
 TTGATGCGGTATACGAGGCCCCACAGGAGGAGAAGGCTGTGTTCTTTCAGGGAATGAATACTGGATCTA  
 CTCAGCCAGCACCTGGAGCGAGGGTACCCCAAGCCACTGACCAGCCTGGGACTGCCCCCTGATGTCCAG  
 CGAGTGGATGCCGCTTTAACTGGAGCAAAAACAAGAAGACATACATCTTTGCTGGAGACAAATTCTGGA  
 GATAAATGAGGTGAAGAAGAAAATGGATCCTGGCTTCCCAAGCTCATCGCAGATGCCTGGAATGCCAT  
 CCCCATAACCTGGATGCCGCTGTGGACCTGCAGGGCGGGGTCACAGCTACTTCTCAAGGGTGCCTAT  
 TACCTGAAGCTGGAGAACCAAGTCTGAAGAGCGTGAAGTTTGAAGCATCAAATCCGACTGGCTAGGCT  
 GC

**ACGCGT**ACGCGGCCGCTCGAGCAGAAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT  
 ACAAGGATGACGACGATAAGGTTTAA

**Protein Sequence:**

>RC238984 representing NM\_001302510  
 Red=Cloning site Green=Tags(s)

MQKFFGLPQTGDLQNTIETMRKPRCGNPDVANYNFFPRPKPKWDKNQITYRIIGYTPDLDPETVDDAFAR  
 AFQVWSDVTPLRFSRIHDGEADIMINFGRWEHGDGYFPDGDGLLAHAFAPGTG VGGDSHFDDDELWTLG  
 EGQVVRVKYGNADGEYCKFPFLFNGKEYNSCTDTGRSDGFLWCSTTYNFEKDGYGFCPHEALFTMGGNA  
 EGQPCKFPFRFQGSYDSCTTEGRTDGYRWCGTTEDYDRDKKYGFCPETAMSTVGGNSEGAPCVFPFTFL  
 GNKYESCTSAGRSDGKMWCAATANYDDDRKWGFCPDQYSLFLVAAHEFGHAMGLEHSQDPGALMAPIYT  
 YTKNFRLSQDDIKGIQELYGASPDIDLGTGPTPTLGPVTPEICKQDIVFDGIAQIRGEIFFKDRFIWRT  
 VTPRDKPMGPLL VATFWPELPEKIDAVYEAPQEEKAVFFAGNEYWIYSASTLERGYPKPLTSLGLPPDVQ  
 RVDAAFNWSKNKTYIFAGDKFWRYNEVKKMDPGFPKLIADAWNAIPDNLDAVVDLQGGHSHYFFKGAY  
 YLKLENQSLKSVKFGSIKSDLWGC

**TRTRPLEQKLI**SEEDLAANDILDYKDDDDKV

**Restriction Sites:**

SgfI-MluI



<b>OTI Disclaimer:</b>	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. <a href="#">More info</a>
<b>OTI Annotation:</b>	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
<b>Components:</b>	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
<b>Reconstitution Method:</b>	<ol style="list-style-type: none"><li>1. Centrifuge at 5,000xg for 5min.</li><li>2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.</li><li>3. Close the tube and incubate for 10 minutes at room temperature.</li><li>4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.</li><li>5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.</li></ol>
<b>RefSeq:</b>	<a href="#">NM_001302510.1</a> , <a href="#">NP_001289439.1</a>
<b>RefSeq Size:</b>	3152 bp
<b>RefSeq ORF:</b>	1755 bp
<b>Locus ID:</b>	4313
<b>UniProt ID:</b>	<a href="#">P08253</a>
<b>Cytogenetics:</b>	16q12.2
<b>Protein Families:</b>	Druggable Genome, Protease
<b>Protein Pathways:</b>	Bladder cancer, GnRH signaling pathway, Leukocyte transendothelial migration, Pathways in cancer
<b>MW:</b>	66.2 kDa

**Gene Summary:**

This gene is a member of the matrix metalloproteinase (MMP) gene family, that are zinc-dependent enzymes capable of cleaving components of the extracellular matrix and molecules involved in signal transduction. The protein encoded by this gene is a gelatinase A, type IV collagenase, that contains three fibronectin type II repeats in its catalytic site that allow binding of denatured type IV and V collagen and elastin. Unlike most MMP family members, activation of this protein can occur on the cell membrane. This enzyme can be activated extracellularly by proteases, or, intracellularly by its S-glutathiolation with no requirement for proteolytical removal of the pro-domain. This protein is thought to be involved in multiple pathways including roles in the nervous system, endometrial menstrual breakdown, regulation of vascularization, and metastasis. Mutations in this gene have been associated with Winchester syndrome and Nodulosis-Arthropathy-Osteolysis (NAO) syndrome. Alternative splicing results in multiple transcript variants encoding different isoforms. [provided by RefSeq, Oct 2014]