

Product datasheet for **RC239010**

STIP1 (NM_001282652) Human Tagged ORF Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	STIP1 (NM_001282652) Human Tagged ORF Clone
Tag:	Myc-DDK
Symbol:	STIP1
Synonyms:	HEL-S-94n; HOP; IEF-SSP-3521; P60; STI1; STI1L
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)
Cell Selection:	Neomycin



[View online »](#)

ORF Nucleotide Sequence:

>RC239010 representing NM_001282652
 Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
 GCC**CGGATCGCC**

ATGGAGTCCGGCAGCCCAATGGGAGAGGTGGAATTTCCAGAACGATCAGAACCAATGGGCGGCCAGC
 GCGGCTACGATTGGCAGTGCAAAAGACCAATCCGTGTCGAGAAGTTCGCTCCTCCTCCATTCGTGGAG
 CCTGAGATGGGTCAATGAGCTGAAGGAGAAAGGCAACAAGGCCCTGAGCGTGGGTAAACATCGATGATGCC
 TTACAGTCTACTCCGAAGCTATTAAGCTGGATCCCCACAACCACGTGCTGTACAGCAACCGTTCTGCTG
 CCTATGCCAAGAAAGGAGACTACCAGAAGGCTTATGAGGATGGCTGCAAGACTGTCGACCTAAAGCCTGA
 CTGGGGCAAGGGCTATTCACGAAAAGCAGCAGCTCTAGAGTTCTTAAACCGCTTTGAAGAAGCCAAGCGA
 ACCTATGAGGAGGGCTTAAACACGAGGCAAATAACCCTCAACTGAAAGAGGGTTTACAGAAATATGGAGG
 CCAGGTTGGCAGAGAGAAAATTCATGAACCTTTCAACATGCCTAATCTGTATCAGAAGTTGGAGAGTGA
 TCCAGGACAAGGACACTACTCAGTGATCTACCTACCGGAGCTGATAGAGCAGCTACGAAACAAGCCT
 TCTGACCTGGGCACGAACTACAAGATCCCGGATCATGACCCTCTCAGCGTCCCTTGGGGTTCGATC
 TGGGCAGTATGGATGAGGAGGAAGAGATTGCAACACCTCCACCACCACCCCTCCCAAAAAGGAGACCAA
 GCCAGAGCCAATGGAAGAAGATCTTCCAGAGAATAAGAAGCAGGCACTGAAAGAAAAAGAGCTGGGGAAC
 GATGCCTACAAGAAGAAAGACTTTGACACAGCCTTGAAGCATTACGACAAAGCCAAGGAGCTGGACCCCA
 CTAACATGACTTACATTACCAATCAAGCAGCGGTATACTTTGAAAAGGGCGACTACAATAAGTGCCGGGA
 GCTTTGTGAGAAGGCCATTGAAGTGGGGAGAGAAAACCGAGAAGACTATCGACAGATTGCCAAAGCATAT
 GCTCGAATTGGCAACTCTACTTCAAAGAAGAAAAGTACAAGGATGCCATCCATTTCTATAACAAGTCTC
 TGGCAGAGCACCGAACCCAGATGTGCTCAAGAAATGCCAGCAGCAGAGAAAATCCTGAAGGAGCAAGA
 GCGGCTGGCCTACATAAAACCCCGACCTGGCTTTGGAGGAGAAGAACAAGGCAACGAGTGTTCAGAAA
 GGGGACTATCCCAAGCCATGAAGCATTATACAGAAGCCATCAAAAGGAACCCGAAAAGATGCCAAATAT
 ACAGCAATCGAGCTGCCTGCTACACCAAACCTCTGGAGTTCAGCTGGCACTCAAGGACTGTGAGGAATG
 TATCCAGCTGGAGCCGACCTTCAAGGGTTATACACGAAAGCCGCTGCGCTGGAAGCGATGAAGGAC
 TACACCAAAGCCATGGATGTGTACCAGAAGGCGCTAGACTGGACTCCAGCTGTAAGGAGGCGGCAGACG
 GCTACCAGCGTGTATGATGGCGCAGTACAACCGCACGACAGCCCGAAGATGTGAAGCGACGAGCCAT
 GGCCGACCCTGAGGTGCAGCAGATCATGAGTGACCCAGCCATGCGCCTTATCCTGGAACAGATGCAGAAG
 GACCCCAAGGCACTCAGCGAACACTTAAAGAATCTGTAATAGCACAGAAGATCCAGAAGCTGATGGATG
 TGGGTCTGATTGCAATTCGG

ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT
 ACAAGGATGACGACGATAAGGTTTAA

Protein Sequence:

>RC239010 representing NM_001282652
 Red=Cloning site Green=Tags(s)

MESGSPMGEVEISRTIRTNGRGQRGYDWQCKRPIRVAEVRSSLHSWSLRWVNLKEKGNKALSVGNIDDA
 LQCYSEAIKLDPHNHVLYSNRSAAYAKKGDYQKAYEDGCKTVDLKPDWKGYSRKAALFLNRFEEAKR
 TYEGLKHEANNPQLKEGLQNMEARLAERKFMNPFNMPNLQKLESDPRTRTLSDPTYRELIEQLRNKP
 SDLGTLKLDQPRIMTTLVLLGVDLGSMDDEEEIATPPPPPPPKETKPEPMEEDLPENKKQALKEKELGN
 DAYKKKDFDTALKHYDKAKELDPTNMTYITNQAAVYFEKGDYKNCRELCEKAIIEVGRENREDYRQIAKAY
 ARIGNSYFKEEKYKDAIHFNKSLAEHRTPDVLKCCQQAIEKILKEQERLAYINPDALAEKNKGNCFQK
 GDYPQAMKHYTEAIKRNPKDAKLYSNRAACYTKLLEFQLALKDCEECIQLEPTFIKGYTRKAAALEAMKD
 YTKAMDVYQKALDLSSCKEAADGYQRCMMAQYNRHDSPEVDKRRAMADPEVQQIMSDPAMRLILEQMOK
 DPQALSEHLKNPVIAQKIQKLMVGLIAIR

TRTRPLEQKLISEEDLAANDILDYKDDDDKV

Restriction Sites:

SgfI-MluI


Cloning Scheme:

Cloning sites used for ORF Shutting:



* The last codon before the Stop codon of the ORF

Plasmid Map:



ACCN: NM_001282652

ORF Size: 1770 bp

OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
Reconstitution Method:	<ol style="list-style-type: none"> 1. Centrifuge at 5,000xg for 5min. 2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA. 3. Close the tube and incubate for 10 minutes at room temperature. 4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom. 5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
RefSeq:	NM_001282652.1 , NP_001269581.1
RefSeq Size:	2766 bp
RefSeq ORF:	1773 bp
Locus ID:	10963
UniProt ID:	P31948
Cytogenetics:	11q13.1
Protein Families:	Stem cell - Pluripotency
Protein Pathways:	Prion diseases
MW:	68.5 kDa
Gene Summary:	STIP1 is an adaptor protein that coordinates the functions of HSP70 (see HSPA1A; MIM 140550) and HSP90 (see HSP90AA1; MIM 140571) in protein folding. It is thought to assist in the transfer of proteins from HSP70 to HSP90 by binding both HSP90 and substrate-bound HSP70. STIP1 also stimulates the ATPase activity of HSP70 and inhibits the ATPase activity of HSP90, suggesting that it regulates both the conformations and ATPase cycles of these chaperones (Song and Masison, 2005 [PubMed 16100115]).[supplied by OMIM, Jul 2009]