

Product datasheet for **RC239017**

LOC100996634 (NM_001277339) Human Tagged ORF Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	LOC100996634 (NM_001277339) Human Tagged ORF Clone
Tag:	Myc-DDK
Symbol:	LOC100996634
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)



[View online »](#)

ORF Nucleotide Sequence:

>RC239017 representing NM_001277339
 Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
 GCC**CGCATCGC**C

ATGGACCTCTGCCTTATCACAAGGGCAGACGGCCAGTGCAGACGCTTCTGTATCCCGAGCTCTTGCCC
 AGTGTCGGGAAGAATCGGATCACACGAGGCTTGAAGGATGAATGCAAGGTCCACAGCGTCCAATTCAT
 GGAGTTCTCGAAGTGGAGAACCAGGATGACTACTACAGCACAGAAGCTGGCTGCATCCACACCCAGGAC
 CAGCAAGGGTGCATGTCATGTATGAGGAAGCCTTGAGTGATCTCAAGGAAGTAGAGGCAGAGTTGCTGC
 TTGTGGCTAGCCATTACATTGAGAAAGAAAGGGGTACAAAGTGGGCAGCCAGCCAAAGAGTGGCCAGGA
 GTGGGGATGGGCCATGCAGATGTGGACAGATTGCTGTGCTGTATGACACGTGGACTTGGGAAGCAGCT
 CTTCTGAAAAACAAGCGCCAGCTTCTTGACAGTTACTTTGAAGCCTACCAGCACAGCTTGGACCCTGAGG
 AGAGGTTTACCTTAGCACAAGTGATGACGGACATCGTGCACCGACGCCGAGGTTTGATCTCAGCCACCC
 CTATTTTACCTAAGGCCTACCAAGAGGAATGCACCTGCCTGAGGCTTACCTGCAGTTAGTGAGGGGCATC
 CTAACACGAGATAGAACGTCAGCGGGAGTATGTCCAGAGGCTTGGCGAGAGGATCATCTGGACAATA
 GCAATACATTTGGACTTCCCCTTAATATAGTTTGCAAAACAATTTTCCATCAACAATAGCTGCCCAGC
 TTCAAAAAATGTTTATCTGCTGGAGTTCCACCCTTCCCTCAGCCTGGCAGCACTTATCCCAAAGCCCTA
 GAGCATCTCTTCAAGAGGCCTGCCATGCCTGCAGACCCACATCAGCCAGTGGCCTTGCCTCCCTGGAGC
 GACATGTCCTGCAGCTGGCCCTGGACCTGTGGCTCACCCAGCCACGCTGAGGCCTGGTACAGTACACA
 GCTCCAGAAAGATCTGTTTTAGCTAAAGTTATGGGAGATCCCTTCTCGTGGGGGAGATAGGCCTGCTG
 GCATTCAAGTCTGCAGCAGACGAAGGACAGAAACAAGACCAAGACTCTCACGTGCTGCTCTTGGAGATGT
 TTTCCAACTTCTGGAATCTTGACCTCCGCCACCGCTGATAGAAATGAGCCTGGAGAGTGACACTT
 AGCTCGGCTCTATAAAGAATTGGCATGGGAGATGGGTTTTGAAGAGTTCCATTTGTATCTGAGGCCTGTT
 CACTTTGAATTTGCATCACACGGGGACAGGGCAGACCAGCCACCTCCTGTGTTTATAACCTCTCTGCTGG
 AGGACAGTGGTCGGGTAGACAGATATTCTCCTACTGCTCTTGTCTTAGCCATCTCTGAGGTGGATGATAA
 CCAAAATTGGCAAATTTAGTTTTTACACAAAGGAAGCCATTCTTAAGCTTTTGTGCATTACAGAGTGGAA
 AATATGCAGGTGACTCTTGCATGCCAAGCTGCTCAGAAAAATGCTTTGATGGCGGCTGTCCAGCAGGTGT
 ACTTCTATTGACCCCTCAAGGGAGCTGCTGCTAATGTCAAGGTGAGCCTGCATGACAGAGAAGCCAG
 AGGCATCTCATCCTCCACCAGCTAATGTCAGGCTGCCCTTTCTTCTTCTTTTCCCTTCTCCCT
 TTCTCCCTTTCTTCTCCCTGCCCTACTTTCTTTCTTCTATCTTTATCTTCTTACTTTCTCTTT
 TTCTTCTTTCTTGTTCCTGTTTG

ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT
 ACAAGGATGACGACGATAAGGTTTAA

Protein Sequence:

>RC239017 representing NM_001277339
 Red=Cloning site Green=Tags(s)


MDPLPYHKRRPVQRQPSVSRALAQCPGRIGSHAGLKDECKVHSVQFMEFSEVENQDDYYSTEAGCIHTQD
 QQGVHVMYEEALSDLKELEAELLVASHYIEKERGHKVGSPKSGQEWGAHADVDRAVLVYDTWTWEAA
 LLENKRQLLDYSFEAYQHTLDPEERFTLAQVMTDIVHRRPRFDLSHPYFIKAYQEECTCLRLHLQLVRGI
 LNQQIERQREYVQRLWREDHLDNSNTFGLPLNIVCKQLISINNSCPASKNVYLLFHPSLSLAALIPKAL
 EHLFQEACHACRPTSASGLASLERHVLQLALDLWLTPTAPEAWYSTQLQKDLFSKVMGDPFLVGEIGLL
 AFKSADEGQKQDQDSHVLLLEMFSLLELLTLRHRLIEMSLESAHLARLYKELAWEMGFEEFHLVLRPV
 HFEFASHGDRADQPPPVFITSLLEDSDRVDRYSPTALVLAISEVDNQGKFSFYTKAAILKLLHSGVE
 NMQVTLACQAAQKNALMAAVQVYFYCTPQGSCPANVKVSLHDREARGISSPPPANVRLPLSSFLFSLLP
 FLSFSFPAPYFLFPIISLSSLLSFFLPFLVSCL

TRTRPLE**QKL**ISEEDLAANDILDYKDDDDKV

Restriction Sites:

Sgfl-MluI

Cloning Scheme:



ACCN: NM_001277339

ORF Size: 1776 bp

OTI Disclaimer: The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

OTI Annotation: This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

Reconstitution Method:

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

Note: Plasmids are not sterile. For experiments where strict sterility is required, filtration with 0.22um filter is required.

RefSeq: [NM_001277339.1](#), [NP_001264268.1](#)


RefSeq Size: 1969 bp

RefSeq ORF: 1779 bp

Locus ID: 100996634

Cytogenetics: 6q21
MW: 67.8 kDa

Product images:



Circular map for RC239017